

Ethicon Product Catalog

Sutures • Adhesives • Drains

Table of Contents

General Information.....	4
Needle Index.....	11
True-to-Scale Needle Diagram.....	29
DERMABOND® Portfolio.....	34
STRATAFIX™ Knotless Tissue Control Devices.....	36
Plus Antibacterial Sutures.....	53
Sutures: Absorbable.....	86
Sutures: Nonabsorbable.....	149
Drains.....	223
Primary Suture Storage Carts, Carrying Racks, and Other Wound Closure Devices.....	226
Suture Chart.....	232
Index.....	235

General Information

Explanation of a Typical Selection and Page

The needle type and size/shape are shown only for needle/suture combinations. A double needle image indicates a "double-armed" suture, that is, a suture with a needle attached to each end.

ATRALOC® Surgical Needles

ATRALOC Needles are supplied in a variety of sizes, shapes and chord lengths. Surgeons should select the type of needle that, in their experience, is appropriate for the specific procedure and tissue.

Needle shapes are generally classified according to the degree of curvature of the body – 5/8, 1/2, 3/8, or 1/4 circle and straight – with taper, cutting, blunt points, or TAPERCUT® Surgical Needle tips.

In general, the same size of needle may be made from finer gauge wire for use in soft or delicate tissues and from heavier gauge wire for use in tough or fibrosed tissues (the surgeon's choice).

Needle Type/Needle Shape Symbol

Needle type and needle shape are shown in this subhead.

Material

Suture materials are designated by color and construction characteristics (eg, Undyed Braided, Violet Braided, etc.).

Length

Suture length is shown in both inches and the approximate metric equivalent. Multi-strand (MS) means there is more than 1 strand of suture per package. For example: MS-2 means 2 strands per pack, MS-4 means 4 strands per pack.

The **CR** refers to CONTROL RELEASE™ Needle Sutures.

Product Name

Sutures and other Ethicon products are designated by trademarked and generic names.

Size

Suture size, denoting the diameter of the material, is in accordance with the United States Pharmacopeia (USP).

Needle Silhouette

Needle silhouettes are shown in ACTUAL SIZE. The designated number above the needle, within the Needle section, is the actual length measured in millimeters.

Code

Product codes ending in suffix G are packaged one dozen packets to a box. Product codes ending in suffix T are packaged two dozen packets to a box. Product codes ending in suffix H are packaged three dozen to a box. In the few instances where no alphabetic suffix is in the product code, the number of packets per box is identified by a footnote. The number of strands per multi-strand packet appears on the left of suture length. Primary suture packets do not have suffix G, H, or T.

Needle Length

Needle length is the measurement in millimeters of the length of the needle along the curvature. Ethicon has applied a rounding convention to all needle lengths. Needles less than 10 mm are rounded to the nearest tenth. Needles greater than 10 mm are rounded to the nearest whole number.

Explanation of a Typical Selection and Page

The Major Characteristics of ATRALOC® Surgical Needles are:

- They must be made from high quality stainless steel
- They resist bending but are processed so that they will tend to bend before breaking
- Taper points must be sharp and contoured for easy passage into tissues
- Cutting points or edges must be sharp and free of burrs
- On most ATRALOC Needles, a Super-Smooth finish is provided that permits the needle to penetrate and pass through with minimal resistance or drag

MULTIPASS® Needles


Elements of Needle Design

Needle design involves analyzing a surgical procedure and the density of the tissue involved in great detail. Ethicon engineers work continuously to improve upon their needle line.


The anatomy of the ideal surgical needle includes three factors: Alloy, Geometry, and Coating.

The combination of coating and alloy make up the MULTIPASS Needle Technology, ensuring excellent strength, tissue penetration, and control. MULTIPASS Needles are the highest performing needles that Ethicon offers.


- Ribbed needles – longitudinal ribs are provided on many ATRALOC Needles to increase the stability of the needle in the needleholder
- The attachment of the needle to the suture material must be secure so that the needle will not separate from the suture material under normal use. An exception is the CONTROL RELEASE™ Needle Suture principle, which facilitates fast separation of needle from suture when desirable



CONTROL RELEASE Needle Sutures






Needle is held securely in the holder. Suture is grasped securely just below needle, pulling strand taut.



Needle is released with a straight tug of the needleholder.

Selection and Use of Needleholders



- The needleholder must be an appropriate size instrument for the size needle selected
- It must be made from good quality steel with a secure jaw design
- Needles should be grasped in an area about one-third to one-half the distance from the swaged area to the point. Avoid placement on or near the swaged area
- The needle should be placed securely in the tip of the needleholder jaws
- Do not damage taper points or cutting edges when using the needleholder to pull the needle out through the tissue. Grasp as far back as possible
- When placing the needle in tissue, any force applied should be in the direction following the curve of the needle
- Do not take excessively large bites of tissue with small needles
- Do not force or twist the needle in an effort to bring the point out through the tissue. Withdraw the needle and replace tissue
- Do not force a dull needle through tissue; obtain a new one
- Avoid using the needle to bridge or approximate tissues for suturing
- If the needle is held too tightly in a sharp, hard-jawed, or defective needleholder, the needle may be damaged or notched in such a manner that it will have more of a tendency to bend or break on successive passes through tissue
- When the needle holder with needle is passed to the surgeon, the needle should be pointing in the right direction in which it will start to be used without need for readjustments
- In some patients, the tissues may be tougher or fibrosed more than normal and require the use of a heavier gauge needle
- In deep or confined areas, ideal positioning of the needle may not be possible. The surgeon should proceed with caution or use a heavier gauge needle


Note: The American College of Surgeons recommends universal adoption of blunt-tip suture needles for the closure of fascia and muscle in order to reduce needlestick injuries.

An ATRALOC® Surgical Needle is an instrument that complements the surgeon's skills. The variety of Ethicon needles offers the surgeon a choice of needles that are dependable for consistent quality for practically all surgical applications.


Steps to Create Performance



Conventional Needle



Super Smooth



Needle Swaging

Controlled Manufacturing

Ethicon makes its own needle/suture combinations. Each step, from raw material to finished needle, is carefully controlled for uniformly superior performance.

Strong, Sharp Micrograin Stainless Steel

Ethicon makes its needles from premium grade stainless steel formulated specifically for surgical use.

The Ethicon process controls the orientation of the micrograins so that ductility and strength are closely controlled – for strong needles that start sharp and stay sharp, thus assuring easy passage through tissue.

Super Smooth

Virtually all ATRALOC® Surgical Needles have a micro-thin super-smooth finish that allows the needle to penetrate and pass through the toughest tissue with an almost fluid-like flow, with minimal resistance and drag, giving the surgeon a feeling of control.


Needle Swaging

Needles are permanently attached to the suture materials, eliminating the need for threading and providing a new, sharp needle each time. ATRALOC Needle swaging gives the surgeon a needle with an unbroken surface from end to end; without double strands, it draws through tissue easier with less tissue trauma.

Ribbed for Security and Control

Longitudinal ribbing on the inside and outside curvatures of most needles has a cross-locking action in the needleholder, which minimizes turning, rocking, or twisting, and provides unsurpassed needle control.

Shapes



Types



Conventional Cutting Needle

Two opposing cutting edges, with a third on inside curve. Change in cross-section from a triangle cutting tip to a flattened body.



Reverse Cutting Needle

Cutting edge on outer curve. For tough, difficult-to-penetrate tissues..



Conventional Spatula Needle

Visibility of the point at bottom provides control of depth penetration.



SABRELOC® Spatula Needle

Side-cutting spatula shaped edges. For layers of scleral or corneal tissue. Spatulated needle point is centered for maximum needle stability in thin sclera. Four equidistant and properly defined edges provide greater control.



ETHIGUARD® Blunt Point Needle

Taper body. For blunt dissection and suturing friable tissue.



TAPERCUT® Surgical Needle

Cutting tip, taper body. For tough tissue, like two needles in one.



MICRO-POINT® Reverse Cutting Needle

Cutting edge on outer curve. Extremely smooth. Extremely sharp for ophthalmic surgery.



Taperpoint Needle

For soft, easily penetrated tissues.



MICRO-POINT® Surgical Needle (Spatula)

Thin, flat profile. Specially designed for ophthalmic anterior segment surgery.



ULTIMA® Ophthalmic Needle (Spatula)

Reduced edge-angles provide better penetration. Readily facilitates knot rotation in ophthalmic surgery.



Precision Cosmetic - Conventional Cutting PRIME® Needle

For delicate plastic or cosmetic surgery. Conventional cutting tip and PRIME geometry for increased sharpness in the tip leading to a flattened body.



VISI-BLACK™ Surgical Needle


Slim taper point needles with a black finish for improved visibility.



Precision Point - Reverse Cutting PRIME Needle

For plastic or cosmetic surgery, reverse cutting tip and PRIME geometry for increased sharpness in the tip.

Trademarks

- ATRALOC® Surgical Needle
- CONTROL RELEASE™ Needle or Needle Suture **CR**
- DERMABOND ADVANCED® Topical Skin Adhesive
- DERMABOND® Mini Topical Skin Adhesive
- DERMABOND® PRINEO® Skin Closure System
- ETHALLOY® Needle Alloy
- ETHIBOND® Extra Polyester Suture
- ETHIBOND EXCEL® Polyester Suture
- ETHIGUARD® Blunt Point Needle
- ETHILON® Nylon Suture
- ETHI-PACK™ Pre-cut Sutures
- EVERPOINT® Cardiovascular Needle **EVP**
- HEMO-SEAL™ Needle Suture **HS**
- J-VAC® Closed Wound Drainage System
- J-VAC® Bulb Suction Reservoir
- J-VAC® Reservoir
- J-VAC® Drainage Adapter
- J-VAC® Silicon Drain
- LAPRA-TY® Suture Clip
- LIGAPAK™ Dispensing Reel
- LIGAPAK® Ligature
- MERSILENE® Polyester Fiber Suture
- MERSILENE® Polyester Fiber Strip
- MICRO-POINT® Surgical Needle
- MONOCRYL® (poliglecaprone 25) Suture
- MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture
- MULTIPLASS® Needles 
- NUROLON® Nylon Suture
- P, PS PRIME® Needles
- PC PRIME® Needle
- PDS® (polydioxanone) Suture
- PDS® II (polydioxanone) Suture
- PDS® Plus Antibacterial (polydioxanone) Suture
- PERMA-HAND® Silk Suture
- PROLENE® Polypropylene Suture
- PRONOVA® Poly (Hexafluoropropylene – VDF) Suture
- RELAY Suture Delivery System
- RELAY Suture Package
- SABRELOC® Spatula Needle
- STRATAFIX™ Knotless Tissue Control Devices
- STRATAFIX™ Spiral Knotless Tissue Control Device
- STRATAFIX™ Spiral PDO Knotless Tissue Control Device
- STRATAFIX™ Spiral PGA-PCL Knotless Tissue Control Device
- STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device
- STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device
- SUTUPAK® Pre-cut Sterile Sutures
- TAPERCUT® Surgical Needle
- ULTIMA® Ophthalmic Needles
- VICRYL® (polyglactin 910) Knitted Mesh
- VICRYL® (polyglactin 910) Woven Mesh
- VICRYL® (polyglactin 910) Mesh Bag
- VICRYL® (polyglactin 910) Perirenal Mesh Bag
- VICRYL® (polyglactin 910) Perisplenic Mesh Bag
- VICRYL® (polyglactin 910) Perihepatic Mesh Bag
- VICRYL® (polyglactin 910) Perihepatic Right Lobe Mesh Bag
- VICRYL® (polyglactin 910) Perihepatic Left Lobe Mesh Bag
- Coated VICRYL® (polyglactin 910) Suture
- Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture
- VICRYL RAPIDE™ (polyglactin 910) Suture
- VISI-BLACK™ Surgical Needle **VB**

Needle Index



1/4 Circle

TAPERCUT® Surgical Needle



Needles	Needle Length	Available Products
V75-3	3.5 mm	ETHILON
V100-3	3.5 mm	ETHILON
CIF-4	13 mm	PROLENE

SABRELOC® Spatula Needle



Needles	Needle Length	Available Products
S-4	8 mm	GUT CHROMIC, MERSILENE
S-29	7.6 mm	COATED VICRYL
S-24	8 mm	MERSILENE
S-14	8 mm	COATED VICRYL, GUT CHROMIC, GUT PLAIN, ETHILON, MERSILENE, PERMA HAND SILK

1/4 Circle

ULTIMA® Ophthalmic Needle



Needles	Needle Length	Available Products
CS90-6	7 mm	VICRYL
CTC-6L	13 mm	PROLENE
CTC-6	13 mm	PROLENE

3/8 Circle

Taperpoint



Needles	Needle Length	Available Products
TE	32 mm	TEMPORARY PACING WIRE
BB-1	22 mm	TEMPORARY PACING WIRE
BB	17 mm	STRATAFIX SPIRAL MONOCRYL PLUS, GUT CHROMIC, ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK, TEMPORARY PACING WIRE
C-1	13 mm	PDS PLUS, GUT CHROMIC, PDS II, ETHIBOND, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK
BV	11 mm	PDS PLUS, PDS II, PROLENE, PRONOVA
BV-1	9.3 mm	PDS PLUS, PDS II, ETHIBOND, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK
BV175-8	9.3 mm	PROLENE, PRONOVA
BV175-7	8 mm	PROLENE, PRONOVA

3/8 Circle

Taperpoint



Needles	Needle Length	Available Products
BV175-6	8 mm	PROLENE, PRONOVA
BV130-5	6.5 mm	COATED VICRYL, ETHILON, PROLENE, PRONOVA
BV130-4	5 mm	COATED VICRYL, ETHILON
BV130-3	3.8 mm	ETHILON
BV100-4	5.1 mm	ETHILON, PROLENE
BV100-3	3.8 mm	VICRYL
BV75-4	5.1 mm	ETHILON
BV75-3	3.8 mm	ETHILON, PROLENE
BV50-3	3.8 mm	ETHILON


TAPERCUT® Surgical Needle



Needles	Needle Length	Available Products
V-26	26 mm	STRATAFIX SPIRAL PDO
V-4	17 mm	STRATAFIX SPIRAL PDO, ETHIBOND, PROLENE, PERMA HAND SILK
CC-1	13 mm	ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK
CC	9.3 mm	PROLENE
V130-5	5.8 mm	COATED VICRYL
V130-4	4.2 mm	COATED VICRYL
V130-3	3.8 mm	ETHILON

3/8 Circle

Reverse Cutting



Needles	Needle Length	Available Products
LR	75 mm	ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
FSLX	36 mm	ETHILON, PROLENE, STEEL, PERMA HAND SILK
FSL	30 mm	COATED VICRYL PLUS, COATED VICRYL, ETHILON, NUROLON, PROLENE, PERMA HAND SILK
FS	26 mm	STRATAFIX SPIRAL PDS PLUS, STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, COATED VICRYL, ETHIBOND, ETHILON, MERSILENE, NUROLON, PROLENE, PERMA HAND SILK
FS-1	24 mm	STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
FS-2	19 mm	STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, STEEL, PERMA HAND SILK
M-2	13 mm	GUT CHROMIC
C-3	13 mm	COATED VICRYL, GUT CHROMIC, ETHILON, PERMA HAND SILK
C-2	13 mm	ETHILON, PERMA HAND SILK

3/8 Circle

Conventional Cutting



Needles	Needle Length	Available Products
PC-5	19 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PC-3	16 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PC-1	13 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHILON, NUROLON, PROLENE

Precision Point Reverse Cutting



Needles	Needle Length	Available Products
PSLX	36 mm	ETHILON, PERMA HAND SILK
PSL	30 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, COATED VICRYL, ETHILON, PERMA HAND SILK
PS	26 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, GUT CHROMIC, ETHILON, PERMA HAND SILK
PS-1	24 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHILON, PROLENE, PERMA HAND SILK
PS-2	19 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK

3/8 Circle

Precision Point Reverse Cutting



Needles	Needle Length	Available Products
PS-3	16 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
P-3	13 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, MERSILENE, PROLENE, PERMA HAND SILK
P-1	11 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
P-6	8 mm	ETHILON, PROLENE

ULTIMA® Ophthalmic Needle



Needles	Needle Length	Available Products
CS140-8	6.5 mm	ETHILON
CS140-6	6.5 mm	VICRYL, ETHILON, MERSILENE, PROLENE


MICRO-POINT® Surgical Needle (Spatula)



Needles	Needle Length	Available Products
TG140-8	6.5 mm	COATED VICRYL, VICRYL, GUT CHROMIC, GUT PLAIN, ETHILON, PROLENE, PERMA HAND SILK
TG140-6	6.5 mm	COATED VICRYL, ETHILON
TG100-8	6.5 mm	PROLENE
VAS100-4	5 mm	ETHILON

1/2 Circle

Blunt Point



Needles	Needle Length	Available Products
BP	80 mm	SURGICAL GUT SUTURE - CHROMIC MERSILENE SUTURE
BP-1	65 mm	MERSILENE SUTURE
BP-2	48 mm	SURGICAL GUT SUTURE - CHROMIC
TPB-1	65 mm	COATED VICRYL PLUS SUTURE PDS PLUS SUTURE
CTXB	48 mm	STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, PROLENE
CTB	40 mm	COATED VICRYL PLUS, PDS PLUS
CTB-1	36 mm	STRATAFIX SYMMETRIC PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, PROLENE
CTB-2	26 mm	COATED VICRYL PLUS, PDS PLUS
SHB	26 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS COATED VICRYL, MONOCRYL

1/2 Circle

Taperpoint



Needles	Needle Length	Available Products
XLH	70 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, PROLENE
LH	40 mm	GUT CHROMIC
MH	36 mm	STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL, ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK
SH	26 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, MERSILENE, NUROLON, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK, TEMPORARY PACING WIRE
SH-1	22 mm	STRATAFIX SPIRAL POLYPROPYLENE, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL, PDS II, ETHIBOND, NUROLON, PROLENE, PRONOVA, PERMA HAND SILK
SH-2	20 mm	ETHIBOND, PROLENE HS

1/2 Circle

Taperpoint



Needles	Needle Length	Available Products
RB-1	17 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, NUROLON, PROLENE, PROLENE HS, PRONOVA, PERMA HAND SILK, TEMPORARY PACING WIRE
RB-2	13 mm	COATED VICRYL PLUS, PDS PLUS, PDS II, PROLENE, PROLENE HS, PROVONA, PERMA HAND SILK
UCL	13 mm	PDS PLUS, PDS II, ETHIBOND
TF	13 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, MONOCRYL, PDS II, ETHIBOND, NUROLON, PROLENE, PERMA HAND SILK
RB-3	11 mm	PROLENE
BVH100-3	5 mm	ETHILON

1/2 Circle

TAPERCUT® Surgical Needle



Needles	Needle Length	Available Products
V-56	60 mm	STEEL
V-40	48 mm	ETHIBOND, STEEL
V-37	40 mm	ETHIBOND, STEEL
V-34	36 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, COATED VICRYL, VICRYL RAPIDE, GUT CHROMIC, GUT PLAIN, MONOCRYL, ETHIBOND, PROLENE
V-7	26 mm	STRATAFIX SPIRAL PDS PLUS, GUT PLAIN, ETHIBOND, PROLENE, PRONOVA, PERMA HAND SILK
V-5	17 mm	PDS PLUS, GUT CHROMIC, PDS II, ETHIBOND, PROLENE, PERMA HAND SILK

1/2 Circle

Reverse Cutting



Needles	Needle Length	Available Products
LS-1	45 mm	STEEL
CPX	48 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II
CP	40 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, PDS II, NUROLON
CP-1	36 mm	STRATAFIX SPIRAL MONOCRYL PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, STEEL, PERMA HAND SILK
CP-2	26 mm	STRATAFIX SPIRAL MONOCRYL PLUS, STRATAFIX SPIRAL PDS PLUS, STRATAFIX SPIRAL PDO, STRATAFIX SPIRAL PGA-PCL, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, PDS II, NUROLON
OS-8	40 mm	STRATAFIX SPIRAL PDS PLUS, STRATAFIX SPIRAL PDO, COATED VICRYL PLUS, COATED VICRYL, ETHIBOND, NUROLON
OS-6	36 mm	STRATAFIX SPIRAL PDS PLUS, STRATAFIX SYMMETRIC PDS PLUS, STRATAFIX SPIRAL PDO, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II, ETHIBOND
OS-4	22 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, PDS II, ETHIBOND
OS-2	17 mm	ETHIBOND
X-1	22 mm	COATED VICRYL PLUS, COATED VICRYL, GUT PLAIN, ETHILON, NUROLON, PERMA HAND SILK

1/2 Circle

Reverse Cutting



Needles	Needle Length	Available Products
J-1	17 mm	COATED VICRYL, GUT CHROMIC, PERMA HAND SILK
M-1	13 mm	GUT CHROMIC, GUT PLAIN
PS-4C	16 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC
PS-4	16 mm	STRATAFIX SPIRAL MONOCRYL PLUS, COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II, ETHIBOND, ETHILON, PROLENE, PERMA HAND SILK
PS-5	13 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, PROLENE
PS-6	10 mm	COATED VICRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
P-2	8 mm	COATED VICRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, PDS II
G-2	8 mm	GUT CHROMIC
G-7	8 mm	PERMA HAND SILK

Conventional Cutting



Needles	Needle Length	Available Products
PC-12	16 mm	COATED VICRYL, MONOCRYL, ETHILON

SABRELOC® Spatula Needle



Needles	Needle Length	Available Products
S-28	7.6 mm	COATED VICRYL
S-2	8 mm	GUT CHROMIC, MERSILENE
S-22	8 mm	ETHILON

1/2 Circle

ULTIMA® Ophthalmic Needle



Needles	Needle Length	Available Products
CS175-8	7 mm	GUT CHROMIC, ETHILON
CS160-8	5.5 mm	VICRYL, ETHILON, PROLENE
CS175-6	7 mm	ETHILON
CS160-6	5.5 mm	VICRYL, ETHILON, MERSILENE, PROLENE
CSM-6	4.3 mm	ETHILON

MICRO-POINT® Surgical Needle (Spatula)



Needles	Needle Length	Available Products
TG175-8	7 mm	ETHILON
TG160-8	5.5 mm	COATED VICRYL
TG160-6	5.5 mm	COATED VICRYL, ETHILON
TG160-4	5.3 mm	ETHILON, MERSILENE

5/8 Circle

Taperpoint



Needles	Needle Length	Available Products
UR-4	40 mm	COATED VICRYL, GUT CHROMIC
UR-5	36 mm	COATED VICRYL PLUS, MONOCRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
UR-6	26 mm	STRATAFIX SPIRAL PDS PLUS, COATED VICRYL PLUS, MONOCRYL PLUS, COATED VICRYL, GUT CHROMIC, MONOCRYL
TF-5	8 mm	PERMA HAND SILK

Straight

Taperpoint



Needles	Needle Lengths	Available Products
ST-1	45 mm	PERMA HAND SILK
ST-3	22 mm	COATED VICRYL, ETHIBOND, PDS II
ST-4	19 mm	COATED VICRYL, ETHIBOND, PDS II



ST100-3	3 mm	ETHILON
---------	------	---------

TAPERCUT® Surgical Needle



Needles	Needle Lengths	Available Products
V-47	36 mm	PROLENE

Straight Cutting



Needles	Needle Lengths	Available Products
SKS-3	60 mm	TEMPORARY PACING WIRE



SKS	60 mm	TEMPORARY PACING WIRE
-----	-------	-----------------------



KS	60 mm	COATED VICRYL PLUS, MONOCRYL PLUS, PDS PLUS, COATED VICRYL, GUT CHROMIC, GUT PLAIN, MONOCRYL, PDS II, ETHIBOND, ETHILON, PROLENE, STEEL, PERMA HAND SILK
----	-------	--



TS	48 mm	PROLENE
----	-------	---------



SC-1	13 mm	GUT PLAIN
------	-------	-----------

SABRELOC® Spatula Needle



Needles	Needle Lengths	Available Products
STC-6	16 mm	PROLENE

Trocar

Trocar

Needles	Needle Lengths	Available Products
STP-10	254 mm	PDS II, ETHIBOND

Compound Curve

ULTIMA® Ophthalmic Needle



Needles	Needle Lengths	Available Products
CSC-6	5 mm	ETHILON, PROLENE

Bi-Curve


ULTIMA® Ophthalmic Needle



Needles	Needle Lengths	Available Products
CSB-6	5 mm	VICRYL, ETHILON

Sternum

Conventional Cutting - Sternum



Needles	Needle Lengths	Available Products
CCS-1	60 mm	COATED VICRYL, PDS II, STEEL

CCS	48 mm	ETHIBOND, STEEL
-----	-------	-----------------

Center Point Spatula

SABRELOC® Spatula Needle



Needles	Needle Lengths	Available Products
RD-1	8 mm	COATED VICRYL, ETHILON, MERSILINE

Endoscopy

Taperpoint



Needles	Needle Lengths	Available Products
EN-S	19 mm	PDS
V-38	30 mm	PDS II

True-to-Scale Needle Diagram



True-to-Scale Needle Diagram


Cutting Edge

PRECISION POINT

Reverse Cutting

3/8 Circle

1/2 Circle




REVERSE CUTTING

3/8 Circle



1/2 Circle



PRECISION COSMETIC 

Conventional Cosmetic PRIME®

3/8 Circle

1/2 Circle



MICRO-POINT® 

ULTIMA® 

Reverse Cutting

3/8 Circle

1/2 Circle


Spatula

1/8 Circle

1/4 Circle

3/8 Circle

1/2 Circle



ULTIMA® 

SPATULA 

MICRO-POINT® SPATULA 

Compound Curve



1/4 Circle

3/8 Circle

1/2 Circle

3/8 Circle

1/2 Circle



SABRELOC® 

Conventional Spatula

1/4 Circle

1/2 Circle

Spatula

1/8 Circle

1/4 Circle

1/2 Circle

Center Point Spatula

1/4 Circle



Straight Transchamber

STC-6

STRAIGHT CUTTING 

SC-1

SKS

TS

SKS-3

KS

STERNUM 



Taperpoint




3/8 Circle



5/8 Circle

1/2 Circle



STRAIGHT 



TAPERCUT®



3/8 Circle

1/4 Circle

Straight

1/2 Circle

V130-3
V130-4
V130-5
CC
CC-1

CC175-6
CC175-8
V-4
V-26

V75-3
V100-3
CIF-4

V-47

V-5
V-7
V-34
V-37
V-40
V-56

Trocar

STP-10

Blunt Point



1/2 Circle

ETHIGUARD®

BP-2
BP-1
BP

SHB
CTB-2
CTB-1

CTB
CTXB
TPB-1

Endoscopy

Taperpoint

TAPERCUT®

EN-S
EN-1
EN-2

EN-5
DT-6
EN-3
EN-4
V-38

DERMABOND®

Portfolio

DERMABOND® Portfolio

DERMABOND ADVANCED® Topical Skin Adhesive



Product Code	Description	Adhesive Volume	EA/BX
DNX12	2-octyl cyanoacrylate (2-OCA) topical skin adhesive with an easy-to-use applicator for precise control for wide and fine line applications	0.7 mL	12

DNX6	2-OCA topical skin adhesive with an easy-to-use applicator for precise control for wide and fine line applications	0.7 mL	6
------	--	--------	---

DERMABOND® PRINEO® Skin Closure System



Product Code	Description	Adhesive Volume	EA/BX
CLR222US	2-OCA topical skin adhesive and 22 cm self-adhering mesh designed for optimized healing and high patient satisfaction ¹	3.8 mL	1 x 22 cm mesh patch

CLR422US	2-OCA topical skin adhesive and 42 cm self-adhering mesh designed for optimized healing and high patient satisfaction ¹	3.8 mL	2 x 22 cm mesh patches
----------	--	--------	------------------------

CLR602US	2-OCA topical skin adhesive and 60 cm self-adhering mesh designed for optimized healing and high patient satisfaction ¹	3.8 mL	1 x 60 cm mesh dispenser
----------	--	--------	--------------------------

DERMABOND® Mini Topical Skin Adhesive



Product Code	Description	Adhesive Volume	EA/BX
DHVM12	2-OCA high viscosity formulation for small incisions and lacerations	0.36 mL	12

¹Double-blinded quantitative market research study comparing surgeon experience with DERMABOND PRINEO System and skin staples in total knee arthroplasty. N=83 orthopaedic surgeons. 90% c.i. Fielded June/July 2017.

Reference: 1. PRINEO Claims Research Quant Detail. August 16, 2017. Ethicon, Inc.


STRATAFIX™
Knotless Tissue
Control Devices



STRATAFIX™ Knotless Tissue Control Devices

STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Devices are the only knotless tissue control device that provides strong, secure closure appropriate for high-tension areas such as fascia.¹⁻⁶






STRATAFIX™ Spiral Knotless Tissue Control Devices provide smooth tissue passage and a secure hold that helps control tension and achieve excellent tissue approximation.⁷⁻¹¹



References: **1.** 100326296: Time Zero Tissue Holding - Competitive Claims Comparisons for STRATAFIX Knotless Tissue Control Devices vs Various Products. 2015. Ethicon, Inc. **2.** Ethicon study AST-2011-0210. Study to evaluate the tissue holding performance at time zero of DOLFIN PDS PLUS barbed suture sizes 1 and 2-0 vs dyed PDS II Plus suture sizes 1 and 2-0 in a continuous stitch pattern—Project DOLFIN 11822, version 1. Approved on July 15, 2011. Ethicon, Inc. **3.** Ethicon study AST-2011-0341. Performance testing of DOLFIN PDS PLUS size 3-0 suture—tissue holding 10 cm incision. Approved on August 22, 2011. Ethicon, Inc. **4.** Ethicon study PSE 09-0204, project number 711822. Exploratory histological and biomechanical evaluation of DOLFIN following closure of the ventral abdominal wall in a porcine model at ± 1 days. Approved on July 7, 2010. Ethicon, Inc. **5.** Ethicon study PSE 10-0012, project 7: number 11822. Model development: histological and biomechanical evaluation of 3-0 DOLFIN barbed suture prototypes, 3-0 Quill suture, and 3-0 V-loc suture at 1 days following closure of the ventral abdominal wall in a rabbit model. Approved on August 8, 2011. Ethicon, Inc. **6.** Performance Testing of STRATAFIX SYMMETRIC PDS PLUS Size 0 & 1 Devices - Initiation Strength in Porcine Tissue. AST-2013-0603. Ethicon, Inc. Moran ME, Marsh C, Perrotti M. Bidirectional-barbed sutured knotless running anastomosis v classic Van Velthoven suturing in a model system. *J Endourol.* 2007;21(10):1175-1178. **7.** Vakili JJ, O'Reilly MP, Sutter EG, Mears SC, Belkoff SM, Khanuja HS. Knee arthroscopy repair with a continuous barbed suture: a biomechanical study. *J Arthroplasty.* 2011;26(5):710-713. **8.** Eickmann T, Quane E. Total knee arthroplasty closure with barbed sutures. *J Knee Surg.* 2010;23(3):163-167. **9.** Einarsson JJ, Chavan NR, Suzuki Y, Jonsdottir G, Vellinga TT, Greenberg JA. Use of bidirectional barbed suture in laparoscopic myomectomy: evaluation of perioperative outcomes, safety, and efficacy. *J Minim Invasive Gynecol.* 2011;18(1):92-95. **10.** Levine BR, Ting N, Della Valle CJ. Use of a barbed suture in the closure of hip and knee arthroplasty wounds. *Orthopedics.* 2011;34(9):e473-e475. doi: 10.3928/01477447-20110714-35. **11.** Ethicon study AST-2012-0331. Tissue gapping under tension of porcine cadaveric skin incisions closed with Stratafix Spiral in comparison to Monocryl in both interrupted and continuous stitching patterns. Approved on August 24, 2012. Ethicon, Inc.



STRATAFIX™ Spiral PDO Knotless Tissue Control Device





Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 CT-3 22mm 1/2 Circle	0	Violet	5"/14cm x 5"/14cm	12	SXPD2B404
 SH 26mm 1/2 Circle	2-0	Violet	4"/10cm x 4"/10cm	12	SXPD2B413
	2-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B414
 MH 36mm 1/2 Circle	2-0	Violet	3"/7cm x 3"/7cm	12	SXPD2B406
	2-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B407
	2-0	Violet	9"/23cm x 9"/23cm	12	SXPD2B408
	2-0	Violet	14"/35cm x 14"/35cm	12	SXPD2B412
	0	Violet	5"/14cm x 5"/14cm	12	SXPD2B409
	0	Violet	9"/23cm x 9"/23cm	12	SXPD2B410
	0	Violet	12"/30cm x 12"/30cm	12	SXPD2B411
 CT-1 36mm 1/2 Circle	1	Violet	5"/14cm x 5"/14cm	12	SXPD2B401
	1	Violet	9"/23cm x 9"/23cm	12	SXPD2B402
	1	Violet	12"/30cm x 12"/30cm	12	SXPD2B403

STRATAFIX™ Spiral PDO Knotless Tissue Control Device








Taperpoint					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>MO-4 36mm 1/2 Circle</p>	1	Violet	14"/35cm x 14"/35cm	12	SXPD2B400
 <p>CTX 48mm 1/2 Circle</p>	1	Violet	14"/35cm x 14"/35cm	12	SXPD2B405

TAPERCUT®					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>V-4 17mm 3/8 Circle</p>	4-0	Violet	9"/23cm x 9"/23cm	12	SXPD2B426
 <p>V-26 26mm 3/8 Circle</p>	4-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B425
	4-0	Violet	9"/23cm x 9"/23cm	12	SXPD2B424

STRATAFIX™ Spiral PDO Knotless Tissue Control Device









Cutting Edge Reverse 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 FS-2 19mm 3/8 Circle	5-0	Violet	3"/7cm x 3"/7cm	12	SXPD2B421
	4-0	Violet	3"/7cm x 3"/7cm	12	SXPD2B422
	4-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B423
 FS 26mm 3/8 Circle	3-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B417
	3-0	Violet	9"/23cm x 9"/23cm	12	SXPD2B419
	2-0	Violet	5"/14cm x 5"/14cm	12	SXPD2B418
	2-0	Violet	9"/23cm x 9"/23cm	12	SXPD2B420
 CP-2 26mm 1/2 Circle	0	Violet	9"/23cm x 9"/23cm	12	SXPD2B415
	1	Violet	9"/23cm x 9"/23cm	12	SXPD2B416
 OS-6 36mm 1/2 Circle	1	Violet	9"/23cm x 9"/23cm	12	SXPD2B200
	1	Violet	12"/30cm x 12"/30cm	12	SXPD2B201
 OS-8 40mm 1/2 Circle	1	Violet	14"/35cm x 14"/35cm	12	SXPD2B202

STRATAFIX™ Spiral PGA-PCL Knotless Tissue Control Device



Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	6"/15cm x 6"/15cm	12	SXMD2B403
	3-0	Undyed	6"/15cm x 6"/15cm	12	SXMD2B402
 MH 36mm 1/2 Circle	2-0	Undyed	5"/14cm x 5"/14cm	12	SXMD2B400
	2-0	Undyed	14"/35cm x 14"/35cm	12	SXMD2B401


Cutting Edge Reverse 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 FS-2 19mm 3/8 Circle	4-0	Undyed	3"/7cm x 3"/7cm	12	SXMD2B405
	4-0	Undyed	5"/14cm x 5"/14cm	12	SXMD2B407
	4-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B409
	3-0	Undyed	3"/7cm x 3"/7cm	12	SXMD2B404
	3-0	Undyed	5"/14cm x 5"/14cm	12	SXMD2B406
	3-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B408
 FS-1 24mm 3/8 Circle	4-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B150
	3-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B410

STRATAFIX™ Spiral PGA-PCL Knotless Tissue Control Device




Cutting Edge Reverse ▾					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 FS 26mm 3/8 Circle	4-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B413
	3-0	Undyed	5"/14cm x 5"/14cm	12	SXMD2B411
	3-0	Undyed	12"/30cm x 12"/30cm	12	SXMD2B412
 CP-2 26mm 1/2 Circle	2-0	Undyed	5"/14cm x 5"/14cm	12	SXMD2B414

STRATAFIX™ Spiral Polypropylene Knotless Tissue Control Device



Taperpoint ⊙					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 CT-3 22mm 1/2 Circle	0	Undyed	5"/14cm x 5"/14cm	12	SXPL2B401
 CT-1 36mm 1/2 Circle	1	Undyed	9"/23cm x 9"/23cm	12	SXPL2B400
 SH-1 22mm 1/2 Circle	2-0	Undyed	12"/30cm	12	SXPL1B400




STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device





Taperpoint					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	6"/15cm	12	SXMP1B434
	4-0	Undyed	12"/30cm	12	SXMP1B435
	3-0	Undyed	6"/15cm	12	SXMP1B424
	3-0	Undyed	9"/23cm	12	SMMP1B438
	3-0	Undyed	12"/30cm	12	SXMP1B425
 CT-2 26mm 1/2 Circle	2-0	Undyed	6"/15cm	12	SXMP1B415
	2-0	Undyed	12"/30cm	12	SXMP1B416
	2-0	Undyed	18"/45cm	12	SXMP1B417
 SH 26mm 1/2 Circle	4-0	Undyed	27"/70cm	12	SXMP1B437
	3-0	Undyed	9"/23cm	12	SXMP1B427
	3-0	Undyed	27"/70cm	12	SXMP1B428
	2-0	Undyed	6"/15cm	12	SXMP1B408
	2-0	Undyed	9"/23cm	12	SXMP1B409
	2-0	Undyed	12"/30cm	12	SXMP1B410
 CT-1 36mm 1/2 Circle	3-0	Undyed	18"/45cm	12	SXMP1B429
	2-0	Undyed	6"/15cm	12	SXMP1B411
	2-0	Undyed	12"/30cm	12	SXMP1B412
	2-0	Undyed	18"/45cm	12	SXMP1B413
	2-0	Undyed	27"/70cm	12	SXMP1B414
 ST-4 19mm Straight	3-0	Undyed	12"/30cm	12	SXMP1B433

STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device





Cutting Edge Reverse ▾					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 PSL 30mm 3/8 Circle	2-0	Undyed	27"/70cm	12	SXMP1B421
 CP-2 26mm 1/2 Circle	2-0	Undyed	12"/30cm	12	SXMP1B420
 CP-1 36mm 1/2 Circle	2-0	Undyed	27"/70cm	12	SXMP1B419

PRIME Reverse Cutting ▾					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS-2 19mm 3/8 Circle	4-0	Undyed	12"/30cm	12	SXMP1B117
	4-0	Undyed	18"/45cm	12	SXMP1B118
	4-0	Undyed	27"/70cm	12	SXMP1B119
	3-0	Undyed	6"/15cm	12	SXMP1B105
	3-0	Undyed	12"/30cm	12	SXMP1B106
	3-0	Undyed	18"/45cm	12	SXMP1B107
	3-0	Undyed	24"/60cm	12	SXMP1B108
	3-0	Undyed	27"/70cm	12	SXMP1B109
 PS-1 24mm 3/8 Circle	4-0	Undyed	12"/30cm	12	SXMP1B114
	4-0	Undyed	18"/45cm	12	SXMP1B115
	4-0	Undyed	27"/70cm	12	SXMP1B116
	3-0	Undyed	12"/30cm	12	SXMP1B101
	3-0	Undyed	18"/45cm	12	SXMP1B102
	3-0	Undyed	24"/60cm	12	SXMP1B103
	3-0	Undyed	27"/70cm	12	SXMP1B104




STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device



PRIME Reverse Cutting ▽					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS 26mm 3/8 Circle	3-0	Undyed	12"/30cm	12	SXMP1B110
	3-0	Undyed	18"/45cm	12	SXMP1B111
	3-0	Undyed	27"/70cm	12	SXMP1B113
 PS-4 16mm 1/2 Circle	4-0	Undyed	18"/45cm	12	SXMP1B120









STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device



Taperpoint ⊙					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	4-0	Violet	6"/15cm	12	SXPP1B427
	4-0	Violet	12"/30cm	12	SXPP1B428
	3-0	Violet	9"/23cm	12	SXMP1B438
	3-0	Violet	6"/15cm	12	SXPP1B422
 CT-2 26mm 1/2 Circle	3-0	Violet	27"/70cm	12	SXPP1B425
	2-0	Violet	6"/15cm	12	SXPP1B413
	2-0	Violet	9"/23cm	12	SXPP1B432
	2-0	Violet	27"/70cm	12	SXPP1B414
	0	Violet	12"/30cm	12	SXPP1B405
 SH 26mm 1/2 Circle	4-0	Violet	27"/70cm	12	SXPP1B426
	3-0	Violet	6"/15cm	12	SXPP1B420
	3-0	Violet	27"/70cm	12	SXPP1B421
	2-0	Violet	6"/15cm	12	SXPP1B415
	2-0	Violet	9"/23cm	12	SXPP1B433
	2-0	Violet	12"/30cm	12	SXPP1B416
	2-0	Violet	27"/70cm	12	SXPP1B417




STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device









Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 SH-1 22mm 1/2 Circle	4-0	Violet	9"/23cm	12	SXPP1B454
	3-0	Violet	9"/23cm	12	SXPP1B453
 CT-1 36mm 1/2 Circle	2-0	Violet	6"/15cm	12	SXPP1B409
	2-0	Violet	9"/23cm	12	SXPP1B456
	2-0	Violet	12"/30cm	12	SXPP1B410
	2-0	Violet	18"/45cm	12	SXPP1B411
	2-0	Violet	27"/70cm	12	SXPP1B412
	0	Violet	6"/15cm	12	SXPP1B406
	0	Violet	9"/23cm	12	SXPP1B455
	0	Violet	12"/30cm	12	SXPP1B450
	0	Violet	18"/45cm	12	SXPP1B407
	0	Violet	27"/70cm	12	SXPP1B408
	1	Violet	12"/30cm	12	SXPP1B429
	1	Violet	18"/45cm	12	SXPP1B430
1	Violet	27"/70cm	12	SXPP1B431	
 MO-4 36mm 1/2 Circle	1	Violet	27"/70cm	12	SXPP1B205
 CTX 48mm 1/2 Circle	1	Violet	27"/70cm	12	SXPP1B401
	1	Violet	36"/90cm	12	SXPP1B402
 UR-6 26mm 5/8 Circle	2-0	Violet	12"/30cm	12	SXPP1B419
 EN-S 19mm Ski	0	Violet	12"/30cm	12	SXPP1B452
 ST-4 19mm Straight	3-0	Violet	6"/15cm	12	SXPP1B424

STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device





TAPERCUT® 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 V-7 26mm 1/2 Circle	2-0	Violet	12"/30cm	12	SXPP1B418
 V-34 36mm 1/2 Circle	2-0	Violet	9"/23cm	12	SXPP1B458
	0	Violet	9"/23cm	12	SXPP1B457
	0	Violet	12"/30cm	12	SXPP1B404


Cutting Edge Reverse 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 FS 26mm 3/8 Circle	2-0	Violet	18"/45cm	12	SXPP1B101
	2-0	Violet	27"/70cm	12	SXPP1B102
 CP-2 26mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1B202
	1	Violet	18"/45cm	12	SXPP1B403
 OS-6 36mm 1/2 Circle	1	Violet	18"/45cm	12	SXPP1B201
	1	Violet	27"/70cm	12	SXPP1B203
 OS-8 40mm 1/2 Circle	1	Violet	27"/70cm	12	SXPP1B204
 X-1 22mm 1/2 Circle	3-0	Violet	27"/70cm	12	SXPP1B120




STRATAFIX™ Spiral PDS™ Plus Knotless Tissue Control Device




PRIME Reverse Cutting 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS-2 19mm 3/8 Circle	4-0	Violet	6"/15cm	12	SXPP1B111
	4-0	Violet	12"/30cm	12	SXPP1B112
	4-0	Violet	18"/45cm	12	SXPP1B113
	4-0	Violet	27"/70cm	12	SXPP1B114
	3-0	Violet	6"/15cm	12	SXPP1B103
	3-0	Violet	12"/30cm	12	SXPP1B104
	3-0	Violet	18"/45cm	12	SXPP1B105
	3-0	Violet	27"/70cm	12	SXPP1B106
 PS-1 24mm 3/8 Circle	4-0	Violet	18"/45cm	12	SXPP1B116
	4-0	Violet	23"/60cm	12	SXPP1B117
	3-0	Violet	27"/70cm	12	SXPP1B110
 PS 26mm 3/8 Circle	4-0	Violet	18"/45cm	12	SXPP1B118
	3-0	Violet	12"/30cm	12	SXPP1B107
	3-0	Violet	18"/45cm	12	SXPP1B108
	3-0	Violet	27"/70cm	12	SXPP1B109






STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device




ETHIGUARD® Blunt Point Needle 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 CTB-1 36mm 1/2 Circle	0	Violet	23"/60cm	12	SXPP1A303
	1	Violet	18"/45cm	12	SXPP1A301
	1	Violet	23"/60cm	12	SXPP1A304
 CTXB 48mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1A302
	0	Violet	23"/60cm	12	SXPP1A305
	1	Violet	18"/45cm	12	SXPP1A300
	1	Violet	23"/60cm	12	SXPP1A306






STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device




Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 RB-1 17mm 1/2 Circle	3-0	Violet	12"/30cm	12	SXPP1A427
 SH-1 22mm 1/2 Circle	2-0	Violet	12"/30cm	12	SXPP1A430
 CT-2 26mm 1/2 Circle	2-0	Violet	6"/15cm	12	SXPP1A415
	2-0	Violet	9"/23cm	12	SXPP1A422
	2-0	Violet	18"/45cm	12	SXPP1A408
	0	Violet	9"/23cm	12	SXPP1A446
	0	Violet	18"/45cm	12	SXPP1A407
	1	Violet	6"/15cm	12	SXPP1A419
 SH 26mm 1/2 Circle	3-0	Violet	6"/15cm	12	SXPP1A413
	3-0	Violet	12"/30cm	12	SXPP1A428
	3-0	Violet	18"/45cm	12	SXPP1A410
	2-0	Violet	6"/15cm	12	SXPP1A414
	2-0	Violet	9"/23cm	12	SXPP1A421
	2-0	Violet	12"/30cm	12	SXPP1A431
	2-0	Violet	18"/45cm	12	SXPP1A409




STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device



Taperpoint 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Violet	12"/30cm	12	SXPP1A429
	3-0	Violet	18"/45cm	12	SXPP1A438
	2-0	Violet	6"/15cm	12	SXPP1A416
	2-0	Violet	9"/23cm	12	SXPP1A423
	2-0	Violet	12"/30cm	12	SXPP1A432
	2-0	Violet	18"/45cm	12	SXPP1A403
	2-0	Violet	24"/60cm	12	SXPP1A439
	0	Violet	6"/15cm	12	SXPP1A418
	0	Violet	9"/23cm	12	SXPP1A425
	0	Violet	12"/30cm	12	SXPP1A433
	0	Violet	18"/45cm	12	SXPP1A401
	0	Violet	24"/60cm	12	SXPP1A440
	1	Violet	6"/15cm	12	SXPP1A420
	1	Violet	12"/30cm	12	SXPP1A435
1	Violet	18"/45cm	12	SXPP1A404	
1	Violet	24"/60cm	12	SXPP1A443	
 <p>CT 40mm 1/2 Circle</p>	0	Violet	12"/30cm	12	SXPP1A412
	0	Violet	18"/45cm	12	SXPP1A406
	0	Violet	24"/60cm	12	SXPP1A441
	1	Violet	18"/45cm	12	SXPP1A405
	1	Violet	24"/60cm	12	SXPP1A444
 <p>CTX 48mm 1/2 Circle</p>	0	Violet	12"/30cm	12	SXPP1A434
	0	Violet	18"/45cm	12	SXPP1A402
	0	Violet	24"/60cm	12	SXPP1A442
	1	Violet	12"/30cm	12	SXPP1A436
	1	Violet	18"/45cm	12	SXPP1A400
	1	Violet	24"/60cm	12	SXPP1A445
 <p>UR-6 26mm 5/8 Circle</p>	0	Violet	6"/15cm	12	SXPP1A447
	2-0	Violet	6"/15cm	12	SXPP1A417


STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device






TAPERCUT® 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 V-7 26mm 1/2 Circle	2-0	Violet	9"/23cm	12	SXPP1A424
 V-34 36mm 1/2 Circle	0	Violet	9"/23cm	12	SXPP1A426


Reverse Cutting 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 CP-2 26mm 1/2 Circle	3-0	Violet	18"/45cm	12	SXPP1A437
 OS-6 36mm 1/2 Circle	0	Violet	18"/45cm	12	SXPP1A200
	0	Violet	24"/60cm	12	SXPP1A202
	1	Violet	18"/45cm	12	SXPP1A201
	1	Violet	24"/60cm	12	SXPP1A203
 OS-8 40mm 1/2 Circle	0	Violet	24"/60cm	12	SXPP1A204
	1	Violet	24"/60cm	12	SXPP1A205

STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device



PRIME Reverse Cutting 					
Needle	Suture USP	Color	Length	EA/BX	Product Code
 PS-2 19mm 3/8 Circle	3-0	Violet	18"/45cm	12	SXPP1A101
 PS-1 24mm 3/8 Circle	3-0	Violet	18"/45cm	12	SXPP1A100

Plus Antibacterial Sutures




Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture





VICRYL Plus Suture is indicated for use in general soft tissue approximation and/or ligation, except for ophthalmic, cardiovascular and neurological tissues.¹

- Shown in vitro to inhibit bacterial colonization of the suture for 7 days or more²




Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture














ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB-2 26mm 1/2 Circle	2-0	Undyed	1 x 27"/70cm	36		VCPB269H
	2-0	Violet	8 x 18"/45cm	12	CR	VCPB726D*
	0	Violet	8 x 18"/45cm	12	CR	VCPB727D*
 SHB 26mm 1/2 Circle	3-0	Undyed	1 x 27"/70cm	36		VCPB416H
	3-0	Undyed	8 x 18"/45cm	12	CR	VCPB864D*
	3-0	Violet	8 x 18"/45cm	12	CR	VCPB774D*
	2-0	Undyed	1 x 27"/70cm	36		VCPB417H
 CTB-1 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		VCPB944H
	2-0	Undyed	8 x 18"/45cm	12	CR	VCPB839D*
	2-0	Undyed	1 x 27"/70cm	36		VCPB259H
	2-0	Undyed	1 x 36"/90cm	36		VCPB945H
	2-0	Violet	8 x 18"/45cm	12	CR	VCPB739D*
	0	Undyed	8 x 18"/45cm	12	CR	VCPB840D*
	0	Violet	8 x 18"/45cm	12	CR	VCPB740D*
	0	Undyed	1 x 27"/70cm	36		VCPB260H
	0	Undyed	8 x 27"/70cm	12	CR	VCPB41D*
	0	Violet	1 x 27"/70cm	36		VCPB340H
	0	Violet	8 x 27"/70cm	12	CR	VCPPB31D*
	0	Undyed	1 x 36"/90cm	36		VCPB946H
	1	Undyed	8 x 18"/45cm	12	CR	VCPB841D*
	1	Violet	8 x 18"/45cm	12	CR	VCPB741D*
	1	Violet	1 x 27"/70cm	36		VCPB341H
1	Undyed	1 x 36"/90cm	36		VCPB947H	

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture

Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture








ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB 40mm 1/2 Circle	0	Undyed	1 x 36"/90cm	36		VCPB958H
 CTXB 48mm 1/2 Circle	0	Undyed	8 x 18"/45cm	12		VCPB724D*
	0	Undyed	1 x 36"/90cm	36		VCPB978H
	1	Undyed	8 x 18"/45cm	12		VCPB725D*
	1	Violet	8 x 18"/45cm	12		VCPB765D*
	1	Undyed	1 x 36"/90cm	36		VCPB977H
 TPB-1 65mm 1/2 Circle	1	Violet	4 x 27"/70cm	12		VCPB650G

Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12		VCP743D*
 RB-2 13mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		VCP433H
	4-0	Undyed	1 x 27"/70cm	36		VCP434H


*Safety organizer tray  CONTROL RELEASE™ Needle Suture




Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		VCP213H
	5-0	Violet	1 x 27"/70cm	36		VCP303H
	4-0	Undyed	8 x 18"/45cm	12	CR	VCP714D*
	4-0	Undyed	1 x 27"/70cm	36		VCP214H
	4-0	Violet	1 x 27"/70cm	36		VCP304H
	3-0	Violet	8 x 18"/45cm	12	CR	VCP713D*
	3-0	Undyed	1 x 27"/70cm	36		VCP215H
	3-0	Violet	1 x 27"/70cm	36		VCP305H
	2-0	Violet	1 x 27"/70cm	36		VCP306H
 SH-1 22mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12	CR	VCP771D*
	4-0	Undyed	1 x 27"/70cm	36		VCP218H
	4-0	Violet	1 x 27"/70cm	36		VCP310H
	3-0	Violet	8 x 18"/45cm	12	CR	VCP772D*
	3-0	Undyed	1 x 27"/70cm	36		VCP219H
	3-0	Violet	1 x 27"/70cm	36		VCP311H
 CT-3 22mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		VCP327H
	2-0	Violet	1 x 27"/70cm	36		VCP328H
	0	Violet	1 x 27"/70cm	36		VCP329H















Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-2 26mm 1/2 Circle</p>	3-0	Undyed	1 x 27"/70cm	36		VCP232H
	3-0	Violet	1 x 27"/70cm	36		VCP332H
	2-0	Violet	8 x 18"/45cm	12	CR	VCP726D*
	2-0	Undyed	1 x 27"/70cm	36		VCP269H
	2-0	Violet	1 x 27"/70cm	36		VCP333H
	0	Violet	8 x 18"/45cm	12	CR	VCP727D*
	0	Undyed	1 x 27"/70cm	36		VCP270H
	0	Violet	1 x 27"/70cm	36		VCP334H
	1	Violet	1 x 27"/70cm	36		VCP335H
 <p>SH 26mm 1/2 Circle</p>	4-0	Violet	8 x 18"/45cm	12	CR	VCP773D*
	4-0	Undyed	1 x 27"/70cm	36		VCP415H
	4-0	Violet	1 x 27"/70cm	36		VCP315H
	3-0	Undyed	8 x 18"/45cm	12	CR	VCP864D*
	3-0	Violet	8 x 18"/45cm	12	CR	VCP774D*
	3-0	Undyed	1 x 27"/70cm	36		VCP416H
	3-0	Violet	1 x 27"/70cm	36		VCP316H
	3-0	Violet	8 x 27"/70cm	12	CR	VCP784D*
	2-0	Violet	8 x 18"/45cm	12	CR	VCP775D*
	2-0	Undyed	1 x 27"/70cm	36		VCP417H
	2-0	Violet	1 x 27"/70cm	36		VCP317H
	2-0	Violet	8 x 27"/70cm	12	CR	VCP785D*
	0	Undyed	1 x 27"/70cm	36		VCP418H


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	3-0	Violet	1 x 36"/90cm	36		VCP527H
 MO-6 26mm 1/2 Circle	2-0	Violet	8 x 18"/45cm	12		VCP789D*
 MO-5 31mm 1/2 Circle	0	Violet	8 x 18"/45cm	12		VCP801D*
	0	Violet	8 x 27"/70cm	12		VCPP81D*
	1	Violet	8 x 27"/70cm	12		VCPP80D*
 CT-1 36mm 1/2 Circle	3-0	Undyed	8 x 18"/45cm	12		VCP838D*
	3-0	Violet	8 x 18"/45cm	12		VCP738D*
	3-0	Undyed	1 x 27"/70cm	36		VCP258H
	3-0	Violet	1 x 27"/70cm	36		VCP338H
	3-0	Undyed	1 x 36"/90cm	36		VCP944H
	3-0	Violet	1 x 36"/90cm	36		VCP344H
	2-0	Undyed	8 x 18"/45cm	12		VCP839D*
	2-0	Violet	8 x 18"/45cm	12		VCP739D*
	2-0	Undyed	1 x 27"/70cm	36		VCP259H
	2-0	Undyed	8 x 27"/70cm	12		VCPP42D*
	2-0	Violet	1 x 27"/70cm	36		VCP339H
	2-0	Undyed	1 x 36"/90cm	36		VCP945H
2-0	Violet	1 x 36"/90cm	36		VCP345H	

*Safety organizer tray  CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	0	Undyed	8 x 18"/45cm	12	CR	VCP840D*
	0	Violet	8 x 18"/45cm	12	CR	VCP740D*
	0	Undyed	1 x 27"/70cm	36		VCP260H
	0	Undyed	8 x 27"/70cm	12	CR	VCPP41D*
	0	Violet	1 x 27"/70cm	36		VCP340H
	0	Violet	8 x 27"/70cm	12	CR	VCPP31D*
	0	Undyed	1 x 36"/90cm	36		VCP946H
	0	Violet	1 x 36"/90cm	36		VCP346H
	1	Undyed	8 x 18"/45cm	12	CR	VCP841D*
	1	Violet	8 x 18"/45cm	12	CR	VCP741D*
	1	Undyed	1 x 27"/70cm	36		VCP261H
	1	Undyed	8 x 27"/70cm	12	CR	VCPP40D*
	1	Violet	1 x 27"/70cm	36		VCP341H
	1	Undyed	1 x 36"/90cm	36		VCP947H
1	Violet	1 x 36"/90cm	36		VCP347H	
 <p>MH 36mm 1/2 Circle</p>	3-0	Violet	1 x 27"/70cm	36		VCP322H
	2-0	Violet	1 x 27"/70cm	36		VCP323H
 <p>MO-4 36mm 1/2 Circle</p>	2-0	Violet	8 x 18"/45cm	12	CR	VCP700D*
	0	Violet	8 x 18"/45cm	12	CR	VCP701D*
	0	Violet	1 x 27"/70cm	36		VCP436H
	0	Violet	8 x 27"/70cm	12	CR	VCPP71D*
	1	Violet	8 x 18"/45cm	12	CR	VCP702D*
	1	Violet	1 x 27"/70cm	36		VCP437H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm 1/2 Circle</p>	3-0	Undyed	1 x 36"/90cm	36		VCP956H
	3-0	Violet	1 x 36"/90cm	36		VCP356H
	2-0	Violet	8 x 18"/45cm	12	CR	VCP751D*
	2-0	Undyed	1 x 27"/70cm	36		VCP275H
	2-0	Violet	1 x 27"/70cm	36		VCP351H
	2-0	Undyed	1 x 36"/90cm	36		VCP957H
	2-0	Violet	1 x 36"/90cm	36		VCP357H
	0	Violet	8 x 18"/45cm	12	CR	VCP752D*
	0	Undyed	1 x 27"/70cm	36		VCP280H
	0	Violet	1 x 27"/70cm	36		VCP352H
	0	Undyed	1 x 36"/90cm	36		VCP958H
	0	Violet	1 x 36"/90cm	36		VCP358H
	1	Violet	8 x 18"/45cm	12	CR	VCP753D*
	1	Undyed	1 x 27"/70cm	36		VCP281H
	1	Violet	1 x 27"/70cm	36		VCP353H
	1	Undyed	1 x 36"/90cm	36		VCP959H
1	Violet	1 x 36"/90cm	36		VCP359H	
 <p>CTX 48mm 1/2 Circle</p>	3-0	Undyed	1 x 36"/90cm	36		VCP980H
	2-0	Undyed	8 x 18"/45cm	12	CR	VCP723D*
	2-0	Violet	1 x 27"/70cm	36		VCP363H
	2-0	Undyed	1 x 36"/90cm	36		VCP979H
	2-0	Violet	1 x 36"/90cm	36		VCP369H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	0	Undyed	8 x 18"/45cm	12	CR	VCP724D*
	0	Violet	8 x 18"/45cm	12	CR	VCP764D*
	0	Violet	1 x 27"/70cm	36		VCP364H
	0	Undyed	1 x 36"/90cm	36		VCP978H
	0	Violet	1 x 36"/90cm	36		VCP370H
	1	Violet	8 x 18"/45cm	12	CR	VCP765D*
	1	Violet	1 x 27"/70cm	36		VCP365H
	1	Undyed	1 x 36"/90cm	36		VCP977H
	1	Violet	1 x 36"/90cm	36		VCP371H
 <p>TP-1 65mm 1/2 Circle</p>	2	Undyed	2 x 27"/70cm	12		VCP849G
	2	Violet	4 x 27"/70cm	12		VCP649G
	2	Undyed	1 x 54"/135cm	24		VCP880T
 <p>XLH 70mm 1/2 Circle</p>	2-0	Violet	1 x 27"/70cm	12		VCP581G
	1	Violet	1 x 27"/70cm	12		VCP583G
 <p>UR-6 26mm 5/8 Circle</p>	2-0	Violet	1 x 27"/70cm	36		VCP602H
	0	Violet	1 x 27"/70cm	36		VCP603H
 <p>UR-5 36mm 5/8 Circle</p>	2-0	Violet	1 x 27"/70cm	36		VCP375H
	0	Violet	1 x 27"/70cm	36		VCP376H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture









TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-34 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		VCP916H
	3-0	Violet	1 x 36"/90cm	36		VCP516H
	2-0	Undyed	1 x 36"/90cm	36		VCP917H
	2-0	Violet	1 x 36"/90cm	36		VCP517H
	0	Violet	1 x 36"/90cm	36		VCP518H
	1	Violet	1 x 36"/90cm	36		VCP519H

Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Violet	1 x 18"/45cm	36		VCP391H
	5-0	Undyed	1 x 27"/70cm	36		VCP421H
	4-0	Violet	1 x 18"/45cm	36		VCP392H
	4-0	Undyed	1 x 27"/70cm	36		VCP422H
	3-0	Violet	1 x 18"/45cm	36		VCP393H
	3-0	Undyed	1 x 27"/70cm	36		VCP423H


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture







Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-1 24mm 3/8 Circle	4-0	Undyed	1 x 27"/70cm	36		VCP441H
	4-0	Violet	1 x 27"/70cm	36		VCP451H
	3-0	Undyed	1 x 27"/70cm	36		VCP442H
	3-0	Violet	1 x 27"/70cm	36		VCP452H
	2-0	Undyed	1 x 27"/70cm	36		VCP443H
 FS 26mm 3/8 Circle	3-0	Violet	8 x 18"/45cm	12	CR	VCP787D*
	3-0	Violet	8 x 18"/45cm	12	CR	VCP787Z
 FSL 30mm 3/8 Circle	3-0	Undyed	1 x 27"/70cm	36		VCP588H
	2-0	Undyed	1 x 27"/70cm	36		VCP589H
 PSL 30mm 3/8 Circle	3-0	Undyed	1 x 27"/70cm	36		VCP502H
	2-0	Undyed	1 x 27"/70cm	36		VCP596H
 OS-4 22mm 1/2 Circle	2-0	Undyed	3 x 18"/45cm	24	CR	VCP748T
	2-0	Violet	3 x 18"/45cm	24	CR	VCP706T
	0	Undyed	3 x 18"/45cm	24	CR	VCP749T
	0	Undyed	1 x 27"/70cm	36		VCP694H
	1	Violet	3 x 18"/45cm	24	CR	VCP708T
	1	Undyed	1 x 27"/70cm	36		VCP695H
 X-1 22mm 1/2 Circle	3-0	Violet	8 x 18"/45cm	12	CR	VCP790D*
	3-0	Undyed	1 x 27"/70cm	36		VCP458H
	3-0	Violet	1 x 27"/70cm	36		VCP460H
	2-0	Undyed	1 x 27"/70cm	36		VCP459H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture







Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CP-2 26mm 1/2 Circle</p>	3-0	Undyed	8 x 18"/45cm	12	CR	VCP761D*
	3-0	Undyed	1 x 27"/70cm	36		VCP868H
	2-0	Undyed	8 x 18"/45cm	12	CR	VCP762D*
	2-0	Undyed	1 x 27"/70cm	36		VCP869H
	0	Undyed	1 x 27"/70cm	36		VCP870H
	1	Undyed	1 x 27"/70cm	36		VCP871H
 <p>CP-1 36mm 1/2 Circle</p>	3-0	Undyed	1 x 27"/70cm	36		VCP265H
	2-0	Undyed	1 x 27"/70cm	36		VCP266H
	2-0	Violet	1 x 27"/70cm	36		VCP466H
	2-0	Violet	1 x 36"/90cm	36		VCP472H
	0	Undyed	1 x 27"/70cm	36		VCP267H
	0	Violet	1 x 27"/70cm	36		VCP467H
	1	Undyed	1 x 27"/70cm	36		VCP268H
	1	Violet	1 x 27"/70cm	36		VCP468H
 <p>OS-6 36mm 1/2 Circle</p>	2-0	Undyed	3 x 18"/45cm	24	CR	VCP750T
	2-0	Violet	3 x 18"/45cm	24	CR	VCP710T
	2-0	Undyed	1 x 27"/70cm	36		VCP533H
	0	Undyed	3 x 18"/45cm	24	CR	VCP754T
	0	Violet	3 x 18"/45cm	24	CR	VCP711T
	0	Undyed	1 x 27"/70cm	36		VCP534H
	1	Undyed	1 x 27"/70cm	36		VCP535H
	 <p>CP 40mm 1/2 Circle</p>	2-0	Violet	1 x 27"/70cm	36	
0		Undyed	1 x 27"/70cm	36		VCP197H
1		Undyed	1 x 27"/70cm	36		VCP196H
1		Violet	1 x 27"/70cm	36		VCP480H
1		Violet	1 x 36"/90cm	36		VCP486H
2		Undyed	1 x 27"/70cm	36		VCP195H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture


























Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>OS-8 40mm 1/2 Circle</p>	2-0	Violet	3 x 18"/45cm	24	CR	VCP716T
	0	Undyed	3 x 18"/45cm	24	CR	VCP756T
	0	Violet	3 x 18"/45cm	24	CR	VCP717T
	0	Undyed	1 x 27"/70cm	36		VCP698H
	1	Undyed	3 x 18"/45cm	24	CR	VCP757T
	1	Violet	3 x 18"/45cm	24	CR	VCP718T
	1	Undyed	1 x 27"/70cm	36		VCP699H
	2	Violet	3 x 18"/45cm	24	CR	VCP719T
 <p>CPX 48mm 1/2 Circle</p>	0	Undyed	1 x 27"/70cm	36		VCP568H
	1	Undyed	1 x 27"/70cm	36		VCP569H

Taperpoint ⊙, Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm/ FSL 30mm 1/2 Circle/ 3/8 Circle</p>	3-0	Undyed	1 x 54"/135cm	36		VCP940H
 <p>CT-1, CP-1 36mm 1/2 Circle</p>	3-0	Undyed	1 x 54"/135cm	36		VCP931H


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture






















Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 x 27"/70cm	36		VCP662H
	3-0	Undyed	1 x 27"/70cm	36		VCP663H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP490G
 P-3 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP493G
	5-0	Undyed	1 x 18"/45cm	36		VCP493H
	5-0	Violet	1 x 18"/45cm	12		VCP463G
	4-0	Undyed	1 x 18"/45cm	12		VCP494G
	4-0	Undyed	1 x 18"/45cm	36		VCP494H
	4-0	Violet	1 x 18"/45cm	12		VCP464G
 PS-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP500G
	5-0	Undyed	1 x 18"/45cm	36		VCP500H
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP495G
	5-0	Undyed	1 x 18"/45cm	36		VCP495H
	4-0	Undyed	1 x 18"/45cm	12		VCP496G
	4-0	Undyed	1 x 18"/45cm	36		VCP496H
	4-0	Undyed	1 x 27"/70cm	36		VCP426H
	3-0	Undyed	1 x 18"/45cm	12		VCP497G
	3-0	Undyed	1 x 18"/45cm	36		VCP497H
	3-0	Undyed	1 x 27"/70cm	36		VCP427H
	2-0	Undyed	1 x 27"/70cm	36		VCP428H


Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture












PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12		VCP682G
	4-0	Undyed	1 x 18"/45cm	36		VCP682H
	4-0	Undyed	1 x 27"/70cm	36		VCP935H
	3-0	Undyed	1 x 18"/45cm	12		VCP683G
	3-0	Undyed	1 x 18"/45cm	36		VCP683H
	3-0	Undyed	1 x 27"/70cm	36		VCP936H
 P-2 8mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP503G
	4-0	Undyed	1 x 18"/45cm	12		VCP504G
 PS-6 10mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP511G
 PS-5 13mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP593G
	4-0	Undyed	1 x 18"/45cm	12		VCP594G
 PS-4 16mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP506G
	4-0	Undyed	1 x 18"/45cm	12		VCP507G
 PS-4C 16mm Compound Curve	4-0	Undyed	1 x 18"/45cm	12		VCP656G

PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP834G
	4-0	Undyed	1 x 18"/45cm	12		VCP835G

Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture




PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP844G
	4-0	Undyed	1 x 18"/45cm	12		VCP845G
 PC-5 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VCP822G
	4-0	Undyed	1 x 18"/45cm	12		VCP823G
	4-0	Undyed	1 x 18"/45cm	36		VCP823H
	3-0	Undyed	1 x 18"/45cm	12		VCP824G

LIGAPAK™ Dispensing Reel						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Undyed	1 x 54"/135cm	12		VCP284G
	3-0	Undyed	1 x 54"/135cm	12		VCP285G
	3-0	Violet	1 x 54"/135cm	12		VCP205G
	3-0	Undyed	1 x 144"/366cm	12		VCP885G
	2-0	Undyed	1 x 54"/135cm	12		VCP286G
	0	Undyed	1 x 54"/135cm	12		VCP287G
	0	Violet	1 x 54"/135cm	12		VCP207G
	1	Violet	1 x 54"/135cm	12		VCP208G


SUTUPAK® Pre-Cut Sterile Sutures						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Undyed	12 x 18"/45cm	12		VCP109G
	4-0	Violet	12 x 18"/45cm	12		VCP103G
	3-0	Violet	3 x 18"/45cm	24		VCP634T

Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture



SUTUPAK® Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Undyed	12 x 18"/45cm	12		VCP110G
	3-0	Violet	12 x 18"/45cm	12		VCP104G
	2-0	Undyed	3 x 18"/45cm	24		VCP645T
	2-0	Violet	3 x 18"/45cm	24		VCP635T
	2-0	Undyed	12 x 18"/45cm	12		VCP111G
	2-0	Violet	12 x 18"/45cm	12		VCP105G
	2-0	Undyed	1 x 54"/135cm	36		VCP607H
	2-0	Violet	1 x 54"/135cm	36		VCP615H
	0	Undyed	3 x 18"/45cm	24		VCP646T
	0	Violet	3 x 18"/45cm	24		VCP636T
	0	Undyed	6 x 18"/45cm	12		VCP112G
	0	Violet	6 x 18"/45cm	12		VCP106G
	0	Undyed	1 x 54"/135cm	36		VCP608H
	0	Violet	1 x 54"/135cm	36		VCP616H
	1	Violet	6 x 18"/45cm	12		VCP107G
	1	Violet	1 x 54"/135cm	36		VCP617H
	2	Violet	1 x 54"/135cm	36		VCP618H




MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture


MONOCRYL Plus Sutures are ideal for subcuticular skin closure.




- High initial breaking strength to maintain wound approximation during the critical first few days of skin healing¹
- Shown in vitro to inhibit bacterial colonization of the suture for 7 days or more²
- Monofilament design and polymer properties minimize drag force and elicit only a slight tissue reaction during absorption¹




MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture









ETHIGUARD® Blunt Point Needle						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SHB 26mm 1/2 Circle	3-0	Undyed	1 x 27"/70cm	36		MCPB416H
	3-0	Violet	1 x 27"/70cm	36		MCPB316H
 CTB-1 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		MCPB944H
	2-0	Undyed	1 x 36"/90cm	36		MCPB945H
	0	Undyed	1 x 36"/90cm	36		MCPB946H
	0	Violet	1 x 36"/90cm	36		MCPB346H


Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		MCP213H
	4-0	Undyed	1 x 27"/70cm	36		MCP214H
	4-0	Violet	1 x 27"/70cm	36		MCP304H
 CT-2 26mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		MCP333H
 SH 26mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	36		MCP415H
	4-0	Violet	1 x 27"/70cm	36		MCP315H
	3-0	Undyed	1 x 27"/70cm	36		MCP416H
	3-0	Violet	1 x 27"/70cm	36		MCP316H
	2-0	Undyed	1 x 27"/70cm	36		MCP417H
	2-0	Violet	1 x 27"/70cm	36		MCP317H


MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture























Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		MCP944H
	3-0	Violet	1 x 36"/90cm	36		MCP344H
	2-0	Undyed	1 x 36"/90cm	36		MCP945H
	2-0	Violet	1 x 36"/90cm	36		MCP345H
	2-0	Violet	1 x 27"/70cm	36		MCP339H
	0	Undyed	1 x 36"/90cm	36		MCP946H
	0	Violet	1 x 36"/90cm	36		MCP346H
	0	Violet	1 x 27"/70cm	36		MCP340H
	1	Violet	1 x 36"/90cm	36		MCP347H
 CT 40mm 1/2 Circle	0	Violet	1 x 36"/90cm	36		MCP358H
 UR-6 26mm 5/8 Circle	3-0	Violet	1 x 27"/70cm	36		MCP604H
	2-0	Violet	1 x 27"/70cm	36		MCP605H
 UR-5 36mm 5/8 Circle	2-0	Violet	1 x 27"/70cm	36		MCP377H
 CP-1 36mm 1/2 Circle	2-0	Undyed	1 x 27"/70cm	36		MCP266H

MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture











Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Undyed	1 x 27"/70cm	36		MCP266H


Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Undyed	1 x 27"/70cm	36		MCP523H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		MCP489G
	5-0	Undyed	1 x 18"/45cm	12		MCP490G
 P-3 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		MCP492G
	5-0	Undyed	1 x 18"/45cm	12		MCP493G
	5-0	Violet	1 x 18"/45cm	12		MCP463G
	4-0	Undyed	1 x 18"/45cm	12		MCP494G
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		MCP495G
	4-0	Undyed	1 x 27"/70cm	36		MCP426H
	4-0	Undyed	1 x 18"/45cm	12		MCP496G
	4-0	Violet	1 x 18"/45cm	12		MCP513G
	3-0	Undyed	1 x 27"/70cm	36		MCP427H
	3-0	Undyed	1 x 18"/45cm	12		MCP497G
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 27"/70cm	36		MCP935H
	4-0	Undyed	1 x 18"/45cm	36		MCP682H
	3-0	Undyed	1 x 27"/70cm	36		MCP936H

MONOCRYL® Plus Antibacterial (poliglecaprone 25) Suture



PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		MCP834G
	4-0	Undyed	1 x 18"/45cm	12		MCP835G
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		MCP844G
	4-0	Undyed	1 x 18"/45cm	12		MCP845G
 PC-5 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		MCP822G
	4-0	Undyed	1 x 18"/45cm	12		MCP823G




PDS[®] Plus Antibacterial (polydioxanone) Suture




PDS Plus Sutures combine the features of soft, pliable, monofilament construction with absorbability and extended wound support¹:





- Retains 60% of its original strength for 6 weeks, providing support to the fascia as it slowly heals²
- Shown in vitro to inhibit bacterial colonization of the suture³
- Monofilament design and polymer properties minimize drag force and elicit only a slight tissue reaction during absorption¹




PDS® Plus Antibacterial (polydioxanone) Suture







Blunt Point 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BP-2 48mm 1/2 Circle	0	Violet	1 x 60"/150cm	12		PDPB990G
 BP-1 65mm 1/2 Circle	0	Violet	30"/75cm Loop Length	12	Looped	PDPB991G
	1	Violet	48"/122cm Loop Length	12	Looped	PDPB880G


ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SHB 26mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		PDPB316H
	2-0	Violet	1 x 27"/70cm	36		PDPB317H
 CTB-2 26mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		PDPB332H
	2-0	Violet	1 x 27"/70cm	36		PDPB333H
	0	Violet	1 x 27"/70cm	36		PDPB334H
 CTB-1 36mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		PDPB339H
	2-0	Violet	1 x 36"/90cm	36		PDPB345
	0	Violet	1 x 27"/70cm	36		PDPB340
	0	Violet	1 x 36"/90cm	36		PDPB346
	1	Violet	1 x 36"/90cm	36		PDPB347

PDS® Plus Antibacterial (polydioxanone) Suture




ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB 40mm 1/2 Circle	0	Violet	1 x 36"/90cm	24		PDPB358T
	1	Violet	1 x 36"/90cm	24		PDPB359T
 CTXB 48mm 1/2 Circle	0	Violet	1 x 36"/90cm	24		PDPB370
	1	Violet	1 x 36"/90cm	24		PDPB371
 TPB-1 65mm 1/2 Circle	1	Violet	1 x 54"/135cm	12		PDPB879G








PDS[®] Plus Antibacterial (polydioxanone) Suture



Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	6-0	Violet	1 x 30"/75cm	36		PDP117H
 BV 11mm 3/8 Circle	6-0	Violet	1 x 30"/75cm	36		PDP123H
 C-1 13mm 3/8 Circle	6-0	Violet	1 x 30"/75cm	36		PDP127H
	5-0	Violet	1 x 30"/75cm	36		PDP126H


PDS® Plus Antibacterial (polydioxanone) Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-2 13mm 1/2 Circle	6-0	Violet	1 x 30"/75cm	36		PDP149H
	5-0	Violet	1 x 30"/75cm	36		PDP148H
 TF 13mm 1/2 Circle	6-0	Violet	1 x 27"/70cm	36		PDP432H
	4-0	Violet	8 x 18"/45cm	12	CR	PDP802D*
 UCL 13mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		PDP114H
 RB-1 17mm 1/2 Circle	5-0	Violet	8 x 24"/60cm	12		PDP800G
	5-0	Violet	1 x 27"/70cm	36		PDP303H
	4-0	Violet	8 x 24"/60cm	12	CR	PDP801D*
	4-0	Violet	1 x 27"/70cm	36		PDP304H
	4-0	Violet	8 x 27"/70cm	12	CR	PDP789D*
	3-0	Violet	1 x 27"/70cm	36		PDP305H
 RB-1 17mm 1/2 Circle	5-0	Violet	1 x 30"/75cm	36		PDP320H
 SH-1 22mm 1/2 Circle	4-0	Violet	1 x 27"/70cm	36		PDP310H
	3-0	Violet	8 x 24"/60cm	12	CR	PDP779D*
	3-0	Violet	1 x 27"/70cm	36		PDP311H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


PDS® Plus Antibacterial (polydioxanone) Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-2 26mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		PDP332H
	2-0	Violet	1 x 27"/70cm	36		PDP333H
	0	Violet	1 x 27"/70cm	36		PDP334H
 SH 26mm 1/2 Circle	4-0	Violet	1 x 27"/70cm	36		PDP315H
	3-0	Violet	1 x 27"/70cm	36		PDP316H
	3-0	Undyed	1 x 27"/70cm	36		PDP416H
	2-0	Violet	1 x 27"/70cm	36		PDP317H
 CT-1 36mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36	CR	PDP338H
	2-0	Violet	8 x 18"/45cm	12		PDP739D*
	2-0	Violet	1 x 27"/70cm	36		PDP339H
	2-0	Undyed	1 x 27"/70cm	36		PDP259H
	0	Violet	1 x 27"/70cm	36		PDP340H
	0	Violet	1 x 36"/90cm	36		PDP346H
	1	Violet	1 x 27"/70cm	36		PDP341H
	1	Violet	1 x 36"/90cm	36		PDP347H
 CT 40mm 1/2 Circle	2-0	Violet	1 x 36"/90cm	36		PDP357H
	0	Violet	1 x 27"/70cm	36		PDP352H
	0	Violet	1 x 36"/90cm	24		PDP358T
	1	Violet	1 x 27"/70cm	36		PDP353H
	1	Violet	1 x 36"/90cm	24		PDP359T

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


PDS® Plus Antibacterial (polydioxanone) Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	0	Violet	8 x 18"/45cm	12	CR	PDP776D*
	0	Violet	1 x 36"/90cm	24		PDP370T
	0	Violet	30"/75cm Loop Length	12	Looped	PDP990G
	1	Violet	1 x 36"/90cm	24		PDP371T
 <p>TP-1 65mm 1/2 Circle</p>	0	Violet	30"/75cm Loop Length	12	Looped	PDP991G
	1	Violet	1 x 54"/135cm	12		PDP879G
	1	Violet	48"/122cm Loop Length	12	Looped	PDP880G
 <p>XLH 70mm 1/2 Circle</p>	0	Violet	1 x 27"/70cm	12		PDP582G
	1	Violet	1 x 54"/135cm	12		PDP877G
	1	Violet	48"/122cm Loop Length	12	Looped	PDP881G
 <p>ST-4 19mm Straight</p>	4-0	Violet	1 x 14"/35cm	12		PDP420G

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture

PDS® Plus Antibacterial (polydioxanone) Suture







TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-5 17mm 1/2 Circle	4-0	Violet	1 x 36"/90cm	36		PDP137H


Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Undyed	1 x 27"/70cm	36		PDP421H
	4-0	Undyed	1 x 27"/70cm	36		PDP422H
	3-0	Undyed	1 x 27"/70cm	36		PDP423H
 FS-1 24mm 3/8 Circle	4-0	Undyed	1 x 27"/70cm	36		PDP441H
	3-0	Undyed	1 x 27"/70cm	36		PDP442H
	2-0	Undyed	1 x 27"/70cm	36		PDP443H
 OS-4 22mm 1/2 Circle	1	Violet	1 x 27"/70cm	36		PDP695H
 CP-2 26mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		PDP968H
	2-0	Violet	8 x 18"/45cm	12	CR	PDP762D*
	2-0	Violet	1 x 27"/70cm	36		PDP969H
	0	Violet	1 x 27"/70cm	36		PDP970H








*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture

PDS® Plus Antibacterial (polydioxanone) Suture
















Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		PDP466H
	0	Violet	1 x 27"/70cm	36		PDP467H
	1	Violet	1 x 27"/70cm	36		PDP468H
 OS-6 36mm 1/2 Circle	0	Violet	1 x 27"/70cm	24		PDP534T
	1	Violet	1 x 27"/70cm	36		PDP535H
 CPX 48mm 1/2 Circle	1	Violet	1 x 27"/70cm	24		PDP569T










Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 x 27"/70cm	36		PDP662H
	3-0	Undyed	1 x 27"/70cm	36		PDP663H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Violet	1 x 18"/45cm	12		PDP487G
	6-0	Undyed	1 x 18"/45cm	12		PDP489G
 P-3 13mm 3/8 Circle	5-0	Violet	1 x 18"/45cm	12		PDP463G
	5-0	Undyed	1 x 18"/45cm	12		PDP493G
	4-0	Undyed	1 x 18"/45cm	12		PDP494G


PDS® Plus Antibacterial (polydioxanone) Suture



PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		PDP495G
	4-0	Undyed	1 x 18"/45cm	12		PDP496G
	4-0	Violet	1 x 18"/45cm	12		PDP513G
	3-0	Undyed	1 x 18"/45cm	12		PDP497G
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12		PDP682G
	3-0	Undyed	1 x 18"/45cm	12		PDP683G
 P-2 8mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		PDP503G
	4-0	Undyed	1 x 18"/45cm	12		PDP504G
 PS-4 16mm 1/2 Circle	4-0	Undyed	1 x 18"/45cm	12		PDP507G

PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		PDP833G
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		PDP844G
	4-0	Undyed	1 x 18"/45cm	12		PDP845G
 PC-5 19mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12		PDP823G
	3-0	Undyed	1 x 18"/45cm	12		PDP824G

Sutures: **Absorbable**








VICRYL RAPIDE™ (polyglactin 910) Suture


VICRYL RAPIDE Suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.¹



VICRYL RAPIDE™ (polyglactin 910) Suture


















Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	12		VR214
 SH 26mm 1/2 Circle	3-0	Undyed	1 x 27"/70cm	12		VR416
	2-0	Undyed	1 x 27"/70cm	12		VR417
 CT-1 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	12		VR944
	2-0	Undyed	1 x 36"/90cm	12		VR945
 CTX 48mm 1/2 Circle	0	Undyed	1 x 36"/90cm	12		VR1695

TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-34 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	12		VR916
	2-0	Undyed	1 x 36"/90cm	12		VR917
	0	Undyed	1 x 36"/90cm	36		VR2218


VICRYL RAPIDE™ (polyglactin 910) Suture











Taperpoint  , Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm/ CP-1 36mm 1/2 Circle/ 1/2 Circle	3-0	Undyed	1 x 54"/135cm	12		VR931
	2-0	Undyed	1 x 54"/135cm	12		VR932

PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VR490
 P-3 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		VR493
	4-0	Undyed	1 x 18"/45cm	12		VR494
 PS-2 19mm 3/8 Circle	4-0	Undyed	1 x 27"/70cm	12		VR426
	4-0	Undyed	1 x 18"/45cm	12		VR496
	3-0	Undyed	1 x 18"/45cm	12		VR497
 PS-1 24mm 3/8 Circle	3-0	Undyed	1 x 27"/70cm	12		VR935

VICRYL RAPIDE™ (polyglactin 910) Suture



PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	VR834
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	VR835
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	VR844
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	VR845




Coated VICRYL[®] (polyglactin 910) Suture





VICRYL Suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.¹

VICRYL Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹















Coated VICRYL® (polyglactin 910) Suture




Blunt Point ETHIGUARD® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SHB 26mm 1/2 Circle	3-0	Undyed	8 x 18"/45cm	12	CR	JB864
 CTB-1 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		JB944
	2-0	Undyed	8 x 18"/45cm	12	CR	JB839
	2-0	Undyed	1 x 27"/70cm	36		JB259
	2-0	Undyed	1 x 36"/90cm	36		JB945
	0	Undyed	8 x 18"/45cm	12	CR	JB840
	0	Undyed	1 x 27"/70cm	36		JB260
	0	Violet	1 x 27"/70cm	36		JB340
	0	Undyed	1 x 36"/90cm	36		JB946
	1	Undyed	8 x 18"/45cm	12	CR	JB841
	1	Violet	1 x 27"/70cm	36		JB341
1	Undyed	1 x 36"/90cm	36		JB947	
 CTXB 48mm 1/2 Circle	0	Undyed	8 x 18"/45cm	12	CR	JB724
	0	Undyed	1 x 36"/90cm	36		JB978
	1	Undyed	8 x 18"/45cm	12	CR	JB725
	1	Undyed	1 x 36"/90cm	36		JB977










Coated VICRYL® (polyglactin 910) Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-4 5mm 3/8 Circle	8-0	Violet	1 x 5"/13cm	12		J405G
 BV130-5 6.5mm 3/8 Circle	8-0	Violet	1 x 5"/13cm	12		J401G
 TF 13mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		J433H
	4-0	Violet	8 x 18"/45cm	12		J743D*
	4-0	Undyed	1 x 27"/70cm	36		J434H
 RB-1 17mm 1/2 Circle	6-0	Undyed	1 x 27"/70cm	36		J212H
	5-0	Undyed	1 x 27"/70cm	36		J213H
	5-0	Violet	1 x 27"/70cm	36		J303H
	4-0	Undyed	8 x 18"/45cm	12		J714D*
	4-0	Undyed	1 x 27"/70cm	36		J214H
	4-0	Violet	1 x 27"/70cm	36		J304H
	3-0	Violet	8 x 18"/45cm	12		J713D*
	3-0	Undyed	1 x 27"/70cm	36		J215H
	3-0	Violet	1 x 27"/70cm	36		J305H
	2-0	Violet	1 x 27"/70cm	36		J306H
 CT-3 22mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		J327H
	2-0	Violet	1 x 27"/70cm	36		J328H
	0	Violet	1 x 27"/70cm	36		J329H


Coated VICRYL® (polyglactin 910) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH-1 22mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12		J771D*
	4-0	Undyed	1 x 27"/70cm	36		J218H
	4-0	Violet	1 x 27"/70cm	36		J310H
	3-0	Violet	8 x 18"/45cm	12		J772D*
	3-0	Undyed	1 x 27"/70cm	36		J219H
	3-0	Violet	1 x 27"/70cm	36		J311H
 CT-2 26mm 1/2 Circle	3-0	Undyed	1 x 27"/70cm	36		J232H
	3-0	Violet	1 x 27"/70cm	36		J332H
	2-0	Violet	8 x 18"/45cm	12		J726D*
	2-0	Undyed	1 x 27"/70cm	36		J269H
	2-0	Violet	1 x 27"/70cm	36		J333H
	0	Violet	8 x 18"/45cm	12		J727D*
	0	Undyed	1 x 27"/70cm	36		J270H
	0	Violet	1 x 27"/70cm	36		J334H
	1	Violet	1 x 27"/70cm	36		J335H
 MO-6 26mm 1/2 Circle	2-0	Violet	8 x 18"/45cm	12		J789D*

*Safety organizer tray  CONTROL RELEASE™ Needle Suture


Coated VICRYL® (polyglactin 910) Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12	CR	J773D*
	4-0	Undyed	1 x 27"/70cm	36		J415H
	4-0	Violet	1 x 27"/70cm	36		J315H
	4-0	Violet	8 x 27"/70cm	12	CR	J783G
	3-0	Undyed	8 x 18"/45cm	12	CR	J864D
	3-0	Violet	8 x 18"/45cm	12	CR	J774D*
	3-0	Undyed	1 x 27"/70cm	36		J416H
	3-0	Violet	1 x 27"/70cm	36		J316H
	3-0	Violet	8 x 27"/70cm	12	CR	J784G
	2-0	Violet	8 x 18"/45cm	12	CR	J775D*
	2-0	Undyed	1 x 27"/70cm	36		J417H
	2-0	Violet	1 x 27"/70cm	36		J317H
	2-0	Violet	8 x 27"/70cm	12	CR	J785G
	0	Undyed	1 x 27"/70cm	36		J418H
 SH 26mm 1/2 Circle	3-0	Violet	1 x 36"/90cm	36		J527H
 MO-5 31mm 1/2 Circle	0	Violet	8 x 18"/45cm	12	CR	J801D*
	0	Violet	8 x 27"/70cm	12	CR	JJ81G
	1	Violet	8 x 27"/70cm	12	CR	JJ80G

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL[®] (polyglactin 910) Suture













Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	4-0	Violet	1 x 27"/70cm	36		J337H
	3-0	Undyed	8 x 18"/45cm	12	CR	J838D*
	3-0	Violet	8 x 18"/45cm	12	CR	J738D*
	3-0	Undyed	1 x 27"/70cm	36		J258H
	3-0	Violet	1 x 27"/70cm	36		J338H
	3-0	Undyed	1 x 36"/90cm	36		J944H
	3-0	Violet	1 x 36"/90cm	36		J344H
	2-0	Undyed	8 x 18"/45cm	12	CR	J839D*
	2-0	Violet	8 x 18"/45cm	12	CR	J739D*
	2-0	Undyed	1 x 27"/70cm	36		J259H
	2-0	Violet	1 x 27"/70cm	36		J339H
	2-0	Undyed	8 x 27"/70cm	12	CR	JJ42G
	2-0	Undyed	1 x 36"/90cm	36		J945H
	2-0	Violet	1 x 36"/90cm	36		J345H
	0	Undyed	8 x 18"/45cm	12	CR	J840D*
	0	Violet	8 x 18"/45cm	12	CR	J740D*
	0	Undyed	1 x 27"/70cm	36		J260H
	0	Violet	1 x 27"/70cm	36		J340H
	0	Violet	3 x 27"/70cm	24		J430T
	0	Undyed	8 x 27"/70cm	12	CR	JJ41G
	0	Violet	8 x 27"/70cm	12	CR	JJ31G
	0	Undyed	1 x 36"/90cm	36		J946H
	0	Violet	1 x 36"/90cm	36		J346H
	1	Undyed	8 x 18"/45cm	12	CR	J841D*
	1	Violet	8 x 18"/45cm	12	CR	J741D*
	1	Undyed	1 x 27"/70cm	36		J261H
	1	Violet	1 x 27"/70cm	36		J341H
	1	Undyed	8 x 27"/70cm	12	CR	JJ40G
	1	Undyed	1 x 36"/90cm	36		J947H
	1	Violet	1 x 36"/90cm	36		J347H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL[®] (polyglactin 910) Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>MH 36mm 1/2 Circle</p>	3-0	Violet	1 x 27"/70cm	36		J322H
	2-0	Violet	1 x 27"/70cm	36		J323H
 <p>MO-4 36mm 1/2 Circle</p>	2-0	Violet	8 x 18"/45cm	12		J700D*
	0	Violet	8 x 18"/45cm	12		J701D*
	0	Violet	1 x 27"/70cm	36		J436H
	1	Violet	8 x 18"/45cm	12		J702D*
	1	Violet	1 x 27"/70cm	36		J437H
 <p>CT 40mm 1/2 Circle</p>	3-0	Undyed	1 x 36"/90cm	36		J956H
	3-0	Violet	1 x 36"/90cm	36		J356H
	2-0	Violet	8 x 18"/45cm	12		J751D*
	2-0	Undyed	1 x 27"/70cm	36		J275H
	2-0	Violet	1 x 27"/70cm	36		J351H
	2-0	Undyed	1 x 36"/90cm	36		J957H
	2-0	Violet	1 x 36"/90cm	36		J357H
	0	Violet	8 x 18"/45cm	12		J752D*
	0	Undyed	1 x 27"/70cm	36		J280H
	0	Violet	1 x 27"/70cm	36		J352H
	0	Undyed	1 x 36"/90cm	36		J958H
	0	Violet	1 x 36"/90cm	36		J358H
	1	Violet	8 x 18"/45cm	12		J753D*
	1	Undyed	1 x 27"/70cm	36		J281H
	1	Violet	1 x 27"/70cm	36		J353H
	1	Undyed	1 x 36"/90cm	36		J959H
1	Violet	1 x 36"/90cm	36		J359H	

*Safety organizer tray  CONTROL RELEASE™ Needle Suture


Coated VICRYL® (polyglactin 910) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	3-0	Undyed	1 x 36"/90cm	36		J980H
	2-0	Undyed	8 x 18"/45cm	12	CR	J723D*
	2-0	Violet	1 x 27"/70cm	36		J363H
	2-0	Undyed	1 x 36"/90cm	36		J979H
	2-0	Violet	1 x 36"/90cm	36		J369H
	0	Undyed	8 x 18"/45cm	12	CR	J724D*
	0	Violet	8 x 18"/45cm	12	CR	J764D*
	0	Violet	1 x 27"/70cm	36		J364H
	0	Undyed	1 x 36"/90cm	36		J978H
	0	Violet	1 x 36"/90cm	36		J370H
	1	Violet	8 x 18"/45cm	12	CR	J765D*
	1	Violet	1 x 27"/70cm	36		J365H
	1	Undyed	1 x 36"/90cm	36		J977H
1	Violet	1 x 36"/90cm	36		J371H	
 <p>TP-1 65mm 1/2 Circle</p>	1	Violet	4 x 27"/70cm	12		J650G
	2	Undyed	2 x 27"/70cm	12		J849G
	2	Violet	4 x 27"/70cm	12		J649G
	2	Undyed	1 x 54"/135cm	24		J880T
 <p>XLH 70mm 1/2 Circle</p>	2-0	Violet	1 x 27"/70cm	12		J581G
	0	Violet	1 x 27"/70cm	12		J582G
	1	Violet	1 x 27"/70cm	12		J583G

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® (polyglactin 910) Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 UR-6 26mm 5/8 Circle	2-0	Violet	1 x 27"/70cm	36		J602H
	0	Violet	1 x 27"/70cm	36		J603H
 UR-5 36mm 5/8 Circle	2-0	Violet	1 x 27"/70cm	36		J375H
	0	Violet	1 x 27"/70cm	36		J376H
 UR-4 40mm 5/8 Circle	0	Violet	1 x 27"/70cm	36		J381H


TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V130-4 4.2mm 1/4 Circle	8-0	Violet	1 x 9"/23cm	12		J408G
 V130-5 5.8mm 1/4 Circle	8-0	Violet	1 x 9"/23cm	12		J409G
 V-34 36mm 1/2 Circle	3-0	Undyed	1 x 36"/90cm	36		J916H
	3-0	Violet	1 x 36"/90cm	36		J516H
	2-0	Undyed	1 x 36"/90cm	36		J917H
	2-0	Violet	1 x 36"/90cm	36		J517H
	0	Violet	1 x 36"/90cm	36		J518H
	1	Violet	1 x 36"/90cm	36		J519H







Coated VICRYL[®] (polyglactin 910) Suture



Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-3 13mm 3/8 Circle	5-0	Violet	1 x 18"/45cm	36		J385H
	4-0	Violet	1 x 18"/45cm	36		J386H
 FS-2 19mm 3/8 Circle	5-0	Violet	1 x 18"/45cm	36		J391H
	5-0	Undyed	1 x 27"/70cm	36		J421H
	4-0	Violet	1 x 18"/45cm	36		J392H
	4-0	Undyed	1 x 27"/70cm	36		J422H
	3-0	Violet	1 x 18"/45cm	36		J393H
	3-0	Undyed	1 x 27"/70cm	36		J423H


Coated VICRYL® (polyglactin 910) Suture







Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-1 24mm 3/8 Circle	4-0	Undyed	1 x 27"/70cm	36		J441H
	4-0	Violet	1 x 27"/70cm	36		J451H
	3-0	Undyed	1 x 27"/70cm	36		J442H
	2-0	Undyed	1 x 27"/70cm	36		J443H
 FS 26mm 3/8 Circle	3-0	Violet	8 x 18"/45cm	12	CR	J787D*
 FSL 30mm 3/8 Circle	3-0	Undyed	1 x 27"/70cm	36		J588H
	2-0	Undyed	1 x 27"/70cm	36		J589H
 J-1 17mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	36		J455H
 OS-4 22mm 1/2 Circle	2-0	Undyed	3 x 18"/45cm	24	CR	J748T
	2-0	Undyed	1 x 27"/70cm	36		J693H
	0	Undyed	3 x 18"/45cm	24	CR	J749T
	0	Violet	3 x 18"/45cm	24	CR	J707T
	0	Undyed	1 x 27"/70cm	36		J694H
	1	Violet	3 x 18"/45cm	24	CR	J708T
	1	Undyed	1 x 27"/70cm	36		J695H
 X-1 22mm 1/2 Circle	3-0	Violet	8 x 18"/45cm	12	CR	J790D*
	3-0	Undyed	1 x 27"/70cm	36		J458H
	3-0	Violet	1 x 27"/70cm	36		J460H
	2-0	Undyed	1 x 27"/70cm	36		J459H
	2-0	Violet	1 x 27"/70cm	36		J461H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® (polyglactin 910) Suture





Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CP-2 26mm 1/2 Circle</p>	3-0	Undyed	8 x 18"/45cm	12	CR	J761D*
	3-0	Undyed	1 x 27"/70cm	36		J868H
	2-0	Undyed	8 x 18"/45cm	12	CR	J762D*
	2-0	Undyed	1 x 27"/70cm	36		J869H
	0	Undyed	1 x 27"/70cm	36		J870H
	1	Undyed	1 x 27"/70cm	36		J871H
 <p>CP-1 36mm 1/2 Circle</p>	3-0	Violet	1 x 36"/90cm	36		J471H
	2-0	Undyed	1 x 27"/70cm	36		J266H
	2-0	Violet	1 x 27"/70cm	36		J466H
	2-0	Violet	1 x 36"/90cm	36		J472H
	0	Undyed	1 x 27"/70cm	36		J267H
	0	Violet	1 x 27"/70cm	36		J467H
	0	Violet	1 x 36"/90cm	36		J473H
	1	Undyed	1 x 27"/70cm	36		J268H
	1	Violet	1 x 27"/70cm	36		J468H
	1	Violet	1 x 36"/90cm	36		J474H
 <p>OS-6 36mm 1/2 Circle</p>	2-0	Undyed	3 x 18"/45cm	24	CR	J750T
	2-0	Violet	3 x 18"/45cm	24	CR	J710T
	2-0	Undyed	1 x 27"/70cm	36		J533H
	0	Undyed	3 x 18"/45cm	24	CR	J754T
	0	Violet	3 x 18"/45cm	24	CR	J711T
	0	Undyed	1 x 27"/70cm	36		J534H
	1	Undyed	1 x 27"/70cm	36		J535H
 <p>CP 40mm 1/2 Circle</p>	0	Violet	1 x 27"/70cm	36		J479H
	0	Violet	1 x 36"/90cm	36		J485H
	1	Undyed	1 x 27"/70cm	36		J196H
	1	Violet	1 x 27"/70cm	36		J480H
	1	Violet	1 x 36"/90cm	36		J486H
	2	Undyed	1 x 27"/70cm	36		J195H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


Coated VICRYL® (polyglactin 910) Suture

















Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>OS-8 40mm 1/2 Circle</p>	2-0	Violet	3 x 18"/45cm	24	CR	J716T
	0	Undyed	3 x 18"/45cm	24	CR	J756T
	0	Violet	3 x 18"/45cm	24	CR	J717T
	0	Undyed	1 x 27"/70cm	36		J698H
	1	Undyed	3 x 18"/45cm	24	CR	J757T
	1	Violet	3 x 18"/45cm	24	CR	J718T
	1	Undyed	1 x 27"/70cm	36		J699H
	2	Violet	3 x 18"/45cm	24	CR	J719T
 <p>CPX 48mm 1/2 Circle</p>	0	Undyed	1 x 27"/70cm	36		J568H
	1	Undyed	1 x 27"/70cm	36		J569H

Taperpoint ⊙, Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm/ FSL 30mm 1/2 Circle, 3/8 Circle</p>	3-0	Undyed	1 x 54"/135cm	36		J940H
	 <p>CT-136mm/ CP-136mm 1/2 Circle/ 1/2 Circle</p>	3-0	Violet	1 x 54"/135cm	36	
2-0		Violet	1 x 54"/135cm	36		J532H


Coated VICRYL® (polyglactin 910) Suture































Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 x 27"/70cm	36		J662H
	3-0	Undyed	1 x 27"/70cm	36		J663H
	3-0	Violet	1 x 27"/70cm	36		J523H


PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Violet	1 x 18"/45cm	12		J551G
 P-1 11mm 3/8 Circle	7-0	Undyed	1 x 18"/45cm	12		J488G
	6-0	Undyed	1 x 18"/45cm	12		J489G
	5-0	Undyed	1 x 18"/45cm	12		J490G
 P-3 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		J492G
	6-0	Undyed	1 x 18"/45cm	36		J492H
	5-0	Undyed	1 x 18"/45cm	12		J493G
	5-0	Undyed	1 x 18"/45cm	36		J493H
	5-0	Violet	1 x 18"/45cm	12		J463G
	4-0	Undyed	1 x 18"/45cm	12		J494G
 PSL 30mm 3/8 Circle	3-0	Undyed	1 x 27"/70cm	36		J502H
	2-0	Undyed	1 x 27"/70cm	36		J596H
	0	Violet	1 x 27"/70cm	36		J610H





Coated VICRYL® (polyglactin 910) Suture












PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-3 13mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J494H
	4-0	Violet	1 x 18"/45cm	12	 MultiPass® NEEDLES	J464G
 PS-3 16mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J499G
	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J500G
	5-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J500H
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J495G
	5-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J495H
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J496G
	4-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J496H
	4-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	J426H
	3-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J497G
	3-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J497H
	3-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	J427H
	2-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	J428H
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J682G
	4-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J682H
	4-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	J935H
	3-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	J683G
	3-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	J683H
	3-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	J936H
 P-2 8mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		J503G
	4-0	Undyed	1 x 18"/45cm	12		J504G
 PS-6 10mm 1/2 Circle	6-0	Undyed	1 x 18"/45cm	12		J510G
	5-0	Undyed	1 x 18"/45cm	12		J511G
 PS-5 13mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		J593G
	4-0	Undyed	1 x 18"/45cm	12		J594G

Coated VICRYL® (polyglactin 910) Suture











PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-4 16mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		J506G
	4-0	Undyed	1 x 18"/45cm	12		J507G
 PS-4C 16mm Compound Curve	4-0	Undyed	1 x 18"/45cm	12		J656G
	3-0	Undyed	1 x 18"/45cm	12		J657G




Cutting Edge Conventional ▲						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS-1 60mm 1/2 Circle	3	Undyed	1 x 27"/70cm	24		J808T



PRIME Conventional Cutting ▲						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		J833G
	5-0	Undyed	1 x 18"/45cm	12		J834G
	4-0	Undyed	1 x 18"/45cm	12		J835G
 PC-3 16mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		J843G
	5-0	Undyed	1 x 18"/45cm	12		J844G
	4-0	Undyed	1 x 18"/45cm	12		J845G

Coated VICRYL[®] (polyglactin 910) Suture











PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-5 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	 MultiPass [®] NEEDLES	J822G
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass [®] NEEDLES	J823G
	4-0	Undyed	1 x 18"/45cm	36	 MultiPass [®] NEEDLES	J823H
	3-0	Undyed	1 x 18"/45cm	12	 MultiPass [®] NEEDLES	J824G
 PC-12 16mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		J813G
	4-0	Undyed	1 x 18"/45cm	12		J814G

Spatula Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-29 7.6mm 1/4 Circle	6-0	Violet	1 x 12"/30cm	12		J556G
	6-0	Violet	1 x 18"/45cm	12		J555G
 S-28 7.6mm 1/2 Circle	6-0	Violet	1 x 18"/45cm	12		J562G

SABRELOC [®] Centerpoint Spatula 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RD-1 8mm 1/4 Circle	6-0	Violet	1 x 12"/30cm	12		J446G


Coated VICRYL[®] (polyglactin 910) Suture



Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	6-0	Violet	1 x 18"/45cm	12		J544G
 S-14 8mm 1/4 Circle	6-0	Violet	1 x 8"/20cm	12		J578G
	6-0	Violet	1 x 8"/20cm	12		J590G
	6-0	Violet	1 x 12"/30cm	12		J552G
	6-0	Undyed	1 x 18"/45cm	12		J670G
	6-0	Violet	1 x 18"/45cm	12		J570G
	5-0	Violet	1 x 8"/20cm	12		J579G
	5-0	Violet	1 x 8"/20cm	12		J591G
	5-0	Violet	1 x 12"/30cm	12		J553G
	5-0	Undyed	1 x 18"/45cm	12		J671G
	5-0	Violet	1 x 18"/45cm	12		J571G
 TG140-6 6.5mm 3/8 Circle	8-0	Violet	1 x 12"/30cm	12		J974G
 TG140-8 6.5mm 3/8 Circle	8-0	Violet	1 x 8"/20cm	12		J547G
	8-0	Violet	1 x 12"/30cm	12		J548G
	8-0	Violet	1 x 12"/30cm	12		J599G*
	7-0	Violet	1 x 12"/30cm	12		J566G
	7-0	Violet	1 x 18"/45cm	12		J546G
 TG160-6 5.5mm 1/2 Circle	8-0	Violet	1 x 12"/30cm	12		J975G
 TG160-8 5.5mm 1/2 Circle	8-0	Violet	1 x 8"/20cm	12		J574G
	7-0	Violet	1 x 12"/30cm	12		J575G
	7-0	Violet	1 x 18"/45cm	12		J576G

*Beaded retraction suture


Coated VICRYL® (polyglactin 910) Suture



SUTUPAK® Pre-Cut Sterile Sutures						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Undyed	3 x 18"/45cm	36		J643H
	4-0	Violet	3 x 18"/45cm	36		J633H
	4-0	Undyed	12 x 18"/45cm	24		J109T*
	4-0	Violet	12 x 18"/45cm	24		J103T*
	4-0	Violet	12 x 18"/45cm	12		J103G*
	3-0	Undyed	3 x 18"/45cm	36		J644H
	3-0	Violet	3 x 18"/45cm	36		J634H
	3-0	Undyed	12 x 18"/45cm	24		J110T*
	3-0	Undyed	12 x 18"/45cm	24		J910T
	3-0	Violet	12 x 18"/45cm	12		J104G*
	3-0	Violet	12 x 18"/45cm	24		J104T*
	3-0	Violet	12 x 18"/45cm	24		J904T
	3-0	Undyed	1 x 54"/135cm	36		J606H
	3-0	Violet	1 x 54"/135cm	36		J614H
	2-0	Undyed	3 x 18"/45cm	36		J645H
	2-0	Violet	3 x 18"/45cm	36		J635H
	2-0	Undyed	12 x 18"/45cm	24		J111T*
	2-0	Undyed	12 x 18"/45cm	24		J911T
	2-0	Violet	12 x 18"/45cm	12		J105G*
	2-0	Violet	12 x 18"/45cm	24		J105T*
	2-0	Violet	12 x 18"/45cm	24		J905T
	2-0	Undyed	1 x 54"/135cm	36		J607H
	2-0	Violet	1 x 54"/135cm	36		J615H
	0	Undyed	3 x 18"/45cm	36		J646H
	0	Violet	3 x 18"/45cm	36		J636H
	0	Undyed	6 x 18"/45cm	24		J112T*
	0	Violet	6 x 18"/45cm	12		J106G*
	0	Violet	6 x 18"/45cm	24		J106T*
	0	Undyed	12 x 18"/45cm	12		J912G
	0	Violet	12 x 18"/45cm	12		J906G
	0	Undyed	1 x 54"/135cm	36		J608H
	0	Violet	1 x 54"/135cm	36		J616H


*Single strand delivery

Coated VICRYL[®] (polyglactin 910) Suture




SUTUPAK [®] Pre-Cut Sterile Sutures						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	1	Undyed	6 x 18"/45cm	24		J113T*
	1	Undyed	12 x 18"/45cm	12		J913G
	1	Violet	12 x 18"/45cm	12		J907G
	1	Violet	1 x 54"/135cm	36		J617H
	2	Violet	1 x 54"/135cm	36		J618H



LIGAPAK [™] Dispensing Reel						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Undyed	1 x 54"/135cm	12		J283G
	4-0	Undyed	1 x 54"/135cm	12		J284G
	4-0	Violet	1 x 54"/135cm	12		J204G
	3-0	Undyed	1 x 54"/135cm	12		J285G
	3-0	Violet	1 x 54"/135cm	12		J205G
	3-0	Undyed	1 x 144"/366cm	12		J885G
	2-0	Undyed	1 x 54"/135cm	12		J286G
	2-0	Violet	1 x 54"/135cm	12		J206G
	0	Undyed	1 x 54"/135cm	12		J287G
	0	Violet	1 x 54"/135cm	12		J207G
	1	Violet	1 x 54"/135cm	12		J208G





VICRYL[®] (polyglactin 910) Suture


VICRYL® (polyglactin 910) Suture



Taperpoint 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 BV100-3 3.8mm 3/8 Circle	9-0	Violet/Monofilament	1 x 5"/13cm	12		V402G*

Spatula Reverse Cutting 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 TG140-8 6.5mm 3/8 Circle	9-0	Violet/Braided	1 x 12"/30cm	12		V549G*


ULTIMA® Reverse Cutting 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 CS90-6 7mm 1/4 Circle	10-0	Violet/Monofilament	1 x 4"/10cm	12		V980G*
 CS140-6 6.5mm 3/8 Circle	10-0	Violet/Monofilament	1 x 4"/10cm	12		V960G*
	9-0	Violet/Braided	1 x 4"/10cm	12		V950G*
 CS160-6 5.5mm 1/2 Circle	10-0	Violet/Monofilament	1 x 12"/30cm	12		V448G*
 CS160-8 5.5mm 1/2 Circle	10-0	Violet/Monofilament	1 x 12"/30cm	12		V450G*
	9-0	Violet/Monofilament	1 x 12"/30cm	12		V449G*
 CS-B-6 5mm Compound Curve	10-0	Violet/Monofilament	1 x 4"/10cm	12		V966G*




MONOCRYL[®] (poliglecaprone 25) Suture




MONOCRYL Suture is a monofilament synthetic absorbable surgical suture prepared from a copolymer of glycolide and epsilon-caprolactone.¹



MONOCRYL Sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.¹




MONOCRYL® (poliglecaprone 25) Suture








ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SHB 26mm 1/2 Circle	3-0	Undyed	1 x 27"/70cm	36		YB416
	3-0	Undyed	1 x 36"/90cm	36		YB944
 CTB-1 36mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		YB257
	0	Violet	1 x 27"/70cm	36		YB258
	0	Undyed	1 x 36"/90cm	36		YB946
	1	Undyed	1 x 36"/90cm	36		YB947
	0	Undyed	1 x 36"/90cm	36		YB978
CTXB 48mm 1/2 Circle	0	Undyed	1 x 36"/90cm	36		YB978

Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	6-0	Undyed	1 x 27"/70cm	36		Y432H
	5-0	Undyed	1 x 27"/70cm	36		Y433H
	4-0	Undyed	1 x 27"/70cm	36		Y434H


MONOCRYL® (poliglecaprone 25) Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	5-0	Undyed	1 x 27"/70cm	36		Y213H
	4-0	Undyed	1 x 27"/70cm	36		Y214H
	3-0	Undyed	1 x 27"/70cm	36		Y215H
	5-0	Violet	1 x 27"/70cm	36		Y303H
	4-0	Violet	1 x 27"/70cm	36		Y304H
	3-0	Violet	1 x 27"/70cm	36		Y305H
 SH-1 22mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	36		Y218H
	3-0	Undyed	1 x 27"/70cm	36		Y219H
	3-0	Undyed	8 x 18"/45cm	12	CR	Y872D*
 SH 26mm 1/2 Circle	4-0	Undyed	1 x 27"/70cm	36		Y415H
	4-0	Violet	1 x 27"/70cm	36		Y315H
	3-0	Undyed	1 x 27"/70cm	36		Y416H
	3-0	Violet	1 x 27"/70cm	36		Y316H
	2-0	Undyed	1 x 27"/70cm	36		Y417H
	2-0	Violet	1 x 27"/70cm	36		Y317H
	2-0	Violet	1 x 36"/90cm	36		Y732H
	0	Violet	1 x 27"/70cm	36		Y318H
 SH 26mm 1/2 Circle	3-0	Violet	1 x 36"/90cm	36		Y527H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


MONOCRYL® (poliglecaprone 25) Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-2 26mm 1/2 Circle</p>	2-0	Violet	1 x 27"/70cm	36		Y333H
	2-0	Violet	1 x 36"/90cm	36		Y762H
	0	Violet	1 x 27"/70cm	36		Y334H
	1	Violet	1 x 27"/70cm	36		Y335H
 <p>MO-5 31mm 1/2 Circle</p>	0	Violet	8 x 18"/45cm	12	CR	Y801D*
 <p>MH 36mm 1/2 Circle</p>	2-0	Undyed	1 x 27"/70cm	36		Y227H
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Violet	8 x 18"/45cm	12	CR	Y738D*
	3-0	Violet	1 x 27"/70cm	36		Y338H
	3-0	Undyed	1 x 36"/90cm	36		Y944H
	3-0	Violet	1 x 36"/90cm	36		Y344H
	2-0	Violet	8 x 18"/45cm	12	CR	Y739D*
	2-0	Violet	1 x 27"/70cm	36		Y339H
	2-0	Undyed	1 x 36"/90cm	36		Y945H
	2-0	Violet	1 x 36"/90cm	36		Y345H
	0	Violet	8 x 18"/45cm	12	CR	Y740D*
	0	Violet	1 x 27"/70cm	36		Y340H
	0	Undyed	1 x 36"/90cm	36		Y946H
	0	Violet	1 x 36"/90cm	36		Y346H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


MONOCRYL® (poliglecaprone 25) Suture







Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm 1/2 Circle	0	Violet	4 x 27"/70cm	12		YY31G
	1	Violet	1 x 27"/70cm	36		Y341H
	1	Violet	4 x 27"/70cm	12		YY30G
	1	Undyed	1 x 36"/90cm	36		Y947H
	1	Violet	1 x 36"/90cm	36		Y347H
 MO-4 36mm 1/2 Circle	0	Violet	8 x 18"/45cm	12		Y701D*
	0	Violet	4 x 27"/70cm	12		YY71G
	0	Violet	1 x 27"/70cm	24		Y436H
	1	Violet	1 x 27"/70cm	24		Y437H
	1	Violet	4 x 27"/70cm	12		YY70G
 CT 40mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		Y350H
	2-0	Violet	1 x 27"/70cm	36		Y351H
	2-0	Violet	1 x 36"/90cm	36		Y357H
	0	Violet	1 x 27"/70cm	36		Y352H
	0	Violet	1 x 36"/90cm	36		Y358H
	0	Violet	8 x 18"/45cm	24		Y752D
	1	Violet	1 x 27"/70cm	36		Y353H
	1	Violet	1 x 36"/90cm	36		Y359H
 CTX 48mm 1/2 Circle	2-0	Violet	1 x 36"/90cm	36		Y397H
	0	Violet	1 x 36"/90cm	36		Y398H
	1	Violet	1 x 36"/90cm	36		Y399H
 UR-6 26mm 5/8 Circle	3-0	Violet	1 x 27"/70cm	36		Y604H
	2-0	Violet	1 x 27"/70cm	36		Y605H
	0	Violet	1 x 27"/70cm	36		Y606H
 UR-5 36mm 5/8 Circle	2-0	Violet	1 x 27"/70cm	36		Y377H
	0	Violet	1 x 27"/70cm	36		Y378H




*Safety organizer tray CONTROL RELEASE™ Needle Suture



MONOCRYL® (poliglecaprone 25) Suture




TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-34 36mm 1/2 Circle	2-0	Undyed	1 x 27"/70cm	36		Y817H
	2-0	Undyed	1 x 36"/90cm	36		Y917H























Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Undyed	1 x 27"/70cm	36		Y266H
	0	Violet	1 x 36"/90cm	36		Y967H
	0	Undyed	1 x 27"/70cm	36		Y267H






Taperpoint  Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1, CP-1 36mm 1/2 Circle	3-0	Undyed	1 x 54"/135cm	36		Y931H
	2-0	Undyed	1 x 54"/135cm	36		Y932H

Straight Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Undyed	1 x 27"/70cm	36		Y523H


MONOCRYL® (poliglecaprone 25) Suture











PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y489G
	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y490G
 P-3 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y492G
	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y493G
	5-0	Violet	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y463G
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y494G
	4-0	Violet	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y464G
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y495G
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y496G
	4-0	Violet	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y513G
	4-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	Y426H
	3-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y497G
	3-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	Y427H
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	36	 MultiPass® NEEDLES	Y682H
	4-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	Y935H
	3-0	Undyed	1 x 27"/70cm	36	 MultiPass® NEEDLES	Y936H
 PS-6 10mm 1/2 Circle	6-0	Undyed	1 x 18"/45cm	12		Y510G

PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y833G
	5-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y834G
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	Y835G


MONOCRYL® (poliglecaprone 25) Suture



PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		Y844G
	4-0	Undyed	1 x 18"/45cm	12		Y845G
 PC-5 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		Y822G
	4-0	Undyed	1 x 18"/45cm	12		Y823G
 PC-12 16mm 1/2 Circle	4-0	Undyed	1 x 18"/45cm	12		Y814G

SUTUPAK® Pre-Cut Sterile Sutures						
No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Violet	6 x 18"/45cm	12		Y908G*
	2-0	Violet	6 x 18"/45cm	12		Y909G*
	0	Violet	6 x 18"/45cm	12		Y910G*
	1	Violet	6 x 18"/45cm	12		Y911G*

*Single strand delivery



PDS[®] II (polydioxanone) Suture


PDS II Suture is prepared from polyester, poly (p-dioxanone).¹




PDS II Sutures are indicated for use in soft tissue approximation, including use in pediatric cardiovascular tissue where growth is expected to occur and ophthalmic surgery. PDS II Suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue.¹





These sutures are particularly useful where the combination of an absorbable suture and extended wound support (up to 6 weeks) is desirable.¹




PDS® II (polydioxanone) Suture









ETHIGUARD® Blunt Point Needle 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTB-1 36mm 1/2 Circle</p>	2-0	Violet	1 x 36"/90cm	36		ZB345
	0	Violet	8 x 18"/45cm	12	CR	ZB740
	0	Violet	1 x 27"/70cm	36		ZB340
	0	Violet	1 x 36"/90cm	36		ZB346
	1	Violet	1 x 36"/90cm	36		ZB347
 <p>CTXB 48mm 1/2 Circle</p>	0	Violet	1 x 36"/90cm	24		ZB370
	1	Violet	1 x 36"/90cm	24		ZB371


Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>BV-1 9.3mm 3/8 Circle</p>	7-0	Violet	1 x 24"/60cm	36		Z155H
	7-0	Violet	1 x 30"/75cm	36		Z135H
	6-0	Violet	1 x 30"/75cm	36		Z117H
 <p>BV 11mm 3/8 Circle</p>	6-0	Violet	1 x 30"/75cm	36		Z123H
 <p>C-1 13mm 3/8 Circle</p>	6-0	Violet	1 x 30"/75cm	36		Z127H
	5-0	Violet	1 x 30"/75cm	36		Z126H





PDS® II (polydioxanone) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-2 13mm 1/2 Circle	6-0	Violet	1 x 30"/75cm	36		Z149H
	5-0	Violet	1 x 30"/75cm	36		Z148H
 TF 13mm 1/2 Circle	6-0	Violet	1 x 27"/70cm	36		Z432H
 UCL 13mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		Z114H
 RB-1 17mm 1/2 Circle	5-0	Violet	1 x 27"/70cm	36		Z303H
	4-0	Violet	1 x 27"/70cm	36		Z304H
	3-0	Violet	1 x 27"/70cm	36		Z305H
 RB-1 17mm 1/2 Circle	5-0	Violet	1 x 30"/75cm	36		Z320H
 SH-1 22mm 1/2 Circle	4-0	Violet	1 x 27"/70cm	36		Z310H
	3-0	Violet	1 x 27"/70cm	36		Z311H
 CT-2 26mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		Z332H
	2-0	Violet	1 x 27"/70cm	36		Z333H
	0	Violet	1 x 27"/70cm	36		Z334H





PDS® II (polydioxanone) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	4-0	Violet	8 x 18"/45cm	12	CR	Z773D*
	4-0	Violet	1 x 27"/70cm	36		Z315H
	3-0	Violet	8 x 18"/45cm	12	CR	Z774D*
	3-0	Undyed	1 x 27"/70cm	36		Z416H
	3-0	Violet	1 x 27"/70cm	36		Z316H
	2-0	Violet	8 x 18"/45cm	12	CR	Z775D*
	2-0	Violet	1 x 27"/70cm	36		Z317H
 CT-1 36mm 1/2 Circle	3-0	Violet	1 x 27"/70cm	36		Z338H
	2-0	Violet	8 x 18"/45cm	12	CR	Z739D*
	2-0	Undyed	1 x 27"/70cm	36		Z259H
	2-0	Violet	1 x 27"/70cm	36		Z339H
	0	Violet	8 x 18"/45cm	12	CR	Z740D*
	0	Violet	1 x 27"/70cm	36		Z340H
	0	Violet	1 x 36"/90cm	36		Z346H
	1	Violet	8 x 18"/45cm	12	CR	Z741D*
	1	Violet	1 x 27"/70cm	36		Z341H
	1	Violet	1 x 36"/90cm	36		Z347H
 CT 40mm 1/2 Circle	2-0	Violet	1 x 36"/90cm	36		Z357H
	0	Violet	1 x 27"/70cm	36		Z352H
	0	Violet	1 x 36"/90cm	24		Z358T
	1	Violet	1 x 27"/70cm	36		Z353H
	1	Violet	1 x 36"/90cm	24		Z359T

PDS® II (polydioxanone) Suture




Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTX 48mm 1/2 Circle	0	Violet	8 x 18"/45cm	12	CR	Z776D*
	0	Violet	1 x 36"/90cm	24		Z370T
	0	Violet	30"/75cm Loop Length	12	Looped	Z990G
	1	Violet	8 x 18"/45cm	12	CR	Z765D*
	1	Violet	1 x 36"/90cm	24		Z371T
 TP-1 65mm 1/2 Circle	0	Violet	30"/75cm Loop Length	12	Looped	Z991G
	1	Violet	2 x 27"/70cm	12		Z690G
	1	Violet	1 x 54"/135cm	12		Z879G
	1	Violet	48"/122cm Loop Length	12	Looped	Z880G
 XLH 70mm 1/2 Circle	1	Violet	1 x 27"/70cm	12		Z583G
	1	Violet	1 x 54"/135cm	12		Z877G
	1	Violet	48"/122cm Loop Length	12	Looped	Z881G
 ST-4 19mm Straight	4-0	Violet	1 x 14"/35cm	12		Z420G

TAPERCUT®						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-5 17mm 1/2 Circle	4-0	Violet	1 x 36"/90cm	36		Z137H





*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture

PDS® II (polydioxanone) Suture









TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-38 30mm 1/2 Circle</p>	0	Violet	1 x 18"/45cm	12		Z00JG


Trocar						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
<p>STP-10 254mm Straight</p>	2-0	Violet	1 x 27"/70cm	12		Z997G

Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>FS-2 19mm 3/8 Circle</p>	5-0	Undyed	1 x 27"/70cm	36		Z421H
	4-0	Undyed	1 x 27"/70cm	36		Z422H
	4-0	Violet	1 x 27"/70cm	36		Z397H
	3-0	Violet	1 x 27"/70cm	36		Z398H
	3-0	Undyed	1 x 27"/70cm	36		Z423H
 <p>FS-1 24mm 3/8 Circle</p>	4-0	Undyed	1 x 27"/70cm	36		Z441H
	3-0	Undyed	1 x 27"/70cm	36		Z442H
	3-0	Violet	1 x 27"/70cm	36		Z452G
	3-0	Violet	1 x 27"/70cm	36		Z452H
	2-0	Violet	1 x 27"/70cm	36		Z451G
	2-0	Violet	1 x 27"/70cm	36		Z451H
	2-0	Undyed	1 x 27"/70cm	36		Z443H
 <p>OS-4 22mm 1/2 Circle</p>	1	Violet	1 x 27"/70cm	36		Z695H

PDS® II (polydioxanone) Suture

























Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-2 26mm 1/2 Circle	3-0	Violet	8 x 18"/45cm	12	CR	Z761D*
	3-0	Violet	1 x 27"/70cm	36		Z968H
	2-0	Violet	8 x 18"/45cm	12	CR	Z762D*
	2-0	Violet	1 x 27"/70cm	36		Z969H
	0	Violet	1 x 27"/70cm	36		Z970H
 CP-1 36mm 1/2 Circle	2-0	Violet	1 x 27"/70cm	36		Z466H
	0	Violet	1 x 27"/70cm	36		Z467H
	1	Violet	1 x 27"/70cm	36		Z468H
 OS-6 36mm 1/2 Circle	0	Violet	1 x 27"/70cm	24		Z534T
	1	Violet	1 x 27"/70cm	36		Z535H
 CP 40mm 1/2 Circle	2	Violet	1 x 27"/70cm	24		Z195T
 CPX 48mm 1/2 Circle	1	Violet	1 x 27"/70cm	24		Z569T

Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Undyed	1 x 27"/70cm	36		Z662H
	3-0	Undyed	1 x 27"/70cm	36		Z663H


*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture



PDS® II (polydioxanone) Suture













PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		Z489G
	6-0	Violet	1 x 18"/45cm	12		Z487G
 P-3 13mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		Z493G
	5-0	Violet	1 x 18"/45cm	12		Z463G
	4-0	Undyed	1 x 18"/45cm	12		Z494G
 PS-2 19mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		Z495G
	4-0	Undyed	1 x 18"/45cm	12		Z496G
	4-0	Violet	1 x 18"/45cm	12		Z513G
	3-0	Undyed	1 x 18"/45cm	12		Z497G
 PS-1 24mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12		Z682G
	3-0	Undyed	1 x 18"/45cm	12		Z683G
 P-2 8mm 1/2 Circle	5-0	Undyed	1 x 18"/45cm	12		Z503G
	4-0	Undyed	1 x 18"/45cm	12		Z504G
 PS-4 16mm 1/2 Circle	4-0	Undyed	1 x 18"/45cm	12		Z507G

PDS® II (polydioxanone) Suture



Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS-1 60mm 1/2 Circle	2	Violet	1 x 27"/70cm	24		Z807T


PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Undyed	1 x 18"/45cm	12		Z833G
 PC-3 16mm 3/8 Circle	5-0	Undyed	1 x 18"/45cm	12		Z844G
	4-0	Undyed	1 x 18"/45cm	12		Z845G
 PC-5 19mm 3/8 Circle	4-0	Undyed	1 x 18"/45cm	12		Z823G
	3-0	Undyed	1 x 18"/45cm	12		Z824G




Surgical Gut Suture – Chromic



Surgical Gut Suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. Surgical Gut Sutures are available in plain or chromic.¹


Surgical Gut Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹




Surgical Gut Suture – Chromic





Blunt Point ●						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BP-1 65mm 1/2 Circle	0	Brown	1 x 27"/70cm	24		47T
 BP 80mm 1/2 Circle	1	Brown	1 x 27"/70cm	12		48G


Taperpoint ⊙						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13mm 3/8 Circle	5-0	Brown	1 x 27"/70cm	36		K895H
 BB 17mm 3/8 Circle	4-0	Brown	1 x 27"/70cm	12		4887G
 RB-1 17mm 1/2 Circle	5-0	Brown	1 x 27"/70cm	36		U202H
	4-0	Brown	1 x 27"/70cm	36		U203H
	3-0	Brown	1 x 27"/70cm	36		U204H





Surgical Gut Suture – Chromic




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-3 22mm 1/2 Circle	3-0	Brown	1 x 27"/70cm	36		892H
	2-0	Brown	1 x 27"/70cm	36		893H
 SH-1 22mm 1/2 Circle	4-0	Brown	1 x 27"/70cm	36		G181H
	3-0	Brown	1 x 27"/70cm	36		G182H
 CT-2 26mm 1/2 Circle	3-0	Brown	1 x 27"/70cm	36		882H
	2-0	Brown	1 x 27"/70cm	36		883H
	0	Brown	1 x 27"/70cm	36		884H
	1	Brown	1 x 27"/70cm	36		885H
 SH 26mm 1/2 Circle	4-0	Brown	1 x 27"/70cm	36		G121H
	3-0	Brown	1 x 27"/70cm	36		G122H
	2-0	Brown	1 x 27"/70cm	36		G123H
	0	Brown	1 x 27"/70cm	36		G124H
 SH 26mm 1/2 Circle	3-0	Brown	1 x 36"/90cm	36		G172H
	2-0	Brown	1 x 36"/90cm	36		G173H
 MO-5 31mm 1/2 Circle	0	Brown	8 x 18"/45cm	12	CR	CC04G
	1	Brown	8 x 18"/45cm	12	CR	CC03G
	1	Brown	1 x 27"/70cm	36		4995H






Surgical Gut Suture – Chromic




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	3-0	Brown	1 x 27"/70cm	36		810H
	3-0	Brown	1 x 36"/90cm	36		922H
	2-0	Brown	1 x 27"/70cm	36		811H
	2-0	Brown	3 x 27"/70cm	24		3811T
	2-0	Brown	1 x 36"/90cm	36		923H
	0	Brown	8 x 18"/45cm	12	CR	CC41G
	0	Brown	1 x 27"/70cm	36		812H
	0	Brown	3 x 27"/70cm	24		3812T
	0	Brown	1 x 36"/90cm	36		924H
	1	Brown	8 x 18"/45cm	12	CR	CC40G
	1	Brown	1 x 27"/70cm	36		813H
	1	Brown	3 x 27"/70cm	12		3813G
1	Brown	1 x 36"/90cm	36		925H	
2	Brown	1 x 27"/70cm	36		814H	
 <p>MH 36mm 1/2 Circle</p>	3-0	Brown	1 x 27"/70cm	36		G126H
	2-0	Brown	1 x 27"/70cm	36		G127H
	0	Brown	1 x 27"/70cm	36		G128H
 <p>MO-4 36mm 1/2 Circle</p>	0	Brown	8 x 18"/45cm	12	CR	CC02G
	0	Brown	1 x 27"/70cm	36		4924H
	1	Brown	8 x 18"/45cm	12	CR	CC01G
	1	Brown	1 x 27"/70cm	36		4925H
	1	Brown	3 x 27"/70cm	24		5589T
	1	Brown	8 x 27"/70cm	12	CR	CC38G
	2-0	Brown	1 x 27"/70cm	24		4923H







Surgical Gut Suture – Chromic




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT 40mm 1/2 Circle</p>	3-0	Brown	1 x 27"/70cm	36		800H
	3-0	Brown	1 x 36"/90cm	36		912H
	2-0	Brown	1 x 27"/70cm	36		801H
	2-0	Brown	1 x 36"/90cm	36		913H
	0	Brown	3 x 27"/70cm	24		3802T
	0	Brown	1 x 27"/70cm	36		802H
	0	Brown	1 x 36"/90cm	36		914H
	1	Brown	1 x 27"/70cm	36		803H
	1	Brown	3 x 27"/70cm	24		3803T
	1	Brown	1 x 36"/90cm	36		915H
 <p>CT 40mm 1/2 Circle</p>	3-0	Brown	1 x 54"/135cm	36		990H
	2-0	Brown	1 x 54"/135cm	36		991H
 <p>LH 40mm 1/2 Circle</p>	2-0	Brown	1 x 27"/70cm	36		G133H
 <p>MO-2 40mm 1/2 Circle</p>	1	Brown	1 x 27"/70cm	36		4915H




Surgical Gut Suture – Chromic






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	3-0	Brown	1 x 27"/70cm	36		862H
	2-0	Brown	1 x 27"/70cm	36		863H
	2-0	Brown	1 x 36"/90cm	36		903H
	0	Brown	1 x 27"/70cm	36		864H
	0	Brown	1 x 36"/90cm	36		904H
	1	Brown	1 x 27"/70cm	36		865H
	1	Brown	1 x 36"/90cm	36		905H
 <p>TP-1 65mm 1/2 Circle</p>	1	Brown	4 x 27"/70cm	12		M760G
	1	Brown	30"/75cm Loop Length	12	Looped	GL30G
	2	Brown	30"/75cm Loop Length	12	Looped	GL31G
 <p>XLH 70mm 1/2 Circle</p>	0	Brown	1 x 27"/70cm	24		44T
 <p>UR-6 26mm 5/8 Circle</p>	3-0	Brown	1 x 27"/70cm	36		N877H
	2-0	Brown	1 x 27"/70cm	36		N878H
	0	Brown	1 x 27"/70cm	36		N879H
 <p>UR-5 36mm 5/8 Circle</p>	2-0	Brown	1 x 27"/70cm	36		U245H
	1	Brown	1 X 27"/70cm	24		U247H
	0	Brown	1 x 27"/70cm	36		U246H


Surgical Gut Suture – Chromic











Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>UR-4 40mm 5/8 Circle</p>	2-0	Brown	1 x 27"/70cm	36		U255H
	0	Brown	1 x 27"/70cm	36		U256H
 <p>SH 26mm/ ST-1 45mm 1/2 Circle/ Straight</p>	3-0	Brown	1 x 36"/90cm	36		G152H


TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-5 17mm 1/2 Circle</p>	4-0	Brown	1 x 27"/70cm	36		U233H
	3-0	Brown	1 x 36"/90cm	36		942H
	2-0	Brown	1 x 27"/70cm	36		837H
	2-0	Brown	1 x 36"/90cm	36		943H
 <p>V-34 36mm 1/2 Circle</p>	0	Brown	1 x 27"/70cm	36		838H
	0	Brown	1 x 36"/90cm	36		944H
	1	Brown	1 x 27"/70cm	36		839H
	1	Brown	1 x 36"/90cm	36		945H



Surgical Gut Suture – Chromic




Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 M-2 13mm 1/4 Circle	4-0	Brown	1 x 12"/30cm	12		752G
 M-2 13mm 1/4 Circle	4-0	Brown	1 x 12"/30cm	12		751G
 C-3 13mm 3/8 Circle	5-0	Brown	1 x 12"/30cm	12		723G
 FS-2 19mm 3/8 Circle	5-0	Brown	1 x 27"/70cm	12		634G
	4-0	Brown	1 x 27"/70cm	36		635H
	3-0	Brown	1 x 27"/70cm	36		636H
 FS-1 24mm 3/8 Circle	2-0	Brown	1 x 27"/70cm	36		637H
 M-1 13mm 1/2 Circle	4-0	Brown	1 x 12"/30cm	12		744G
 J-1 17mm 1/2 Circle	4-0	Brown	1 x 18"/45cm	12		724G
	4-0	Brown	1 x 27"/70cm	36		U213H
 CP-2 26mm 1/2 Circle	3-0	Brown	1 x 27"/70cm	36		886H
	2-0	Brown	1 x 27"/70cm	36		887H
	0	Brown	1 x 27"/70cm	24		888H

Surgical Gut Suture – Chromic





















Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-1 36mm 1/2 Circle	2-0	Brown	1 x 27"/70cm	36		816H
	0	Brown	1 x 27"/70cm	36		817H
	1	Brown	1 x 36"/90cm	36		975H
 CP 40mm 1/2 Circle	3-0	Brown	1 x 27"/70cm	36		832H
	2-0	Brown	1 x 27"/70cm	36		833H
	1	Brown	1 x 27"/70cm	36		805H

Taperpoint ⊙, Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1, CP-1 36mm 1/2 Circle	3-0	Brown	1 x 54"/135cm	36		982H
	2-0	Brown	1 x 54"/135cm	36		983H
 CT, FSL 40mm 1/2 Circle	2-0	Brown	1 x 54"/135cm	36		977H


Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Brown	1 x 27"/70cm	36		654H











Surgical Gut Suture – Chromic




PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Brown	1 x 18"/45cm	12		794G
	6-0	Blue	1 x 18"/45cm	12		795G
 G-1 11mm 3/8 Circle	6-0	Brown	1 x 18"/45cm	12		796G
 G-1 11mm 3/8 Circle	6-0	Brown	1 x 18"/45cm	12		790G
 G-3 13mm 3/8 Circle	5-0	Brown	1 x 18"/45cm	12		792G
	4-0	Brown	1 x 18"/45cm	12		793G
 P-3 13mm 3/8 Circle	5-0	Brown	1 x 18"/45cm	12		687G
	4-0	Brown	1 x 18"/45cm	12		1654G
 PS-3 16mm 3/8 Circle	6-0	Brown	1 x 18"/45cm	12		1635G
	5-0	Brown	1 x 18"/45cm	12		1636G
 PS-2 19mm 3/8 Circle	4-0	Brown	1 x 18"/45cm	12		1637G
	3-0	Brown	1 x 27"/70cm	36		1638G
	2-0	Brown	1 x 27"/70cm	36		1621H
 PS 26mm 3/8 Circle	0	Brown	1 x 27"/70cm	36		598H







Surgical Gut Suture – Chromic





PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-2 8mm 1/2 Circle	4-0	Brown	1 x 18"/45cm	12		798G
 P-2 8mm 1/2 Circle	5-0	Brown	1 x 18"/45cm	12		1658G
	4-0	Brown	1 x 18"/45cm	12		1659G
 PS-6 10mm 1/2 Circle	6-0	Brown	1 x 18"/45cm	12		1816G
 PS-5 13mm 1/2 Circle	5-0	Brown	1 x 18"/45cm	12		1652G
	4-0	Brown	1 x 18"/45cm	12		1653G
 PS-4 16mm 1/2 Circle	6-0	Brown	1 x 18"/45cm	12		1641G
	5-0	Brown	1 x 18"/45cm	12		1642G
	4-0	Brown	1 x 18"/45cm	12		1643G
 PS-4C 16mm Compound Curve	4-0	Brown	1 x 18"/45cm	12		2644G
	3-0	Brown	1 x 18"/45cm	12		2643G


Surgical Gut Suture – Chromic



Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	7-0	Brown	1 x 18"/45cm	12		1744G
 S-14 8mm 1/4 Circle	6-0	Blue	1 x 18"/45cm	12		1791G
	5-0	Brown	1 x 12"/30cm	12		1766G
	5-0	Brown	1 x 18"/45cm	12		1792G
 S-4 8mm 1/4 Circle	4-0	Brown	1 x 18"/45cm	12		1793G
 TG140-8 6.5mm 3/8 Circle	10-0	Blue	1 x 6"/15cm	12		1757G
	10-0	Blue	1 x 12"/30cm	12		1771G
	9-0	Blue	1 x 6"/15cm	12		1754G
	7-0	Blue	1 x 18"/45cm	12		1797G
	6-0	Blue	1 x 18"/45cm	12		1731G
 S-2 8mm 1/2 Circle	4-0	Brown	1 x 18"/45cm	12		1798G

ULTIMA® Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CS175-8 7mm 1/2 Circle	7-0	Brown	1 x 18"/45cm	12		1745G

Surgical Gut Suture – Chromic



LIGAPAK™ Dispensing Reel


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Brown	1 x 54"/135cm	12		L111G
	3-0	Brown	1 x 54"/135cm	12		L112G
	2-0	Brown	1 x 54"/135cm	12		L113G
	0	Brown	1 x 54"/135cm	12		L114G
	1	Brown	1 x 54"/135cm	12		L115G

Pre-Cut Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Brown	1 x 54"/135cm	36		S111H
	3-0	Brown	1 x 54"/135cm	36		S112H
	2-0	Brown	1 x 54"/135cm	36		S113H
	0	Brown	1 x 54"/135cm	36		S114H
	1	Brown	1 x 54"/135cm	36		S115H
	2	Brown	1 x 54"/135cm	36		S116H
	3	Brown	1 x 54"/135cm	36		S117H

SUTUPAK® Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	2-0	Brown	12 x 18"/45cm	24		SG13T
	0	Brown	12 x 18"/45cm	24		SG14T
	1	Brown	12 x 18"/45cm	24		SG15T




Surgical Gut Suture - Plain









Surgical Gut Suture is an absorbable, sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. Surgical Gut Sutures are available in plain or chromic.¹

Surgical Gut Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological tissues.¹




Surgical Gut Suture – Plain






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RB-1 17mm 1/2 Circle	4-0	Yellowish Tan	1 x 27"/70cm	36		U207H
 CT-3 22mm 1/2 Circle	3-0	Yellowish Tan	1 x 27"/70cm	36		N862H
	2-0	Yellowish Tan	1 x 27"/70cm	36		N863H
	0	Yellowish Tan	1 x 27"/70cm	36		N864H
 SH 26mm 1/2 Circle	4-0	Yellowish Tan	1 x 27"/70cm	36		G321H
	3-0	Yellowish Tan	1 x 27"/70cm	36		G322H
 CT-1 36mm 1/2 Circle	3-0	Yellowish Tan	3 x 18"/45cm	24		5503T
	3-0	Yellowish Tan	1 x 27"/70cm	36		842H
	2-0	Yellowish Tan	1 x 27"/70cm	36		843H
	0	Yellowish Tan	1 x 27"/70cm	36		844H
 CT 40mm 1/2 Circle	3-0	Yellowish Tan	8 x 18"/45cm	12	CR	CP53G
	3-0	Yellowish Tan	1 x 27"/70cm	36		852H
	2-0	Yellowish Tan	1 x 27"/70cm	36		853H
 CTX 48mm 1/2 Circle	3-0	Yellowish Tan	1 x 27"/70cm	36		873H
	2-0	Yellowish Tan	1 x 27"/70cm	36		872H
 XLH 70mm 1/2 Circle	3-0	Yellowish Tan	1 x 27"/70cm	24		52T
	2-0	Yellowish Tan	1 x 27"/70cm	24		53T


Surgical Gut Suture – Plain






TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-7 26mm 1/2 Circle	3-0	Yellowish Tan	1 x 18"/45cm	12		849G
 V-34 36mm 1/2 Circle	3-0	Yellowish Tan	1 x 27"/70cm	36		846H


Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Yellowish Tan	1 x 27"/70cm	12		H820G
	4-0	Yellowish Tan	1 x 27"/70cm	36		H821H
	3-0	Yellowish Tan	1 x 27"/70cm	36		H822H
 FS-1 24mm 3/8 Circle	3-0	Yellowish Tan	1 x 27"/70cm	36		H810H
	2-0	Yellowish Tan	1 x 27"/70cm	36		H823H
 M-1 13mm 1/2 Circle	4-0	Yellowish Tan	1 x 12"/30cm	12		742G
 X-1 22mm 1/2 Circle	3-0	Brown	1 x 18"/45cm	12		612G
 CP-1 36mm 1/2 Circle	2-0	Yellowish Tan	1 x 27"/70cm	36		898H

Surgical Gut Suture – Plain



















Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-1 11mm 3/8 Circle	6-0	Yellowish Tan	1 x 18"/45cm	12		774G
 G-3 13mm 3/8 Circle	5-0	Yellowish Tan	1 X 18"/45cm	12		772G
	4-0	Yellowish Tan	1 X 18"/45cm	12		773G

Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SC-1 13mm Straight	4-0	Yellowish Tan	1 x 18"/45cm	36		1824H
 SC-1 13mm Straight	4-0	Yellowish Tan	1 x 18"/45cm	36		1828H
 KS 60mm Straight	3-0	Yellowish Tan	1 x 27"/70cm	36		665H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-6 8mm 1/4 Circle	6-0	Blue	1 x 18"/45cm	12		775G
 G-1 11mm 3/8 Circle	6-0	Yellowish Tan	1 x 18"/45cm	12		770G
 P-1 11mm 3/8 Circle	6-0	Yellowish Tan	1 x 18"/45cm	12		1657G


Surgical Gut Suture – Plain



PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-3 13mm 3/8 Circle	5-0	Yellowish Tan	1 x 18"/45cm	12		686G
	4-0	Yellowish Tan	1 x 18"/45cm	12		1644G
 PS-3 16mm 3/8 Circle	5-0	Yellowish Tan	1 x 18"/45cm	12		1626G
 PS-2 19mm 3/8 Circle	4-0	Yellowish Tan	1 x 18"/45cm	36		1627G
	3-0	Yellowish Tan	1 x 27"/70cm	36		1630H
 PS-1 24mm 3/8 Circle	3-0	Yellowish Tan	1 x 27"/70cm	36		1640H
 P-2 8mm 1/2 Circle	5-0	Yellowish Tan	1 x 18"/45cm	12		658G
 PS-4 16mm 1/2 Circle	5-0	Yellowish Tan	1 x 18"/45cm	12		1632G
	4-0	Yellowish Tan	1 x 18"/45cm	12		1633G

PRIME Conventional Cutting  Fast-Absorbing Sutures						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 5-0 Yellowish Tan 1 x 18"/45cm 12 	5-0	Yellowish Tan	1 x 18"/45cm	12		1915G

Surgical Gut Suture – Plain



Spatula Reverse Cutting

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	5-0	Yellowish Tan	1 x 12"/30cm	12		1775G
	5-0	Yellowish Tan	1 x 18"/45cm	12		1772G
 TG140-8 6.5mm 3/8 Circle	6-0	Yellowish Tan	1 x 18"/45cm	12		1735G

LIGAPAK™ Dispensing Reel


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Yellowish Tan	1 x 54"/135cm	12		L101G
	3-0	Yellowish Tan	1 x 54"/135cm	12		L102G
	2-0	Yellowish Tan	1 x 54"/135cm	12		L103G
	0	Yellowish Tan	1 x 54"/135cm	12		L104G

Pre-Cut Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Yellowish Tan	12 x 18"/45cm	24		SG2T
	3-0	Yellowish Tan	1 x 54"/135cm	36		S102H
	2-0	Yellowish Tan	1 x 54"/135cm	36		S103H
	0	Yellowish Tan	1 x 54"/135cm	36		S104H
	1	Yellowish Tan	1 x 54"/135cm	36		S105H

Sutures:

Nonabsorbable




PRONOVA[®] Poly (Hexafluoropropylene - VDF) Suture
















PRONOVA Suture (clear or pigmented) is a nonabsorbable, sterile surgical suture made from polymer blend of poly (vinylidene fluoride) and poly (vinylidene fluoride-co-hexafluoropropylene). The suture is pigmented blue to enhance visibility.¹

PRONOVA Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular and neurological procedures.¹




PRONOVA[®] Poly (Hexafluoropropylene - VDF) Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	2 x 24"/60cm	12		M3739*
 BV175-6 8mm 3/8 Circle	8-0	Blue	2 x 24"/60cm	12		M3743*
	7-0	Blue	4 x 24"/60cm	12		M3735*
	7-0	Blue	2 x 24"/60cm	12		M3737*
	7-0	Blue	1 x 24"/60cm	36		3735H*
 BV175-7 8mm 3/8 Circle	8-0	Blue	2 x 24"/60cm	12		M3767*
	7-0	Blue	4 x 30"/75cm	12		M3770*
	7-0	Blue	2 x 24"/60cm	12		M3766*
	7-0	Blue	2 x 18"/45cm	12		M3768*
 BV175-8 9.3mm 3/8 Circle	7-0	Blue	2 x 24"/60cm	12		M3747*

*Straight pack

PRONOVA® Poly (Hexafluoropropylene - VDF) Suture






Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	7-0	Blue	4 x 30"/75cm	12		M3703*
	7-0	Blue	1 x 30"/75cm	36		3703H
	7-0	Blue	4 x 24"/60cm	12		M3702*
	7-0	Blue	1 x 24"/60cm	36		3702H
	7-0	Blue	1 x 24"/60cm	36	VB	3304H
	7-0	Blue	1 x 18"/45cm	36		3701H
	6-0	Blue	4 x 30"/75cm	12		M3709*
	6-0	Blue	1 x 30"/75cm	36		3709H
	6-0	Blue	4 x 24"/60cm	12		M3805*
	6-0	Blue	1 x 24"/60cm	36		3805H
	6-0	Blue	1 x 24"/60cm	36	VB	3305H
	6-0	Blue	1 x 18"/45cm	36		3806H
 BV 11mm 3/8 Circle	6-0	Blue	1 x 30"/75cm	36		3776H
 C-1 13mm 3/8 Circle	6-0	Blue	4 x 30"/75cm	12		M3706*
	6-0	Blue	1 x 30"/75cm	36		3706H
	6-0	Blue	1 x 30"/75cm	36	VB	3307H
	6-0	Blue	1 x 24"/60cm	36		3726H
	5-0	Blue	4 x 36"/90cm	12		M3720*
	5-0	Blue	1 x 36"/90cm	36		3720H
	5-0	Blue	1 x 36"/90cm	36	VB	3321H
	5-0	Blue	4 x 24"/60cm	12		M3725*
	5-0	Blue	1 x 24"/60cm	36		3725H
	5-0	Blue	1 x 18"/45cm	36		3717H

*Straight pack


VB VISI-BLACK™ Surgical Needle






PRONOVA® Poly (Hexafluoropropylene - VDF) Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BB 17mm 3/8 Circle	5-0	Blue	1 x 36"/90cm	36		3580H
	4-0	Blue	1 x 36"/90cm	36		3581H
 RB-2 13mm 1/2 Circle	6-0	Blue	1 x 30"/75cm	36		3711H
	5-0	Blue	1 x 30"/75cm	36		3710H
 RB-1 17mm 1/2 Circle	5-0	Blue	1 x 36"/90cm	36		3556H
	5-0	Blue	1 x 24"/60cm	36		3555H
	4-0	Blue	1 x 36"/90cm	36		3557H
	3-0	Blue	1 x 36"/90cm	36		3558H
 SH-1 22mm 1/2 Circle	4-0	Blue	1 x 30"/75cm	36		3526H
	3-0	Blue	1 x 30"/75cm	36		3762H
 SH 26mm 1/2 Circle	4-0	Blue	1 x 36"/90cm	36		3521H
	3-0	Blue	1 x 36"/90cm	36		3522H
	2-0	Blue	1 x 36"/90cm	36		3523H

PRONOVA® Poly (Hexafluoropropylene - VDF) Suture



TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MH 36mm 1/2 Circle	3-0	Blue	1 x 36"/90cm	36		3842H
 CC-1 13mm 3/8 Circle	6-0	Blue	1 x 30"/75cm	36		3707H
 V-7 26mm 1/2 Circle	4-0 3-0	Blue Blue	1 x 36"/90cm 1 x 36"/90cm	36 36		3975H 3976H



ETHIBOND EXCEL® Polyester Suture

ETHIBOND EXCEL Suture is a nonabsorbable, braided, sterile, surgical suture composed of Poly (ethylene terephthalate). It is prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component.¹

The suture is uniformly coated with polybutylate or poly {oxy-1, 4 butanediyl oxy (1, 6-dioxo-1, 6 hexanediyl)}. The highly adherent coating is a relatively nonreactive nonabsorbable compound which acts as a lubricant to mechanically improve the physical properties of the uncoated suture by improving handling qualities as contrasted to the braided, uncoated fiber.¹






ETHIBOND EXCEL Suture is braided for optimal handling properties, and dyed green for enhanced visibility in the surgical field.¹

ETHIBOND EXCEL Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹



ETHIBOND EXCEL® Polyester Suture










Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13mm 3/8 Circle	6-0	Green	1 x 30"/75cm	36		X889H
	5-0	Green	1 x 30"/75cm	36		X890H
 C-1 13mm 3/8 Circle	6-0	Green	1 x 18"/45cm	36		X706H
	5-0	Green	1 x 18"/45cm	36		X717H
	5-0	Green	1 x 30"/75cm	36		X727H
	4-0	Green	1 x 30"/75cm	36	Pledget*	X726H
 BB 17mm 3/8 Circle	4-0	Green	1 x 30"/75cm	36		X881H
 BB 17mm 3/8 Circle	4-0	Green	1 x 30"/75cm	36		X581H
	3-0	Green	1 x 36"/90cm	36		X588H
 TF 13mm 1/2 Circle	5-0	Green	1 x 24"/60cm	36		X205H
	4-0	Green	1 x 24"/60cm	36		X204H
 UCL 13mm 1/2 Circle	2-0	Green	1 x 27"/70cm	36		X114H

*7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

ETHIBOND EXCEL® Polyester Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>RB-1 17mm 1/2 Circle</p>	5-0	Green	1 x 30"/75cm	36		X550H
	5-0	Green	1 x 36"/90cm	36		X556H
	4-0	Green	1 x 30"/75cm	36		X551H
	4-0	Green	4 x 30"/75cm	12		MX551
	4-0	Green	1 x 36"/90cm	36		X557H
	3-0	Green	4 x 30"/75cm	12		MX552
	3-0	Green	1 x 36"/90cm	36		X558H
	2-0	Green	4 x 30"/75cm	12		MX553
	2-0	White	4 x 30"/75cm	12		MX523
	2-0	Grn/Wht	8 x 30"/75cm	12	Pledget*	PXX41
	2-0	Grn/Wht	8 x 30"/75cm	12	Pledget†	PXX43
	2-0	Grn/Wht	10 x 30"/75cm	12		10X42
2-0	Grn/Wht	10 x 30"/75cm	12	Pledget‡	SXX43	
 <p>RB-1 17mm 1/2 Circle</p>	5-0	Green	1 x 30"/75cm	36		X870H
	4-0	Green	1 x 30"/75cm	36		X871H
	3-0	Green	1 x 30"/75cm	36		X872H
	2-0	Green	1 x 30"/75cm	36		X873H
 <p>SH-2 20mm 1/2 Circle</p>	2-0	Green	1 x 30"/75cm	36	Pledget*	PX82H
	2-0	Green	1 x 30"/75cm	36		X582H
	2-0	White	1 x 30"/75cm	36	Pledget*	PX83H
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX80
	2-0	Grn/Wht	10 x 30"/75cm	12		PXX82
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget§	SXP82
	2-0	Green	1 x 36"/90cm	36		X559H
 <p>MO-7 22mm 1/2 Circle</p>	2-0	Green	8 x 18"/45cm	12		CX42D
	0	Green	8 x 18"/45cm	12		CX41D

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget

†3mm x 3mm x 1.5mm Firm TFE Polymer Pledget






‡3mm x 3mm x 1.5mm Soft TFE Polymer Pledget

§7mm x 3mm x 1.5mm Soft TFE Polymer Pledget

^{||}Safety organizer tray

ETHIBOND EXCEL® Polyester Suture







Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SH-1 22mm 1/2 Circle</p>	4-0	Green	1 x 36"/90cm	36		X761H
	3-0	Green	1 x 36"/90cm	36		X762H
	2-0	Green	4 x 30"/75cm	12		MX763
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX92
	2-0	Green	1 x 36"/90cm	36		X763H
 <p>CT-2 26mm 1/2 Circle</p>	2-0	Green	8 x 18"/45cm	12	CR	CX26D [†]
	2-0	Green	1 x 30"/75cm	36		X411H
	0	Green	8 x 18"/45cm	12	CR	CX27D [†]
	0	Green	1 x 30"/75cm	36		X412H
 <p>MO-6 26mm 1/2 Circle</p>	2-0	Green	8 x 18"/45cm	12	CR	CX46D [†]
	0	Green	8 x 18"/45cm	12	CR	CX45D [†]
 <p>SH 26mm 1/2 Circle</p>	4-0	Green	8 x 18"/45cm	12	CR	CX14D [†]
	4-0	Green	1 x 30"/75cm	36		X831H
	3-0	Green	8 x 18"/45cm	12	CR	CX13D [†]
	3-0	Green	1 x 30"/75cm	36		X832H
	2-0	Green	8 x 18"/45cm	12	CR	CX12D [†]
	2-0	Green	1 x 30"/75cm	36		X833H
	2-0	Green	4 x 30"/75cm	12		MX833
	2-0	White	4 x 30"/75cm	12		MX823
	0	Green	1 x 30"/75cm	36		X834H

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
[†]Safety organizer tray

ETHIBOND EXCEL® Polyester Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>SH 26mm 1/2 Circle</p>	4-0	Green	1 x 36"/90cm	36		X521H
	3-0	Green	1 x 30"/75cm	36		X562H
	3-0	Green	1 x 36"/90cm	36		X522H
	2-0	Green	1 x 30"/75cm	36		X563H
	2-0	Green	4 x 30"/75cm	12		MX563
	2-0	Grn/Wht	10 x 30"/75cm	12		10X82
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX76
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	SXP76
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget†	SXX86
	2-0	Green	1 x 36"/90cm	12	Pledget*	PX514
	2-0	Green	1 x 36"/90cm	36	Pledget*	SXP02
	2-0	Green	1 x 36"/90cm	36		X523H
	2-0	White	1 x 36"/90cm	36		X513H
	2-0	Grn/Wht	6 x 36"/90cm	12	Pledget*	PXX85
	2-0	Grn/Wht	10 x 36"/90cm	12		10X86
	2-0	Grn/Wht	10 x 36"/90cm	12	Pledget*	PXX86
	2-0	Grn/Wht	10 x 36"/90cm	12	Pledget*	SXP86
	0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX87
0	Green	1 x 36"/90cm	36		X524H	
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Green	8 x 18"/45cm	12	CR	CX22D‡
	2-0	Green	1 x 30"/75cm	36		X423H
	0	Green	8 x 18"/45cm	12	CR	CX21D‡
	0	Green	1 x 30"/75cm	36		X424H
	1	Green	1 x 30"/75cm	36		X425H
 <p>MH 36mm 1/2 Circle</p>	2-0	Green	1 x 36"/90cm	36		X843H
	0	Green	1 x 36"/90cm	36		X844H




*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
 †3mm x 3mm x 1.5mm Firm TFE Polymer Pledget





‡Safety organizer tray

CR CONTROL RELEASE™ Needle Suture

ETHIBOND EXCEL® Polyester Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	0	Green	8 x 18"/45cm	12	CR	CX31D*
	0	Green	1 x 30"/75cm	36		X864H
	1	Green	8 x 18"/45cm	12	CR	CX30D*
	1	Green	1 x 30"/75cm	36		X865H
 <p>XLH 70mm 1/2 Circle</p>	1	Green	1 x 30"/75cm	12		X190G

TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CC-1 13mm 3/8 Circle</p>	5-0	Green	1 x 30"/75cm	36		X728H
 <p>V-4 17mm 3/8 Circle</p>	4-0	Green	1 x 30"/75cm	36		X942H

ETHIBOND EXCEL® Polyester Suture






TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-4 17mm 3/8 Circle	4-0	Green	1 x 30"/75cm	36		X925H
	3-0	Green	1 x 30"/75cm	36		X926H
	2-0	Green	1 x 30"/75cm	36		X927H
	2-0	Green	4 x 30"/75cm	12		MX927
	2-0	White	4 x 30"/75cm	12		MX944
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX22
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	SXP22
 V-5 17mm 1/2 Circle	5-0	Green	1 X 30"/75cm	36		X951H
 V-7 26mm 1/2 Circle	3-0	Green	1 x 30"/75cm	36		X963H
	2-0	Green	1 x 30"/75cm	36		X964H
 V-34 36mm 1/2 Circle	0	Green	1 x 30"/75cm	36		X444H
 V-40 48mm 1/2 Circle	5	Green	4 x 30"/75cm	12		MB46G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget

ETHIBOND EXCEL® Polyester Suture



TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-5 17mm 1/2 Circle</p>	5-0	Green	1 x 30"/75cm	36		X934H
	4-0	Green	1 x 30"/75cm	36		X935H
	3-0	Green	1 x 30"/75cm	36		X936H
	3-0	White	1 x 30"/75cm	36		X916H
	3-0	Grn/Wht	8 x 30"/75cm	12	Pledget†	SXX53
	3-0	Grn/Wht	10 x 30"/75cm	12		10X53
	3-0	Green	1 x 36"/90cm	36		X806H
	2-0	Green	1 x 30"/75cm	36	Pledget*	PX937
	2-0	Green	1 x 30"/75cm	36		X937H
	2-0	White	1 x 30"/75cm	36		X917H
	2-0	Green	4 x 30"/75cm	12	Pledget*	MPX37
	2-0	White	4 x 30"/75cm	12	Pledget*	MPX17
	2-0	Grn/Wht	8 x 30"/75cm	12	Pledget*	PXX50
	2-0	Grn/Wht	8 x 30"/75cm	12	Pledget†	SXX50
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX52
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	SXP52
	2-0	Grn/Wht	10 x 30"/75cm	12		10X52
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget†	SXP54
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget†	SXX54
	2-0	Green	1 x 36"/90cm	36	Pledget*	SX54H
	2-0	Green	1 x 36"/90cm	36	Pledget*	PX54H
	2-0	Green	1 x 36"/90cm	36		X807H
	2-0	White	1 x 36"/90cm	36	Pledget†	SX17H
2-0	White	1 x 36"/90cm	36	Pledget*	PX17H	
 <p>V-7 26mm 1/2 Circle</p>	4-0	Green	1 x 30"/75cm	36		X975H
	3-0	Green	1 x 30"/75cm	36		X976H
	2-0	Green	1 x 30"/75cm	36	Pledget*	PX977
	2-0	Green	1 x 30"/75cm	36	Pledget*	SXP07
	2-0	Green	1 x 30"/75cm	36		X977H
	2-0	White	1 x 30"/75cm	36		X987H
	2-0	Green	4 x 30"/75cm	12		MX977



*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget
 †3mm x 3mm x 1.5mm Firm TFE Polymer Pledget

ETHIBOND EXCEL® Polyester Suture



TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-7 26mm 1/2 Circle	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	PXX77
	2-0	Grn/Wht	10 x 30"/75cm	12	Pledget*	SXP77
	2-0	Grn/Wht	10 x 30"/75cm	12		10X72
	2-0	Green	1 x 36"/90cm	36	Pledget*	PX62H
	2-0	Green	1 x 36"/90cm	36		X966H
	2-0	White	1 x 36"/90cm	36	Pledget*	PX64H
	2-0	Grn/Wht	8 x 36"/90cm	12	Pledget*	PXX79
	2-0	Grn/Wht	10 x 36"/90cm	12		10X66
	0	Green	1 x 36"/90cm	36		X905H
	0	White	1 x 36"/90cm	36	Pledget*	PX65H
 V-37 40mm 1/2 Circle	2	Green	4 x 30"/75cm	12		MX69G
	5	Green	4 x 30"/75cm	12		MB66G








Trocar						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
STP-10 250mm Straight	2-0	Green	1 x 27"/70cm	12		X997G

Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	4-0	Green	1 x 18"/45cm	12		X683G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget


ETHIBOND EXCEL® Polyester Suture




Reverse Cutting						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-1 24mm 3/8 Circle	3-0	Green	1 x 18"/45cm	36		X663H
 FS 26mm 3/8 Circle	2-0	Green	1 x 18"/45cm	36		X664H
 LR 75mm 3/8 Circle	2	Green	1 x 20"/50cm	24		X406T
	5	Green	1 x 20"/50cm	24		B409T
 LR 75mm 3/8 Circle	2	Green	1 x 30"/75cm	24		X496T
	5	Green	1 x 30"/75cm	24		B499T
 OS-2 17mm 1/2 Circle	0	Green	3 x 18"/45cm	24		CX74T
	1	Green	3 x 18"/45cm	24		CX75T
 OS-4 22mm 1/2 Circle	0	Green	1 x 30"/75cm	36		X517H
	1	Green	1 x 30"/75cm	36		X518H
	2	Green	1 x 30"/75cm	36		X519H

ETHIBOND EXCEL® Polyester Suture



Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 OS-6 36mm 1/2 Circle	1	Green	1 x 30"/75cm	36		X538H
 OS-8 40mm 1/2 Circle	1	Green	1 x 30"/75cm	36		X548H

Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	3-0	Green	1 x 30"/75cm	36		X622H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-3 13mm 3/8 Circle	5-0	Green	1 x 18"/45cm	12		X698G
 PS-3 16mm 3/8 Circle	5-0	Green	1 x 18"/45cm	12		X681G
 PS-2 19mm 3/8 Circle	4-0	White	1 x 18"/45cm	12		X692G
 PS-4 16mm 1/2 Circle	4-0	White	1 x 18"/45cm	12		X695G

ETHIBOND EXCEL® Polyester Suture




Cutting Edge Conventional

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CCS 48mm 1/2 Circle</p>	5	Green	4 x 30"/75cm	12		MB47G

SUTUPAK® Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Green	12 x 18"/45cm	36		X182H
	4-0	Green	12 x 18"/45cm	36		X183H
	4-0	Green	12 x 30"/75cm	36		X303H
	3-0	Green	12 x 18"/45cm	36		X184H
	3-0	Green	12 x 30"/75cm	36		X304H
	2-0	Green	12 x 18"/45cm	36		X185H
	2-0	Green	12 x 30"/75cm	36		X305H
	0	Green	6 x 18"/45cm	36		X186H
	0	Green	6 x 30"/75cm	36		X306H




ETHILON[®] Nylon Suture


















ETHILON Suture is nonabsorbable, sterile surgical monofilament suture composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6,6. ETHILON Suture is dyed black or green to enhance visibility in tissue.¹

ETHILON Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹




ETHILON® Nylon Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-3 3.8mm 3/8 Circle	10-0	Black	1 x 5"/13cm	12		2820G
	9-0	Black	1 x 5"/13cm	12		2819G
	8-0	Black	1 x 5"/13cm	12		2822G
 BV50-3 3.8mm 3/8 Circle	11-0	Black	1 x 5"/13cm	12		2881G
 BV75-3 3.8mm 3/8 Circle	10-0	Black	1 x 5"/13cm	12		2870G
 BV130-4 5mm 3/8 Circle	9-0	Black	1 x 5"/13cm	12		2813G
	8-0	Black	1 x 5"/13cm	12		2815G
 BV100-4 5.1mm 3/8 Circle	10-0	Black	1 x 5"/13cm	12		2830G
	9-0	Black	1 x 5"/13cm	12		2829G
 BV75-4 5.1mm 3/8 Circle	10-0	Black	1 x 5"/13cm	12		2850G
 BV130-5 6.5mm 3/8 Circle	10-0	Black	1 x 5"/13cm	12		2810G
	9-0	Black	1 x 5"/13cm	12		2809G
	8-0	Black	1 x 5"/13cm	12		2808G
 BV130-5 6.5mm 3/8 Circle	9-0	Black	1 x 10"/25cm	12		2800G
 BVH100-3 5mm 1/2 Circle	10-0	Black	1 x 5"/13cm	12		2860G
 CT 40mm 1/2 Circle	0	Black	48"/122cm Loop Length	24	Looped	L886T

ETHILON® Nylon Suture











Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTX 48mm 1/2 Circle	1	Black	1 x 30"/75cm	36		830H
 TP-1 65mm 1/2 Circle	0	Black	48"/122cm Loop Length	12	Looped	L880G
	1	Black	1 x 60"/150cm	12		824G
	2	Black	1 x 60"/150cm	12		825G
 ST100-3 3mm Straight	10-0	Black	1 x 5"/13cm	12		2849G

TAPERCUT®						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V100-3 3.5mm 1/4 Circle	9-0	Black	1 x 5"/13cm	12		2894G
 V130-3 3.8mm 1/4 Circle	9-0	Black	1 x 5"/13cm	12		2898G
	8-0	Black	1 x 5"/13cm	12		2897G
 V75-3 3.5mm 1/4 Circle	10-0	Black	1 x 5"/13cm	12		2892G


Reverse Cutting						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-2 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		667G





ETHILON® Nylon Suture




Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-3 13mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		668G
 FS-2 19mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		661G
	5-0	Black	1 x 18"/45cm	36		661H
	4-0	Black	1 x 18"/45cm	12		662G
	4-0	Black	1 x 18"/45cm	36		662H
 FS-1 24mm 3/8 Circle	4-0	Black	1 x 18"/45cm	36		1629H
	3-0	Black	1 x 18"/45cm	12		663G
	3-0	Black	1 x 18"/45cm	36		663H
	3-0	Black	1 x 30"/75cm	36		669H
 FS 26mm 3/8 Circle	2-0	Black	1 x 18"/45cm	12		664G
	2-0	Black	1 x 18"/45cm	36		664H
 FSL 30mm 3/8 Circle	4-0	Black	1 x 30"/75cm	36		1670H
	3-0	Black	1 x 30"/75cm	36		1671H
 PSL 30mm 3/8 Circle	4-0	Black	1 x 30"/75cm	36		562H
	3-0	Black	1 x 30"/75cm	36		1691H
 FSLX 36mm 3/8 Circle	3-0	Black	1 x 30"/75cm	36		1673H
	2-0	Black	1 x 30"/75cm	36		1674H

ETHILON® Nylon Suture





































Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 LR 75mm 3/8 Circle	2	Black	1 x 20"/50cm	12		470G*
	2	Black	1 x 20"/50cm	24		460T
 LR 75mm 3/8 Circle	1	Black	1 x 30"/75cm	24		489T
	2	Black	1 x 30"/75cm	24		490T
 X-1 22mm 1/2 Circle	3-0	Black	1 x 18"/45cm	12		642G
 PSLX 36mm 3/8 Circle	3-0	Black	1 x 30"/75cm	36		1683H
	2-0	Black	1 x 30"/75cm	36		1697H

Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Black	1 x 30"/75cm	36		626H
	3-0	Black	1 x 30"/75cm	36		627H
	2-0	Black	1 x 30"/75cm	36		628H


*Lead shot and bolster

ETHILON® Nylon Suture

















PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-6 8mm 1/4 Circle	7-0	Black	1 x 18"/45cm	12		1647G
 P-1 11mm 3/8 Circle	7-0	Black	1 x 18"/45cm	12		1696G
	6-0	Undyed	1 x 18"/45cm	12		689G
	6-0	Black	1 x 18"/45cm	12		697G
	6-0	Black	1 x 18"/45cm	36		697H
	6-0	Green	1 x 18"/45cm	12		G697G
	5-0	Green	1 x 18"/45cm	12		G695G
 P-3 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		1698G
	6-0	Green	1 x 18"/45cm	12		G696G
	5-0	Undyed	1 x 18"/45cm	12		690G
	5-0	Black	1 x 18"/45cm	12		698G
	5-0	Black	1 x 18"/45cm	36		698H
	5-0	Green	1 x 18"/45cm	12		G698G
	4-0	Undyed	1 x 18"/45cm	12		691G
	4-0	Black	1 x 18"/45cm	12		699G
	4-0	Black	1 x 18"/45cm	36		699H
 PS-3 16mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		1665G
	5-0	Black	1 x 18"/45cm	12		1668G
	5-0	Black	1 x 18"/45cm	36		1668H
 PS-2 19mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		1666G
	5-0	Black	1 x 18"/45cm	36		1666H
	5-0	Green	1 x 18"/45cm	12		G666G
	4-0	Undyed	1 x 18"/45cm	12		1611G
	4-0	Black	1 x 18"/45cm	12		1667G
	4-0	Black	1 x 18"/45cm	36		1667H
	4-0	Green	1 x 18"/45cm	12		G667G
	3-0	Black	1 x 18"/45cm	36		1669H
	2-0	Black	1 x 18"/45cm	36		593H

ETHILON® Nylon Suture





PRIME Reverse Cutting ▽						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PS-1 24mm 3/8 Circle	3-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1663G
	3-0	Black	1 x 18"/45cm	36	 MultiPass® NEEDLES	1663H
 PS 26mm 3/8 Circle	2-0	Black	1 x 18"/45cm	36	 MultiPass® NEEDLES	585H
 PS-4 16mm 1/2 Circle	6-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1660G
	5-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1661G
	4-0	Undyed	1 x 18"/45cm	12	 MultiPass® NEEDLES	1603G
	4-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1662G

PRIME Conventional Cutting ▲						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1856G*
	6-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1956G*
	5-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1855G*
	5-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1955G*
	4-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1854G*
	4-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1954G*
 PC-3 16mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1866G*
	6-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1966G*
	5-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1865G*
	5-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1965G*
	4-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1864G*
	4-0	Black	1 x 18"/45cm	12	 MultiPass® NEEDLES	1964G*


*Pliabilized

ETHILON® Nylon Suture




PRIME Conventional Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-5 19mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		1895G*
	5-0	Black	1 x 18"/45cm	12		1995G*
	4-0	Black	1 x 18"/45cm	12		1894G*
	4-0	Black	1 x 18"/45cm	12		1994G*
	3-0	Black	1 x 18"/45cm	12		1893G*
	3-0	Black	1 x 18"/45cm	12		1993G*
 PC-12 16mm 1/2 Circle	5-0	Black	1 x 18"/45cm	12		1845G
	5-0	Black	1 x 18"/45cm	12		1945G


SABRELOC® Centerpoint Spatula 








Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 RD-1 8mm 1/4 Circle	5-0	Black	1 x 18"/45cm	12		749G

Spatula Reverse Cutting 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG100-8 6.5mm 1/4 Circle	8-0	Black	1 x 12"/30cm	12		1714G
 S-14 8mm 1/4 Circle	5-0	Black	1 x 18"/45cm	12		7731G
 VAS100-4 5mm 3/8 Circle	9-0	Black	1 x 6"/15cm	12		2890G


ETHILON® Nylon Suture












Spatula Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TG140-6 6.5mm 3/8 Circle	10-0	Black	1 x 12"/30cm	12		7770G
 TG140-8 6.5mm 3/8 Circle	10-0	Black	1 x 12"/30cm	12		7718G
	9-0	Black	1 x 12"/30cm	12		7717G
 TG160-4 5.3mm 1/2 Circle	10-0	Black	1 x 12"/30cm	12		7707G
 TG160-4-3M 4.3mm 1/2 Circle	10-0	Black	1 x 8"/20cm	12		7711G
 TG160-6 5.5mm 1/2 Circle	10-0	Black	1 x 12"/30cm	12		7756G
	10-0	Black	1 x 6"/15cm	12		7757G
	9-0	Black	1 x 12"/30cm	12		7760G
 TG175-8 7mm 1/2 Circle	8-0	Black	1 x 12"/30cm	12		1716G
 S-24 8mm 1/2 Circle	5-0	Black	1 x 18"/45cm	12		1740G

ULTIMA® Reverse Cutting 

Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CS44-6 6.5mm 1/8 Circle	10-0	Black	1 x 4"/10cm	12		9065G
 CS140-6 6.5mm 3/8 Circle	10-0	Black	1 x 4"/10cm	12		9060G

ETHILON[®] Nylon Suture

ULTIMA [®] Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CSB-6 5mm Compound Curve	10-0	Black	1 x 12"/30cm	12		9007G
 CS140-6 6.5mm 3/8 Circle	9-0	Black	1 x 12"/30cm	12		9013G
	10-0	Black	1 x 12"/30cm	12		9003G
 CS140-8 6.5mm 3/8 Circle	10-0	Black	1 x 12"/30cm	12		9033G
 CSM-6 4.3mm 1/2 Circle	10-0	Black	1 X 8"/20cm	12		9006G
	10-0	Black	1 x 12"/30cm	12		9005G
 CS160-6 5.5mm 1/2 Circle	10-0	Black	1 x 4"/10cm	12		9061G
 CS160-6 5.5mm 1/2 Circle	11-0	Black	1 X 8"/20cm	12		9016G
	10-0	Black	1 X 6"/15cm	12		9001G
	10-0	Black	1 x 12"/30cm	12		9000G
	9-0	Black	1 x 12"/30cm	12		9011G
 CS160-8 5.5mm 1/2 Circle	10-0	Black	1 x 6"/15cm	12		9031G
	10-0	Black	1 x 12"/30cm	12		9030G
 CS175-6 7mm 1/2 Circle	10-0	Black	1 x 12"/30cm	12		9002G
 CS175-8 7mm 1/2 Circle	10-0	Black	1 x 12"/30cm	12		9032G
	9-0	Black	1 x 12"/30cm	12		9036G
	8-0	Black	1 x 12"/30cm	12		9038G
 CS-C-6 5mm Compound Curve	10-0	Black	1 x 6"/15cm	12		9010G



MERSILENE® Polyester Fiber Suture




MERSILENE Suture is a nonabsorbable, braided, sterile surgical suture composed of Poly (ethylene terephthalate). It is prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component.¹

MERSILENE Sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹



MERSILENE® Polyester Fiber Suture









Taperpoint 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 SH 26mm 1/2 Circle	3-0	Green/Braided	1 x 30"/75cm	36		R832H
	2-0	Green/Braided	1 x 30"/75cm	36		R833H
 CT-1 36mm 1/2 Circle	0	Green/Braided	1 x 30"/75cm	36		R424H
	0	Green/Braided	8 x 18"/45cm	24		MR21T

Reverse Cutting 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	4-0	White/Braided	1 x 18"/45cm	36		R633H
 FS-1 24mm 3/8 Circle	3-0	White/Braided	1 x 18"/45cm	36		R647H
 FS 26mm 3/8 Circle	2-0	White/Braided	1 x 18"/45cm	36		R665H

MERSILENE® Polyester Fiber Suture



Reverse Cutting 					
Needle	Suture USP	Color/Construction	Length	EA/BX	Product Code
 <p>LR 75mm 3/8 Circle</p>	1	Green/Braided	1 x 30"/75cm	24	R495T







PRIME Reverse Cutting 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>P-3 13mm 3/8 Circle</p>	5-0	White/Braided	1 x 18"/45cm	12		R690G
	4-0	White/Braided	1 x 18"/45cm	12		R691G

SABRELOC® Centerpoint Spatula 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>RD-1 8mm 1/4 Circle</p>	5-0	Green/Braided	1 x 18"/45cm	12		740G
	5-0	White/Braided	1 x 18"/45cm	12		746G



MERSILENE® Polyester Fiber Suture



Spatula Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	6-0	White/Braided	1 x 18"/45cm	12		1722G
	5-0	White/Braided	1 x 18"/45cm	12		1760G
	5-0	Green/Braided	1 x 18"/45cm	12		1761G
	5-0	White/Braided	1 x 18"/45cm	12		1764G
 S-24 8mm 1/4 Circle	4-0	Green/Braided	1 x 18"/45cm	12		1762G
 S-4 8mm 1/4 Circle	4-0	White/Braided	1 x 18"/45cm	12		1778G
 TG160-4 5.3mm 1/2 Circle	11-0	Green/ Monofilament	1 x 12"/30cm	12		R750G
 TG160-4-3M 4.3mm 1/2 Circle	11-0	Green/Monofilament	1 x 8"/20cm	12		R749G
 S-2 8mm 1/2 Circle	4-0	White/Braided	1 x 18"/45cm	12		1779G


ULTIMA® Reverse Cutting

Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 CS140-6 6.5mm 3/8 Circle	10-0	Green/Monofilament	1 x 12"/30cm	12		R770G
 CS160-6 5.5mm 1/2 Circle	10-0	Green/Monofilament	1 x 12"/30cm	12		R756G

MERSILENE® Polyester Fiber Strip



Taperpoint 						
Needle	Suture USP	Color/Construction	Length	EA/BX	Special	Product Code
 <p>CTX 48mm 1/2 Circle</p>	5mm Tape	White	1 x 12"/30cm	6		RS22
 <p>MO-4 36mm 1/2 Circle</p>	5mm Tape	White	1 x 12"/30cm	6		RS23




NUROLON[®] Nylon Suture























NUROLON Suture is a nonabsorbable sterile surgical braided suture composed of the long-chain aliphatic polymer Nylon 6,6. It is dyed black to enhance visibility in tissue. The suture is also available undyed (clear).¹

NUROLON Suture is indicated for use in general soft tissue approximating and/or ligation, including use in cardiovascular, ophthalmic procedures, and neurological procedures.¹




NUROLON® Nylon Suture





Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	4-0	Black	8 x 18"/45cm	12		C584D*
	4-0	Black	8 x 18"/45cm	24		N104T
 RB-1 17mm 1/2 Circle	4-0	Black	8 x 18"/45cm	12		C554D*
	4-0	Black	8 x 18"/45cm	24		N124T
	3-0	Black	8 x 18"/45cm	12		C553D*
 SH-1 22mm 1/2 Circle	3-0	Black	8 x 18"/45cm	12		C503D*
	2-0	Black	8 x 18"/45cm	12		C502D*
 MO-7 22mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12		C542D*
	0	Black	8 x 18"/45cm	12		C541D*
 SH 26mm 1/2 Circle	4-0	Black	8 x 18"/45cm	12		C514D*
	3-0	Black	8 x 18"/45cm	12		C513D*
	2-0	Black	8 x 18"/45cm	12		C512D*
 CT-2 26mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12		C526D*
	0	Black	8 x 18"/45cm	12		C527D*
 MO-6 26mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12		C546D*
	0	Black	8 x 18"/45cm	12		C545D*

*Safety organizer tray  CONTROL RELEASE™ Needle Suture

NUROLON[®] Nylon Suture




Taperpoint						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Black	8 x 18"/45cm	12		C522D*
	0	Black	1 x 30"/75cm	36		5424H
	0	Black	8 x 18"/45cm	12		C521D*
	1	Black	1 x 30"/75cm	36		5425H
	1	Black	8 x 18"/45cm	12		C520D*
 <p>CTX 48mm 1/2 Circle</p>	0	Black	8 x 18"/45cm	12		C551D*
	1	Black	8 x 18"/45cm	12		C550D*




Reverse Cutting						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>FS 26mm 3/8 Circle</p>	4-0	Black	8 x 18"/45cm	12		C570D*
 <p>FSL 30mm 3/8 Circle</p>	3-0	Black	8 x 18"/45cm	12		C566D*
 <p>X-1 22mm 1/2 Circle</p>	3-0	Black	8 x 18"/45cm	12		C563D*


*Safety organizer tray CONTROL RELEASE™ Needle Suture

NUROLON® Nylon Suture



Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CP-2 26mm 1/2 Circle	3-0	Black	8 x 18"/45cm	12	CR	C576D*
	2-0	Black	8 x 18"/45cm	12	CR	C575D*
 CP 40mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12	CR	C572D*
 OS-8 40mm 1/2 Circle	1	Black	3 x 18"/45cm	24	CR	C595T

PRIME Conventional Cutting ▲						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12	 MultiPass NEEDLES	5666G
	5-0	Black	1 x 18"/45cm	12	 MultiPass NEEDLES	5665G



PERMA-HAND® Silk Suture


PERMA-HAND Suture is a nonabsorbable, sterile, surgical suture composed of an organic protein called fibroin. This protein is derived from the domesticated species *Bombyx mori* (B. mori) of the family Bombycidae.¹









The sutures are processed to remove the natural waxes and gums, dyed black and coated with a special wax mixture. PERMA-HAND Suture is also available in its natural color, also available is PERMA-HAND Virgin Silk Suture in which the sericin gum is not removed and serves to hold the filaments together.¹

PERMA-HAND Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.¹













PERMA-HAND[®] Silk Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV-1 9.3mm 3/8 Circle	6-0	Black	1 x 18"/45cm	36		K802H
 BV-1 9.3mm 3/8 Circle	7-0	Black	1 x 18"/45cm	36	Siliconized silk	K809H
	7-0	Black	1 x 18"/45cm	36		Q809H
	6-0	Black	1 x 18"/45cm	36		K801H
 C-1 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		707G
	6-0	Black	1 x 30"/75cm	36		K889H
	5-0	Black	1 x 30"/75cm	36		K890H
	4-0	Black	1 x 30"/75cm	36		K891H
 C-1 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12	Siliconized silk	706G
	6-0	Black	1 x 18"/45cm	36		Q706H
	4-0	Black	1 x 18"/45cm	36		709H
 BB 17mm 3/8 Circle	5-0	Black	1 x 30"/75cm	36		K880H
	4-0	Black	1 x 30"/75cm	36		K881H
	3-0	Black	1 x 30"/75cm	36		K882H
	2-0	Black	1 x 30"/75cm	36		K883H
 BB 17mm 3/8 Circle	5-0	Black	1 x 24"/60cm	36		K580H
	4-0	Black	1 x 24"/60cm	36		K581H
 RB-2 13mm 1/2 Circle	4-0	Black	1 x 18"/45cm	36		N267H


PERMA-HAND[®] Silk Suture








Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	5-0	Black	1 x 18"/45cm	36		N266H
	4-0	Black	8 x 18"/45cm	12	CR	C084D*
	4-0	Black	8 x 18"/45cm	24		M104T
 RB-1 17mm 1/2 Circle	5-0	Black	1 x 30"/75cm	36		K870H
	4-0	Black	8 x 18"/45cm	12	CR	C054D*
	4-0	Black	1 x 30"/75cm	36		K871H
	3-0	Black	8 x 18"/45cm	12	CR	C053D*
	3-0	Black	1 x 30"/75cm	36		K872H
 RB-1 17mm 1/2 Circle	5-0	Black	1 x 24"/60cm	36		K570H
	4-0	Black	1 x 24"/60cm	36		K571H
 RB-1 17mm 1/2 Circle	3-0	Black	1 x 24"/60cm	36		K572H
	2-0	Black	1 x 30"/75cm	36		K873H
 CT-3 22mm 1/2 Circle	2-0	Black	1 x 20"/50cm	36		413H
 MO-7 22mm 1/2 Circle	0	Black	8 x 18"/45cm	12	CR	C041D*
 SH-1 22mm 1/2 Circle	4-0	Black	8 x 18"/45cm	12	CR	C004D*
	3-0	Black	8 x 18"/45cm	12	CR	C003D*
	3-0	Black	8 x 18"/45cm	24		M3T
 CT-2 26mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12	CR	C026D*
	0	Black	8 x 18"/45cm	12	CR	C027D*

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


PERMA-HAND[®] Silk Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MO-6 26mm 1/2 Circle	0	Black	8 x 18"/45cm	12	CR	C045D*
 SH 26mm 1/2 Circle	4-0	Black	5 x 18"/45cm	24	CR	C0145
	4-0	Black	8 x 18"/45cm	12	CR	C014D*
	4-0	Black	1 x 30"/75cm	36		K831H
	4-0	Black	8 x 30"/75cm	12	CR	C018D*
	3-0	Black	5 x 18"/45cm	12	CR	C0135
	3-0	Black	8 x 18"/45cm	12	CR	C013D*
	3-0	Black	8 x 18"/45cm	24		M13T
	3-0	Black	1 x 30"/75cm	36		K832H
	3-0	Black	8 x 30"/75cm	12	CR	C017D*
	2-0	Black	5 x 18"/45cm	24	CR	C0125
	2-0	Black	8 x 18"/45cm	12	CR	C012D*
	2-0	Black	1 x 30"/75cm	36		K833H
	2-0	Black	8 x 30"/75cm	12	CR	C016D*
	0	Black	1 x 30"/75cm	36		K834H
	1	Black	1 x 30"/75cm	36		K835H
 SH 26mm 1/2 Circle	2-0	Black	1 x 24"/60cm	36		K533H
	0	Black	1 x 24"/60cm	36		K534H
 CT-1 36mm 1/2 Circle	2-0	Black	8 x 18"/45cm	12	CR	C022D*
	2-0	Black	1 x 30"/75cm	36		423H
	0	Black	8 x 18"/45cm	12	CR	C021D*
	0	Black	1 x 30"/75cm	36		424H

*Safety organizer tray **CR** CONTROL RELEASE™ Needle Suture


PERMA-HAND® Silk Suture















Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 MH 36mm 1/2 Circle	2-0	Black	8 x 24"/60cm	24		M72T
	2-0	Black	1 x 30"/75cm	36		K843H
	0	Black	1 x 30"/75cm	36		K844H
	1	Black	1 x 30"/75cm	36		K845H
	2	Black	1 x 30"/75cm	36		K846H
 TF-5 8mm 5/8 Circle	5-0	Black	1 x 18"/45cm	36		N271H
	4-0	Black	1 x 18"/45cm	36		N272H
 ST-1 45mm Straight	3-0	Black	1 x 30"/75cm	36		K852H


TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 V-4 17mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		K945G
 V-5 17mm 1/2 Circle	4-0	Black	1 x 30"/75cm	36		K952H
 V-7 26mm 1/2 Circle	2-0	Black	1 x 30"/75cm	36		K964H








PERMA-HAND[®] Silk Suture




Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13 mm 3/8 Circle	4-0	Black	1 x 12"/30cm	12		7734G
 C-2 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		714G
 C-3 13mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		736G
 C-3 13mm 3/8 Circle	4-0	Black	1 x 12"/30cm	12		735G
 CC-1 13mm 3/8 Circle	6-0	Black	1 x 12"/30cm	12		7729G
 FS-2 19mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		682G
	4-0	Black	1 x 18"/45cm	12		683G
 FS-1 24mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		629G
	3-0	Black	1 x 18"/45cm	12		684G
	3-0	Black	1 x 18"/45cm	12		684H
 FS 26mm 3/8 Circle	2-0	Black	1 x 18"/45cm	12		685G
	2-0	Black	1 x 18"/45cm	12		685H
 FSL 30mm 3/8 Circle	3-0	Black	1 x 30"/75cm	36		675H
	2-0	Black	1 x 18"/45cm	12		677G
	2-0	Black	1 x 30"/75cm	36		679H
	0	Black	1 x 18"/45cm	12		678G
	0	Black	1 x 18"/45cm	12		678H
	0	Black	1 x 30"/75cm	36		680H
 FSLX 36mm 3/8 Circle	2-0	Black	1 x 30"/75cm	36		1689H

PERMA-HAND® Silk Suture














Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PSLX 36mm 3/8 Circle	2-0	Black	1 x 30"/75cm	36		591H
	2	Black	1 x 30"/75cm	24		482T
	5	Black	1 x 30"/75cm	24		485T
 J-1 17mm 1/2 Circle	4-0	Black	1 x 18"/45cm	12		734G
 X-1 22mm 1/2 Circle	3-0	Black	1 x 18"/45cm	12		632G
	2-0	Black	1 x 18"/45cm	12		737G
 CP-1 36mm 1/2 Circle	2-0	Black	1 x 30"/75cm	36		443H
 G-7 8mm 1/2 Circle	7-0	Black	1 x 18"/45cm	12		766G
	6-0	Black	1 x 18"/45cm	12		765G
 G-6 8mm 1/4 Circle	7-0	Black	1 x 18"/45cm	12		768G
	6-0	Black	1 x 18"/45cm	12		769G


PERMA-HAND® Silk Suture











Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 G-1 11mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		786G
 G-1 11mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		780G
 G-3 13mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		789G
 G-3 13mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		783G




Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 KS 60mm Straight	4-0	Black	1 x 30"/75cm	36		621H
	3-0	Black	1 x 30"/75cm	36		622H
	2-0	Black	1 x 30"/75cm	36		623H
	0	Black	1 x 30"/75cm	36		624H

PRIME Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-1 11mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		639G
 P-3 13mm 3/8 Circle	6-0	Black	1 x 18"/45cm	12		1639G
	5-0	Black	1 x 18"/45cm	12		640G
	4-0	Black	1 x 18"/45cm	12		641G
 PS-3 16mm 3/8 Circle	6-0	Black	1 x 18"/45cm	36		1675H
	5-0	Black	1 x 18"/45cm	36		1678H
 PS-2 19mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		1676G
	4-0	Black	1 x 18"/45cm	12		1677G
	3-0	Black	1 x 18"/45cm	36		1679H
	2-0	Black	1 x 18"/45cm	36		583H
 PS-1 24mm 3/8 Circle	3-0	Black	1 x 18"/45cm	12		1684G


PERMA-HAND[®] Silk Suture

PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PSL 30mm 3/8 Circle	2-0	Black	1 x 18"/45cm	36		673H
	2-0	Black	1 x 30"/75cm	36		486H
	0	Black	1 x 18"/45cm	36		580H
	0	Black	1 x 30"/75cm	36		590H
 PS 26mm 3/8 Circle	2-0	Black	1 x 18"/45cm	36		1588H
 PS-4 16mm 1/2 Circle	4-0	Black	1 x 18"/45cm	12		1687G

PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-3 16mm 3/8 Circle	5-0	Black	1 x 18"/45cm	12		1925G
 PC-5 19mm 3/8 Circle	4-0	Black	1 x 18"/45cm	12		1984G

Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 S-14 8mm 1/4 Circle	6-0	Black	1 x 18"/45cm	12		1780G
 TG140-8 6.5mm 3/8 Circle	7-0	Black	1 x 18"/45cm	12		7733G
	6-0	Black	1 x 18"/45cm	12		1732G

PERMA-HAND® Silk Suture




LIGAPAK™ Dispensing Reel

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	4-0	Black	1 x 144"/365cm	12		LA53G
	3-0	Black	1 x 144"/365cm	12		LA54G
	2-0	Black	1 x 144"/365cm	12		LA55G

SUTUPAK® Pre-Cut Sterile Sutures


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Black	2 x 60"/150cm	12		SA11G
	4-0	Black	17 x 18"/45cm	36		SA63H
	4-0	Black	3 x 24"/60cm	36		SA73H
	4-0	Black	10 x 30"/75cm	36		SA83H
	4-0	Black	2 x 60"/150cm	12		SA10G
	3-0	Black	17 x 18"/45cm	36		SA64H
	3-0	Black	3 x 24"/60cm	36		SA74H
	3-0	Black	10 x 30"/75cm	36		SA84H
	3-0	Black	2 x 60"/150cm	12		SA9G
	2-0	Black	17 x 18"/45cm	36		SA65H
	2-0	Black	3 x 24"/60cm	36		SA75H
	2-0	Black	10 x 30"/75cm	36		SA85H
	2-0	Black	2 x 60"/150cm	36		SA5H
	2-0	Black	2 x 60"/150cm	36		SA8H
	1	Black	13 x 24"/60cm	12		SA77G
	1	Black	10 x 30"/75cm	12		SA87G
	1	Black	2 x 60"/150cm	36		SA7H
	0	Black	17 x 18"/45cm	12		SA66G
	0	Black	3 x 24"/60cm	12		SA76G
	0	Black	10 x 30"/75cm	12		SA86G
	0	Black	2 x 60"/150cm	36		SA6H

PERMA-HAND[®] Silk Suture



SUTUPAK Pre-Cut Sutures in Labyrinth Package


No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Black	12 x 18"/45cm	36		A182H
	5-0	Black	12 x 30"/75cm	36		A302H
	4-0	Black	12 x 18"/45cm	36		A183H
	4-0	Black	12 x 30"/75cm	36		A303H
	3-0	Black	12 x 18"/45cm	36		A184H
	3-0	Black	12 x 30"/75cm	36		A304H
	2-0	Black	12 x 18"/45cm	36		A185H
	2-0	Black	12 x 30"/75cm	36		A305H
	1	Black	6 x 18"/45cm	36		A187H
	1	Black	6 x 30"/75cm	36		A307H
	0	Black	6 x 18"/45cm	36		A186H
	0	Black	6 x 30"/75cm	36		A306H




PROLENE® Polypropylene Suture



Widely used in general, cardiovascular, plastic, and orthopedic surgery, PROLENE Sutures do not adhere to tissue and are therefore efficacious as a pull-out suture. PROLENE Sutures are relatively biologically inert, offering proven strength, reliability, and versatility. PROLENE Sutures are recommended for use where minimal suture reaction is desired, such as in contaminated and infected wounds to minimize later sinus formation and suture extrusion.¹

















PROLENE suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.¹




PROLENE® Polypropylene Suture







































ETHIGUARD® Blunt Point Needle						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTB-1 36mm 1/2 Circle	0	Blue	1 x 30"/75cm	36		PB21H
	1	Blue	1 x 30"/75cm	36		PB20H
 CTXB 48mm 1/2 Circle	1	Blue	1 x 30"/75cm	36		PB30H


Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV75-3 3.8mm 3/8 Circle	10-0	Blue	1 x 5"/13cm	12		2790G
 BV100-4 5.1mm 3/8 Circle	10-0	Blue	1 x 5"/13cm	12		2794G
	9-0	Blue	1 x 5"/13cm	12		2793G
 BV130-4 5.1mm 3/8 Circle	10-0	Blue	1 x 2"/5cm	12		2795G
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	1 x 5"/13cm	12		2775G
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	1 x 18"/45cm	36		8730H
	8-0	Blue	1 x 18"/45cm	36	 	EP8730H
	8-0	Blue	2 x 18"/45cm	12		M8738
	8-0	Blue	2 x 18"/45cm	12	 	EPM8738
	8-0	Blue	4 x 18"/45cm	12		M8730
	8-0	Blue	4 x 18"/45cm	12	 	EPM8730









































PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV130-5 6.5mm 3/8 Circle	8-0	Blue	1 x 24"/60cm	36	 MultiPass® NEEDLES	8732H
	8-0	Blue	1 x 24"/60cm	36	 EVP  MultiPass® NEEDLES	EP8732H
	8-0	Blue	2 x 24"/60cm	12	 MultiPass® NEEDLES	M8739
	8-0	Blue	2 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8739
	8-0	Blue	4 x 24"/60cm	12	 MultiPass® NEEDLES	M8732
	8-0	Blue	4 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8732
 BV175-6 8mm 3/8 Circle	8-0	Blue	1 x 18"/45cm	36	 MultiPass® NEEDLES	8740H
	8-0	Blue	1 x 18"/45cm	36	 EVP  MultiPass® NEEDLES	EP8740H
	8-0	Blue	2 x 18"/45cm	12	 MultiPass® NEEDLES	M8742
	8-0	Blue	2 x 18"/45cm	12	 EVP  MultiPass® NEEDLES	EPM8742
	8-0	Blue	1 x 24"/60cm	36	 MultiPass® NEEDLES	8741H
	8-0	Blue	1 x 24"/60cm	36	 EVP  MultiPass® NEEDLES	EP8741H
	8-0	Blue	2 x 24"/60cm	12	 MultiPass® NEEDLES	M8743
	8-0	Blue	2 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8743
	8-0	Blue	4 x 24"/60cm	12	 MultiPass® NEEDLES	M8741
	8-0	Blue	4 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8741
	7-0	Blue	1 x 18"/45cm	36	 MultiPass® NEEDLES	8734H
	7-0	Blue	1 x 18"/45cm	36	 EVP  MultiPass® NEEDLES	EP8734H
	7-0	Blue	2 x 18"/45cm	12	 MultiPass® NEEDLES	M8736
	7-0	Blue	2 x 18"/45cm	12	 EVP  MultiPass® NEEDLES	EPM8736
	7-0	Blue	4 x 18"/45cm	12	 MultiPass® NEEDLES	M8734
	7-0	Blue	4 x 18"/45cm	12	 EVP  MultiPass® NEEDLES	EPM8734






























PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV175-6 8mm 3/8 Circle	7-0	Blue	1 x 24"/60cm	36		8735H
	7-0	Blue	1 x 24"/60cm	36	 	EP8735H
	7-0	Blue	2 x 24"/60cm	12		M8737
	7-0	Blue	2 x 24"/60cm	12	 	EPM8737
	7-0	Blue	4 x 24"/60cm	12		M8735
	7-0	Blue	4 x 24"/60cm	12	 	EPM8735
	7-0	Blue	2 x 30"/75cm	12		M8733
	7-0	Blue	2 x 30"/75cm	12	 	EPM8733
	7-0	Blue	4 x 30"/75cm	12		M8731
	7-0	Blue	4 x 30"/75cm	12	 	EPM8731
 BV175-7 8mm 3/8 Circle	8-0	Blue	2 x 24"/60cm	12		M8767
	8-0	Blue	2 x 24"/60cm	12	 	EPM8767
	7-0	Blue	2 x 18"/45cm	12		M8765
	7-0	Blue	2 x 18"/45cm	12	 	EPM8765
	7-0	Blue	4 x 18"/45cm	12		M8764
	7-0	Blue	1 x 24"/60cm	36		8766H
	7-0	Blue	1 x 24"/60cm	36	 	EP8766H
	7-0	Blue	2 x 24"/60cm	12		M8766
	7-0	Blue	2 x 24"/60cm	12	 	EPM8766
	7-0	Blue	4 x 24"/60cm	12		M8775
	7-0	Blue	4 x 24"/60cm	12	 	EPM8775
	7-0	Blue	2 x 30"/75cm	12		M8768
	7-0	Blue	2 x 30"/75cm	12	 	EPM8768
	7-0	Blue	4 x 30"/75cm	12		M8785
	7-0	Blue	4 x 30"/75cm	12	 	EPM8785

















PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV175-8 9.3mm 3/8 Circle	8-0	Blue	1 x 18"/45cm	36		8752H
	8-0	Blue	1 x 18"/45cm	36	 	EP8752H
	8-0	Blue	2 x 18"/45cm	12		M8752
	8-0	Blue	2 x 18"/45cm	12	 	EPM8752
	8-0	Blue	1 x 24"/60cm	36		8753H
	8-0	Blue	1 x 24"/60cm	36	 	EP8753H
	8-0	Blue	2 x 24"/60cm	12		M8753
	8-0	Blue	2 x 24"/60cm	12	 	EPM8753
	8-0	Blue	1 x 30"/75cm	36		8754H
	8-0	Blue	1 x 30"/75cm	36	 	EP8754H
	8-0	Blue	2 x 30"/75cm	12		M8754
	8-0	Blue	2 x 30"/75cm	12	 	EPM8754
	7-0	Blue	1 x 18"/45cm	36		8746H
	7-0	Blue	1 x 18"/45cm	36	 	EP8746H
	7-0	Blue	2 x 18"/45cm	12		M8746
	7-0	Blue	2 x 18"/45cm	12	 	EPM8746
	7-0	Blue	1 x 24"/60cm	36		8747H
	7-0	Blue	1 x 24"/60cm	36	 	EP8747H



PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV175-8 9.3mm 3/8 Circle	7-0	Blue	2 x 24"/60cm	12		M8747
	7-0	Blue	2 x 24"/60cm	12	 	EPM8747
	7-0	Blue	2 x 30"/75cm	12		M8755
	7-0	Blue	2 x 30"/75cm	12	 	EPM8755
	7-0	Blue	4 x 24"/60cm	12		M8745
	7-0	Blue	4 x 24"/60cm	12	 	EPM8745
	7-0	Blue	4 x 30"/75cm	12	 	EPM8751
	7-0	Blue	1 x 30"/75cm	36		8755H
	7-0	Blue	1 x 30"/75cm	36	 	EP8755H





PROLENE® Polypropylene Suture




Taperpoint 							
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code	
 BV-1 9.3mm 3/8 Circle	7-0	Blue	1 x 18"/45cm	36		8701H	
	7-0	Blue	1 x 18"/45cm	36	VB	8301H	
	7-0	Blue	4 x 18"/45cm	12		M8701	
	7-0	Blue	4 x 18"/45cm	12	VB	M8301	
	7-0	Blue	1 x 24"/60cm	36		8702H	
	7-0	Blue	1 x 24"/60cm	36	VB	8304H	
	7-0	Blue	4 x 24"/60cm	12		M8702	
	7-0	Blue	4 x 24"/60cm	12	VB	M8304	
	7-0	Blue	4 x 24"/60cm	12	EVP	EPM8702	
	7-0	Blue	1 x 30"/75cm	36		8703H	
	7-0	Blue	4 x 18"/45cm	12		EPM8701	
	7-0	Blue	2 x 30"/75cm	12		M8203	
	7-0	Blue	2 x 30"/75cm	12	EVP	EPM8203	
	7-0	Blue	4 x 30"/75cm	12		M8703	
	7-0	Blue	4 x 30"/75cm	12	VB	M8303	
	7-0	Blue	4 x 30"/75cm	12	EVP	EPM8703	
	6-0	Blue	1 x 18"/45cm	36		8806H	
	6-0	Blue	1 x 18"/45cm	36	EVP	EP8806H	
	6-0	Blue	4 x 18"/45cm	12		M8806	
	6-0	Blue	4 x 18"/45cm	12	EVP	EPM8806	
	6-0	Blue	1 x 24"/60cm	36		8805H	
	6-0	Blue	1 x 24"/60cm	36	EVP	EP8805H	
	6-0	Blue	1 x 24"/60cm	36	VB	8305H	
	6-0	Blue	2 x 24"/60cm	12		M8205	
	6-0	Blue	2 x 24"/60cm	12	EVP	EPM8205	
	6-0	Blue	4 x 24"/60cm	12		M8805	
	6-0	Blue	4 x 24"/60cm	12	VB	M8305	
	6-0	Blue	4 x 24"/60cm	12	EVP	EPM8805	
	6-0	Blue	1 x 30"/75cm	36		8709H	
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8709H	
	6-0	Blue	4 x 30"/75cm	12		M8709	
	6-0	Blue	4 x 30"/75cm	12	EVP	EPM8709	
	5-0	Blue	1 x 24"/60cm	36		HS	9702H
	5-0	Blue	1 x 24"/60cm	36		HS EVP	EPH9702H


PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 BV 11mm 3/8 Circle	7-0	Blue	4 x 24"/60cm	12		M8600
	6-0	Blue	1 x 24"/60cm	36		8610H
	6-0	Blue	1 x 24"/60cm	36	EVP	EP8610H
	6-0	Blue	4 x 24"/60cm	12		M8610
	6-0	Blue	4 x 24"/60cm	12	EVP	EPM8610
	6-0	Blue	1 x 30"/75cm	36		8776H
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8776H
	6-0	Blue	4 x 30"/75cm	12		M8776
	6-0	Blue	4 x 30"/75cm	12	EVP	EPM8776
 C-1 13mm 3/8 Circle	6-0	Blue	1 x 30"/75cm	36		8889H
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8889H
	5-0	Blue	1 x 30"/75cm	36		8890H
	5-0	Blue	1 x 30"/75cm	36	HS EVP	EPH8890H
 C-1 13mm 3/8 Circle	7-0	Blue	1 x 18"/45cm	36		8800H
	7-0	Blue	1 x 18"/45cm	36	EVP	EP8800H
	7-0	Blue	4 x 18"/45cm	12		M8800
	7-0	Blue	4 x 18"/45cm	12	EVP	EPM8800
	7-0	Blue	4 x 24"/60cm	12		M8727
	7-0	Blue	4 x 24"/60cm	12	EVP	EPM8727
	7-0	Blue	1 x 30"/75cm	36		8708H
	7-0	Blue	1 x 30"/75cm	36	EVP	EP8708H







PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13mm 3/8 Circle	6-0	Blue	1 x 18"/45cm	36		8718H
	6-0	Blue	1 x 18"/45cm	36	VB	8318H
	6-0	Blue	1 x 18"/45cm	36	EVP	EP8718H
	6-0	Blue	4 x 18"/45cm	12		M8718
	6-0	Blue	4 x 18"/45cm	12	VB	M8318
	6-0	Blue	4 x 18"/45cm	12	EVP	EPM8718
	6-0	Blue	1 x 24"/60cm	36		8726H
	6-0	Blue	1 x 24"/60cm	36	EVP	EP8726H
	6-0	Blue	4 x 24"/60cm	12		M8726
	6-0	Blue	4 x 24"/60cm	12	EVP	EPM8726
	6-0	Blue	1 x 30"/75cm	36		8706H
	6-0	Blue	1 x 30"/75cm	36	VB	8307H
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8706H
	6-0	Blue	2 x 30"/75cm	12		M8206
	6-0	Blue	2 x 30"/75cm	12	EVP	EPM8206
	6-0	Blue	4 x 30"/75cm	12		M8706
	6-0	Blue	4 x 30"/75cm	12	VB	M8307
	6-0	Blue	4 x 30"/75cm	12	EVP	EPM8706
	5-0	Blue	1 x 18"/45cm	36		8717H
	5-0	Blue	4 x 18"/45cm	12		M8717
	5-0	Blue	1 x 24"/60cm	36		8725H
	5-0	Blue	1 x 24"/60cm	36	VB	8325H
	5-0	Blue	1 x 24"/60cm	36	HS EVP	EPH8725H
	5-0	Blue	4 x 24"/60cm	12		M8725
	5-0	Blue	4 x 24"/60cm	12	VB	M8325
	5-0	Blue	1 x 30"/75cm	36	HS EVP	EPH7477H
	5-0	Blue	1 x 36"/90cm	36		8720H
	5-0	Blue	1 x 36"/90cm	36	VB	8321H
	5-0	Blue	1 x 36"/90cm	36	HS EVP	EPH8720H
	5-0	Blue	2 x 36"/90cm	12		M8220
	5-0	Blue	4 x 36"/90cm	12		M8720
	5-0	Blue	4 x 36"/90cm	12	VB	M8321







PROLENE® Polypropylene Suture




Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 C-1 13mm 3/8 Circle	4-0	Blue	1 x 30"/75cm	36	HS	9706H
	4-0	Blue	1 x 36"/90cm	36		8519H
 BB 17mm 3/8 Circle	4-0	Blue	1 x 30"/75cm	36		8881H
 BB 17mm 3/8 Circle	5-0	Blue	1 x 24"/60cm	36		8860H
	5-0	Blue	1 x 36"/90cm	36		8580H
	4-0	Blue	1 x 24"/60cm	36		8861H
	4-0	Blue	1 x 36"/90cm	36		8581H
	4-0	Blue	1 x 36"/90cm	36	VB	8381H
	3-0	Blue	1 x 24"/60cm	36		8862H
 RB-3 11mm 1/2 Circle	7-0	Blue	1 x 24"/60cm	36		8206H
 RB-2 13mm 1/2 Circle	6-0	Blue	1 x 18"/45cm	36		8714H
	6-0	Blue	1 x 18"/45cm	36	EVP	EP8714H
	6-0	Blue	1 x 30"/75cm	36		8711H
	6-0	Blue	1 x 30"/75cm	36	EVP	EP8711H
	6-0	Blue	4 x 30"/75cm	12		M8711
	6-0	Blue	4 x 30"/75cm	12	EVP	EPM8711
	5-0	Blue	1 x 18"/45cm	36		8713H
	5-0	Blue	1 x 24"/60cm	36		8554H
	5-0	Blue	1 x 30"/75cm	36		8710H
	5-0	Blue	1 x 30"/75cm	36	HS EVP	EPH8710H
	5-0	Blue	4 x 30"/75cm	12		M8710
	5-0	Blue	1 x 36"/90cm	36	HS EVP	EPH8716H







PROLENE® Polypropylene Suture



Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TF 13mm 1/2 Circle	5-0	Blue	1 x 24"/60cm	36		8205H
	4-0	Blue	1 x 24"/60cm	36		8204H
 RB-1 17mm 1/2 Circle	5-0	Blue	1 x 30"/75cm	36		8870H
	4-0	Blue	1 x 30"/75cm	36		8871H
	3-0	Blue	1 x 30"/75cm	36		8872H
 RB-1 17mm 1/2 Circle	5-0	Blue	1 x 18"/45cm	36		8756H
	5-0	Blue	1 x 24"/60cm	36		8555H
	5-0	Blue	1 x 30"/75cm	36	HS	HS6856H
	5-0	Blue	1 x 36"/90cm	36		8556H
	5-0	Blue	1 x 36"/90cm	36	VB	8356H
	5-0	Blue	1 x 36"/90cm	36	HS	HS8556H
	5-0	Blue	4 x 36"/90cm	12		M8556
	4-0	Blue	1 x 18"/45cm	36		8757H
	4-0	Blue	1 x 36"/90cm	36		8557H
	4-0	Blue	1 x 36"/90cm	36	VB	8357H
	4-0	Blue	1 x 36"/90cm	36	HS	HS6857H
	4-0	Blue	4 x 36"/90cm	12		M8557
	3-0	Blue	1 x 36"/90cm	36		8558H
	3-0	Blue	1 x 36"/90cm	36	HS	9556H
	2-0	Blue	1 x 36"/90cm	36		8559H
 SH-2 20mm 1/2 Circle	4-0	Blue	1 x 36"/90cm	36	HS	HS6861H
 MO-7 22mm 1/2 Circle	0	Blue	8 x 18"/45cm	12	CR	C841G


PROLENE® Polypropylene Suture









Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SH-1 22mm 1/2 Circle	4-0	Blue	1 x 30"/75cm	36		8526H
	3-0	Blue	1 x 30"/75cm	36		8762H
 CT-2 26mm 1/2 Circle	3-0	Blue	1 x 30"/75cm	36		8410H
	2-0	Blue	8 x 18"/45cm	12	CR	C826G
	2-0	Blue	1 x 30"/75cm	36		8411H
	0	Blue	8 x 18"/45cm	12	CR	C827G
	0	Blue	1 x 30"/75cm	36		8412H
 MO-6 26mm 1/2 Circle	2-0	Blue	1 x 30"/75cm	36		8417H
	0	Blue	8 x 18"/45cm	12	CR	C845G
	0	Blue	1 x 30"/75cm	36		8418H
 SH 26mm 1/2 Circle	4-0	Blue	1 x 30"/75cm	36		8831H
	3-0	Blue	1 x 30"/75cm	36		8832H
	2-0	Blue	1 x 30"/75cm	36		8833H
	0	Blue	1 x 30"/75cm	36		8834H
 SH 26mm 1/2 Circle	5-0	Blue	1 x 36"/90cm	36	HS	HS6822H
	4-0	Blue	1 x 36"/90cm	36		8521H
	4-0	Blue	1 x 36"/90cm	36	HS	HS6855H
	3-0	Blue	1 x 36"/90cm	36		8522H
	3-0	Blue	1 x 36"/90cm	36	Pledget*	P8522
	3-0	Blue	1 x 48"/122cm	36		8534H
	2-0	Blue	1 x 36"/90cm	36		8523H
	2-0	Blue	1 x 48"/122cm	36		8533H


*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget




PROLENE® Polypropylene Suture










Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CT-1 36mm 1/2 Circle	3-0	Blue	1 x 30"/75cm	36		8422H
	2-0	Blue	8 x 18"/45cm	12	CR	C822G
	2-0	Blue	1 x 30"/75cm	36		8423H
	0	Blue	8 x 18"/45cm	12	CR	C821G
	0	Blue	1 x 30"/75cm	36		8424H
	1	Blue	1 x 30"/75cm	36		8425H
 MH 36mm 1/2 Circle	3-0	Blue	1 x 36"/90cm	36		8842H
	2-0	Blue	1 x 36"/90cm	36		8843H
	2-0	Blue	1 x 48"/122cm	36		8853H
 CT 40mm 1/2 Circle	2-0	Blue	1 x 30"/75cm	36		8433H
	0	Blue	1 x 30"/75cm	36		8434H
	1	Blue	1 x 30"/75cm	36		8435H
 CTX 48mm 1/2 Circle	0	Blue	8 x 18"/45cm	12	CR	C811G
	0	Blue	1 x 30"/75cm	36		8454H
	1	Blue	1 x 30"/75cm	36		8455H
 TP-1 65mm 1/2 Circle	1	Blue	1 x 60"/150cm	12		8824G
	2	Blue	1 x 60"/150cm	12		8825G


PROLENE® Polypropylene Suture
















Taperpoint 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 XLH 70mm 1/2 Circle	0	Blue	3 x 24"/60cm	12		3844G
	1	Blue	1 x 60"/150cm	12		8845G
 XLH 70mm 1/2 Circle	1	Blue	1 x 30"/75cm	12		8190G


TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CIF-4 13mm 1/4 Circle	10-0	Blue	1 x 8"/20cm	12		788G
 CC 9.3mm 3/8 Circle	7-0	Blue	1 x 18"/45cm	36		8705H
	7-0	Blue	4 x 18"/45cm	12		M8705
	7-0	Blue	4 x 18"/45cm	12	 	EPM8705








PROLENE® Polypropylene Suture





TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CC 9.3mm 3/8 Circle</p>	7-0	Blue	1 x 24"/60cm	36	 MultiPass® NEEDLES	8704H
	7-0	Blue	4 x 24"/60cm	12	 MultiPass® NEEDLES	M8704
	7-0	Blue	4 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8704
	6-0	Blue	1 x 24"/60cm	36	 MultiPass® NEEDLES	8807H
	6-0	Blue	1 x 24"/60cm	36	 EVP  MultiPass® NEEDLES	EP8807H
	6-0	Blue	4 x 24"/60cm	12	 MultiPass® NEEDLES	M8807
	6-0	Blue	4 x 24"/60cm	12	 EVP  MultiPass® NEEDLES	EPM8807
 <p>CC-1 13mm 3/8 Circle</p>	6-0	Blue	1 x 18"/45cm	36	 MultiPass® NEEDLES	8719H
	6-0	Blue	1 x 18"/45cm	36	 EVP  MultiPass® NEEDLES	EP8719H
	6-0	Blue	4 x 18"/45cm	12	 MultiPass® NEEDLES	M8719
	6-0	Blue	4 x 18"/45cm	12	 EVP  MultiPass® NEEDLES	EPM8719
	6-0	Blue	1 x 30"/75cm	36	 MultiPass® NEEDLES	8707H
	6-0	Blue	1 x 30"/75cm	36	 EVP  MultiPass® NEEDLES	EP8707H
	6-0	Blue	4 x 30"/75cm	12	 MultiPass® NEEDLES	M8707
	6-0	Blue	4 x 30"/75cm	12	 EVP  MultiPass® NEEDLES	EPM8707
	5-0	Blue	1 x 36"/90cm	36	 MultiPass® NEEDLES	8721H
	5-0	Blue	4 x 36"/90cm	12	 MultiPass® NEEDLES	M8721
 <p>V-4 17mm 3/8 Circle</p>	4-0	Blue	1 x 36"/90cm	36		8925H
	3-0	Blue	1 x 36"/90cm	36		8926H

PROLENE® Polypropylene Suture










TAPERCUT® 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
CC175-8 9.3mm 3/8 Circle	7-0	Blue	4 x 24"/60cm	12		EPM8704SL
	7-0	Blue	4 x 30"/75cm	12		EPM8812SL
 V-5 17mm 1/2 Circle	5-0	Blue	1 x 36"/90cm	36	Pledget*	8934H
	4-0	Blue	1 x 36"/90cm	36		8935H
	4-0	Blue	1 x 36"/90cm	36		P8935
	3-0	Blue	1 x 36"/90cm	36		8936H
	2-0	Blue	1 x 36"/90cm	36		8937H
 V-7 26mm 1/2 Circle	4-0	Blue	1 x 36"/90cm	36		8975H
	3-0	Blue	1 x 36"/90cm	36		8976H
	2-0	Blue	1 x 36"/90cm	36		8977H
	0	Blue	1 x 36"/90cm	36		8978H
 V-34 36mm 1/2 Circle	0	Blue	1 x 30"/75cm	36		8444H
 V-47 36mm Straight	4-0	Blue	1 x 14"/35cm	12		8592G



Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-2 19mm 3/8 Circle	5-0	Blue	1 x 18"/45cm	12		8661G
	4-0	Blue	1 x 18"/45cm	12		8683G
	3-0	Blue	1 x 18"/45cm	12		8665G

*7mm x 3mm x 1.5mm Firm TFE Polymer Pledget


PROLENE® Polypropylene Suture

























Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 FS-1 24mm 3/8 Circle	4-0	Blue	1 x 18"/45cm	36		8629H
	3-0	Blue	1 x 18"/45cm	12		8684G
 FS 26mm 3/8 Circle	2-0	Blue	1 x 18"/45cm	36		8685H
 FSL 30mm 3/8 Circle	3-0	Blue	1 x 30"/75cm	36		8675H
 FSLX 36mm 3/8 Circle	3-0	Blue	1 x 30"/75cm	36		8649H
	2-0	Blue	1 x 30"/75cm	36		8689H
	0	Blue	1 x 30"/75cm	36		8690H
 LR 75mm 3/8 Circle	2	Blue	3 x 20"/50cm	24		3846T
 STC-6 16mm Straight	10-0	Blue	1 x 8"/20cm	12		1713G


Straight Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 TS 48mm Straight	2-0	Blue	1 x 14"/35cm	12		8571G
 KS 60mm Straight	4-0	Blue	1 x 30"/75cm	36		8621H
	3-0	Blue	1 x 30"/75cm	36		8622H
	2-0	Blue	1 x 30"/75cm	36		8623H


PROLENE® Polypropylene Suture





PRIME Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 P-6 8mm 1/4 Circle	7-0	Blue	1 x 18"/45cm	12		8648G
 P-1 11mm 3/8 Circle	7-0	Blue	1 x 18"/45cm	12		8696G
	6-0	Undyed	1 x 18"/45cm	12		8606G
	6-0	Blue	1 x 18"/45cm	12		8697G
 P-3 13mm 3/8 Circle	6-0	Blue	1 x 18"/45cm	12		8695G
	5-0	Undyed	1 x 18"/45cm	12		8605G
	5-0	Blue	1 x 18"/45cm	12		8698G
	4-0	Undyed	1 x 18"/45cm	12		8604G
	4-0	Blue	1 x 18"/45cm	36		8699G
 PS-3 16mm 3/8 Circle	6-0	Blue	1 x 18"/45cm	12		8680G
	5-0	Blue	1 x 18"/45cm	12		8681G
 PS-2 19mm 3/8 Circle	5-0	Blue	1 x 18"/45cm	12		8686G
	4-0	Blue	1 x 18"/45cm	12		8682G
	4-0	Blue	1 x 18"/45cm	36		8682H
	3-0	Blue	1 x 18"/45cm	36		8687H
 PS-1 24mm 3/8 Circle	3-0	Blue	1 x 18"/45cm	12		8663G
 PS-5 13mm 1/2 Circle	5-0	Blue	1 x 18"/45cm	12		8655G
	4-0	Blue	1 x 18"/45cm	12		8656G
 PS-4 16mm 1/2 Circle	4-0	Undyed	1 x 18"/45cm	12		8603G

PROLENE® Polypropylene Suture










PRIME Conventional Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 PC-1 13mm 3/8 Circle	6-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8617G
	5-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8618G
	4-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8619G
 PC-3 16mm 3/8 Circle	6-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8636G
	5-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8635G
	4-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8634G
 PC-5 19mm 3/8 Circle	5-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8630G
	4-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8631G
	3-0	Blue	1 x 18"/45cm	12	 MultiPass® NEEDLES	8632G

Spatula Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 SM-1 6mm 1/8 Circle	5-0	Blue	1 x 12"/30cm	12		7740G
 TG140-8 6.5mm 3/8 Circle	10-0	Blue	1 x 6"/15cm	12		1757G
	10-0	Blue	1 x 12"/30cm	12		1771G
	9-0	Blue	1 x 6"/15cm	12		1754G

ULTIMA® Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTC-6L 13mm 1/4 Circle	10-0	Blue	1 x 4"/10cm	12		9091G

PROLENE® Polypropylene Suture



ULTIMA® Reverse Cutting 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CTC-6 13mm/ CS160-6 5.5mm 1/4 Circle/1/2 Circle	10-0	Blue	1 x 4"/10cm	12		9090G
 CS140-6 6.5mm 3/8 Circle	10-0	Blue	1 x 12"/30cm	12		1710G
 CS160-6 5.5mm 1/2 Circle	10-0	Blue	1 x 12"/30cm	12		7761G
 CS160-8 5.5mm 1/2 Circle	10-0	Blue	1 x 6"/15cm	12		1755G
 CSC-6 5mm Compound Curve	10-0	Blue	1 x 12"/30cm	12		7796G


Surgical Stainless Steel Suture

Surgical stainless steel suture is a nonabsorbable, surgical suture composed of 316L stainless steel. The suture is available as a monofilament and multifilament suture, and it comes in sterile and non-sterile strands.¹

Surgical stainless steel suture is indicated for use in abdominal wound closure, hernia repair, sternal closure and orthopaedic procedures including cerclage and tendon repair.¹






Surgical Stainless Steel Suture


Taperpoint 




Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>CT-1 36mm 1/2 Circle</p>	2-0	Silver	8 x 18"/45cm	12		M682G
 <p>CTX 48mm 1/2 Circle</p>	2-0	Silver	8 x 24"/60cm	12		M414G

TAPERCUT® 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>V-40 48mm 1/2 Circle</p>	4	Silver	4 x 18"/45cm	12		M651G
	5	Silver	4 x 18"/45cm	12		M650G
	5	Silver	2 x 18"/45cm	12		M656G
	6	Silver	4 x 18"/45cm	12		M649G
 <p>V-56 60mm 1/2 Circle</p>	5	Silver	4 x 30"/75cm	12		M435G
 <p>V-37 40mm 1/2 Circle</p>	1	Silver	4 x 18"/45cm	12		M660G

Surgical Stainless Steel Suture








Reverse Cutting ▾						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>FS-2 19mm 3/8 Circle</p>	4-0	Silver	1 x 18"/45cm	12		603G
	3-0	Silver	1 x 18"/45cm	12		607G
 <p>FS-1 24mm 3/8 Circle</p>	3-0	Silver	1 x 18"/45cm	12		631G
 <p>FSLX 36mm 3/8 Circle</p>	2-0	Silver	1 x 18"/45cm	12		618G
 <p>LS-1 45mm 3/8 Circle</p>	2	Silver	4 x 20"/50cm	12		M400G
 <p>LR 75mm 3/8 Circle</p>	2	Silver	1 x 30"/75cm	12		467G
 <p>CP-1 36mm 1/2 Circle</p>	3-0	Silver	1 x 20"/50cm	12		497G

Surgical Stainless Steel Suture




Cutting Edge Conventional 						
Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 CCS 48mm 1/2 Circle	4	Silver	4 x 18"/45cm	12		M652G
	5	Silver	2 x 18"/45cm	12		M657G
	5	Silver	4 x 18"/45cm	12		M653G
	6	Silver	4 x 18"/45cm	12		M654G
	7	Silver	4 x 18"/45cm	12		M655G
 CCS-1 60mm 1/2 Circle	5	Silver	2 x 30"/75cm	12		M436G



Taperpoint  , Straight Cutting Temporary Pacing Wire 						
Needle	Suture USP	Color	Length	EA/BX	Product Code	
 BB 17mm/ SKS 60mm 3/8 Circle/Straight	2-0	Dark Blue	1 x 24"/60cm	12	TPW20*	
 BB-1 22mm/ SKS-3 90mm 3/8 Circle/Straight	2-0	Dark Blue	1 x 24"/60cm	12	TPW40*†	
	2-0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW42*†‡	
 RB-1 17mm/ SKS 60mm 1/2 Circle/Straight	2-0	Dark Blue	1 x 24"/60cm	12	TPW10*	

*One dozen packets per box †1.5" long score marks; others .875" ‡Multistrand 2"- 60 cm


Surgical Stainless Steel Suture



Taperpoint , Straight Cutting Temporary Pacing Wire 

Needle	Suture USP	Color	Length	EA/BX	Product Code
 <p>SH 26mm/ SKS-3 90mm 1/2 Circle/Straight</p>	2-0	Dark Blue	1 x 24"/60cm	12	TPW30*†
	2-0	Dark Blue	1 x 24"/60cm	12	TPW90*†§
	2-0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW32*†‡
	2-0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW72*†
	2-0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW92*†§
	0	Dark Blue	1 x 24"/60cm	12	TPW50*
	0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW52*†
 <p>TE 32mm/ SKS-3 90mm 3/8 Circle/Straight</p>	2-0	Dark Blue	1 x 24"/60cm	12	TPW80*†§
	0	Light Blue/ Dark Blue	2 x 24"/60cm	12	TPW62*†

Straight Cutting 


Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
 <p>KS 60mm Straight</p>	5-0	Silver	8 x 16"/40cm	12		M645G

SUTUPAK® Pre-Cut Sterile Sutures

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	3-0	Silver	17 x 18"/45cm	12		SS30G
	2-0	Silver	17 x 18"/45cm	12		SS28G
	0	Silver	17 x 18"/45cm	12		SS26G

*One dozen packets per box †1.5" long score marks; others .875" ‡Multistrand 2" - 60 cm §TP wave codes

Surgical Stainless Steel Suture



ETHI-PACK™ Pre-cut Sutures (Non-sterile)

No Needle	Suture USP	Color	Length	EA/BX	Special	Product Code
	5-0	Silver	1 x 18"/45cm	12		DS35
	4-0	Silver	1 x 18"/45cm	12		DS32
	3-0	Silver	1 x 18"/45cm	12		DS30
	2-0	Silver	1 x 18"/45cm	12		DS28
	0	Silver	1 x 18"/45cm	12		DS26
	1	Silver	1 x 18"/45cm	12		DS25
	2	Silver	12 x 18"/45cm	12		DS24
	3	Silver	12 x 18"/45cm	12		DS23
	4	Silver	12 x 18"/45cm	12		DS22
	5	Silver	12 x 18"/45cm	12		DS20
	7	Silver	12 x 18"/45cm	12		DS18


Drains

BLAKE® Silicone Drains


BLAKE® Cardio Connectors - Sterile



Product Code	Description
BCC1	BLAKE Cardio Connector 1:1



BCC2	BLAKE Cardio Connector 2:1
------	----------------------------



BCC3	BLAKE Cardio Connector 3:1
------	----------------------------

BLAKE® Silicone Drains (Flat, Full, or 3/4-Fluted) - Sterile



Product Code	Description
2210	7mm Flat, 3/4 Fluted
2211	7mm Flat, Full Fluted
2212	7mm Flat, Full Fluted, 3/16" Trocar
2213	10mm Flat, 3/4 Fluted
2214	10mm Flat, Full Fluted
2215	10mm Flat, Full Fluted, 3/16" Trocar
2216	7mm Flat, 3/4 Fluted, 3/16" Trocar
2217	10mm Flat, Full Fluted, 3/16" Trocar

BLAKE® Silicone Drains (Round Hubless) - Sterile



Product Code	Description
2226	Hubless, 10 FR Round
2227	Hubless, 10 FR Round, 1/8" Trocar
2228	Hubless, 15 FR Round
2229	Hubless, 15 FR Round, 3/16" Trocar
2230	Hubless, 19 FR Round
2231	Hubless, 19 FR Round, 1/4" Trocar
2232	Hubless, 19 FR Round, 1/4" Bendable Trocar
2233	Hubless, 15 FR Round, 3/16" Bendable Trocar
2234	Hubless, 24 FR Round

J-VAC® Reservoirs

J-VAC Reservoirs



Product Code	Description
2161	150ml Closed Wound Drainage Reservoir
2162	450ml Closed Wound Drainage Reservoir
2163	300ml Closed Wound Drainage Reservoir

J-VAC Bulb Suction Reservoir - Sterile



Product Code	Description
2160	100 cc Bulb Suction Reservoir

J-VAC® Drainage Adapters




Product Code	Description
2199	J-VAC Drainage Adapter 1/4"
2209	JJ-VAC Drainage Adapter 1/8"
2298	J-VAC Drainage Adapter 3/16"

Primary Suture
Storage Carts,
Carrying Racks,
and Other Wound
Closure Devices


Other Wound Closure Devices

ENDOSUTURE ENDO-HOLDER™




Product Code	Description	EA/BX
EC11	EN-3 Needle - O, 42" ETHIBOND Suture	12
EC12	EN-3 Needle - 2-O, 42" ETHIBOND Suture	12
ESS15	ENDO-HOLDER Cannula (for use with 5mm diameter trocar)	1

ENDOKNOT® Pre-Tied Suture



Product Code	Description	EA/BX
EX10G	ETHIBOND® Extra Polyester Suture, O, 42"	12
JK10G	Coated VICRYL® (polyglactin 910) Suture, ST-3 needle, O, 42"	12
ZK10G	PDS® II (polydioxanone) Suture, ST-3 needle, O, 42"	12

ENDOLOOP® Ligature made from PDS® II (polydioxanone) Suture Material



Product Code	Description	EA/BX
EZ10G	ENDOLOOP Ligature PDS II Suture, O, 18"	12

ENDOLOOP® Ligature made from VICRYL® (polyglactin 910) Suture Material



Product Code	Description	EA/BX
EJ10G	ENDOLOOP Ligature VICRYL Suture, O, 18"	12

ENDOPATH® Needle Holder



Product Code	Description	EA/BX
E705R	ENDOPATH Needle Holder, reusable, 5mm, 30cm length	1
E705RL	ENDOPATH Needle Holder, reusable, 5mm, 45cm length	1
SRNHI	Self Righting Needle Holder, reusable, 5mm, 30cm length	1
SRNHIL	Self Righting Needle Holder, reusable, 5mm, 45cm length	1


Other Wound Closure Devices

LAPRA-TY® Suture Clip



Product Code	Description	EA/BX
KA200	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	1
XC200	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	6

Suture Assistant



Product Code	Needle Type	EA/BX
SW100	N/A	4
SW110	6" ETHIBOND EXCEL® Polyester Suture (size 0) and EEN (canoe) needle	10
SW112	6" ETHIBOND EXCEL Suture (size 2-0) and EEN (canoe) needle	10
SW120	6" ETHIBOND EXCEL Suture (size 0) and SH needle	10
SW122	6" ETHIBOND EXCEL Suture (size 2-0) and SH Needle	10
SW222	6" Coated VICRYL® (polyglactin 910) Suture (size 2-0) and SH Needle	10

PROXIMATE® Rotating Head Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PRR35	35	3.9	0.53	5.7	Regular	6
PRW35	35	3.9	0.58	6.9	White	6

PROXIMATE® Plus MD Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PMR35	35	3.9	0.53	5.7	Regular	6
PMW35	35	3.9	0.58	6.9	Wide	6

Other Wound Closure Devices

PROXIMATE® PX Fixed Head Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
PXR35	35	3.9	0.53	5.7	Regular	6
PXW35	35	3.9	0.58	6.9	Wide	6

Staple Extractor



Product Code	Number of Staples	EA/BX
PSX	Squeeze handle staple extractor	12

Retention Suture Bolster



Product Code	Description
450G	6 Bolsters: 2 1/2" Long 3/16" Diameter 1/32" Wall. Sterile Latex Tubing. 6 per packet. 12 packets per box.

Bone Wax



Product Code	Description
W31G	2.5 Grams, 1 Dozen per Box

TFE Polymer Pledgets



Product Code	Description
PCP20	TFE Pledget, White, Firm. 7.0mm x 3.5mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP30	TFE Pledget, White, Firm. 9.5mm x 4.8mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP40	TFE Pledget, White, Soft. 9.5mm x 4.8mm x 1.5mm. 36 envelopes/box (216 pledgets).
PCP50	TFE Pledget, White, Soft. 7.0mm x 3.5mm x 1.5mm. 36 envelopes/box (216 pledgets).


Other Wound Closure Devices

Suture Cart




Product Code	Description
CART	38"L x 18"W x 38"H. No suture racks included. Brackets to attach racks included. Some assembly required.

3-Wide Lazy Susan




Product Code	Description
LS3	15 3/8"L X 11 1/2"W X 20 3/4"H. Includes 6 modular racks. Can be expanded with additional modular racks.

4-Wide Lazy Susan




Product Code	Description
LS4	20 1/2"L X 11 1/2"W X 20 3/4"H. Includes 8 modular racks. Can be expanded with additional modular racks.

Dispenser Rack



Product Code	Description
MR20	25 3/4"L X 5 3/4"W X 27"H. Includes 20 modules, 10 back plates (each with 4 fasteners), 20 module clips, 30 divider panels.


Dispenser Rack



Product Code	Description
MR4	10 1/4"L X 5 3/4"W X 13 1/2"H. Includes 4 modules, 3 back plates (each with 4 fasteners), 4 module clips, 8 divider panels.

Other Wound Closure Devices

Portable Carrying Rack




Product Code	Description
PCR1	15 3/8"L x 6 1/8"W x 17 5/8"H. Includes 6 modular racks.

2-Sided Portable Carrying Rack



Product Code	Description
PCR2	15 3/8"L x 12 1/4"W x 16 3/4"H. Includes 12 modular racks.

I.V. Pole with Rack




Product Code	Description
SPR1	18 1/4"L x 16 1/2"W x 56 1/2"H. Includes 18 modules racks, 15 module clips, 30 divider panels. Includes 18 modules racks, 15 module clips, 30 divider panels.

Suture Chart



Absorbable Sutures: The Strength You Need for the Choices You Make



Nonabsorbable Sutures: The Strength You Need for the Choices You Make

	SUTURE BRAND	MATERIAL	SIZES	FREQUENT USES	BSR [†] PROFILE	
NATURAL	PERMA-HAND [®] Silk Suture	Braided Black	5, 2 through 7/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Gradual loss of all tensile strength over time	
	SURGICAL STAINLESS STEEL Suture	Monofilament Silver Metallic	7 through 5/0	Soft tissue approximation: n Abdominal n Sternal Closure n Orthopedic procedures including cerclage and tendon repair	Indefinite	
SYNTHETIC	NUROLON [®] Nylon Suture	Braided Black	1 through 6/0	General soft tissue approximation and/or ligation: n Neurological n Cardiovascular n Ophthalmic	Gradual loss of tensile strength over time	
	ETHILON [®] Nylon Suture	Monofilament Black/Green/Undyed (Clear)	2 through 11/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological n Skin Closure	Gradual loss of tensile strength over time	
	MERSILENE [®] Polyester Fiber Suture	Braided Green / Undyed (white) Monofilament Green	1 through 6/0 10/0 through 11/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	ETHIBOND EXCEL [®] Polyester Suture	Braided Green / Undyed (white)	5, 2 through 7/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	PROLENE [®] Polypropylene Suture	Monofilament Blue / Clear	2 through 10/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	
	PRONOVA [®] Poly (Hexafluoropropylene-VDF) Suture	Monofilament Blue	2/0 through 8/0	General soft tissue approximation and/or ligation: n Cardiovascular n Ophthalmic n Neurological	Indefinite	

† BSR=breaking strength retention (tensile strength of suture in vivo)

Contact and More Information252

Explanation of a Typical Selection and Page5-6

General Information 4-10

Selection and Use of Needleholders 7

Shapes 8

Steps to Create Performance 8

Table of Contents 3

Trademarks 10

Types 9

DERMABOND® Portfolio 34-35

DERMABOND ADVANCED® Topical Skin Adhesive 35

DNX12 35

DNX6 35

DERMABOND® Mini Topical Skin Adhesive 35

DHVM12 35

DERMABOND® PRINEO® Skin Closure System 35

CLR222US 35

CLR422US 35

CLR602US 35

Drains 223-225

BLAKE® Silicone Drains 224

2210 224

2211 224

2212 224

2213 224

2214 224

2215 224

2216 224

2217 224

2226 224

2227 224

2228 224

2229 224

2230 224

2231 224

2232 224

2233 224

2234 224

BCC1 224

BCC2 224

BCC3 224

J-VAC® Reservoirs 225

2160 225

2161 225

2162 225

2163 225

2199 225

2209 225

2298 225

Needle Index 11-28

1/2 Circle 18-24

1/4 Circle 12-13

3/8 Circle 13-17

5/8 Circle 24

Center Point Spatula 28

Compound Curve 27

Endoscopy 28

Sternum 27

Straight 25

Trocar 26

Plus Antibacterial Sutures 53-85

Coated VICRYL® Plus Antibacterial (polyglactin 910) Suture 54-70

VCP103G 69

VCP104G 70

VCP105G 70

VCP106G 70

VCP107G 70

VCP109G 69

VCP110G 70

VCP111G 70

VCP112G 70

VCP195H 65

VCP196H 65

VCP197H 65

VCP205G 69

VCP207G 69

VCP208G 69

VCP213H 57

VCP214H 57

VCP215H 57

VCP218H 57

VCP219H 57

VCP232H 58

VCP258H 59

VCP259H 59

VCP260H 60

VCP261H 60

VCP265H 65

VCP266H 65

VCP267H 65

VCP268H 65

VCP269H 58

VCP270H 58

VCP275H 61

VCP280H 61

VCP281H 61

VCP284G 69

VCP285G 69

VCP286G 69

VCP287G 69

VCP303H 57

VCP304H 57

VCP305H 57

VCP306H 57

VCP310H 57

VCP311H 57

VCP315H 58

VCP316H 58

VCP317H 58

VCP322H 60

VCP323H 60

VCP327H 57

VCP328H 57

VCP329H 57

VCP332H 58

VCP333H 58

VCP334H 58

VCP335H 58

VCP338H 59

VCP339H 59

VCP340H 60

VCP341H 60

VCP344H 59

VCP345H 59

VCP346H 60

VCP347H 60

VCP351H 61

VCP352H 61

VCP353H 61

VCP356H 61

VCP357H 61

VCP358H 61

VCP359H 61

VCP363H 61

VCP364H 62

VCP365H 62

VCP369H 61

VCP370H 62

VCP371H 62

VCP375H.....	62	VCP517H.....	63	VCP727D.....	58
VCP376H.....	62	VCP518H.....	63	VCP738D.....	59
VCP391H.....	63	VCP519H.....	63	VCP739D.....	59
VCP392H.....	63	VCP527H.....	59	VCP740D.....	60
VCP393H.....	63	VCP533H.....	65	VCP741D.....	60
VCP415H.....	58	VCP534H.....	65	VCP743D.....	56
VCP416H.....	58	VCP535H.....	65	VCP748T.....	64
VCP417H.....	58	VCP568H.....	66	VCP749T.....	64
VCP418H.....	58	VCP569H.....	66	VCP750T.....	65
VCP421H.....	63	VCP581G.....	62	VCP751D.....	61
VCP422H.....	63	VCP583G.....	62	VCP752D.....	61
VCP423H.....	63	VCP588H.....	64	VCP753D.....	61
VCP426H.....	67	VCP589H.....	64	VCP754T.....	65
VCP427H.....	67	VCP593G.....	68	VCP756T.....	66
VCP428H.....	67	VCP594G.....	68	VCP757T.....	66
VCP433H.....	56	VCP596H.....	64	VCP761D.....	65
VCP434H.....	56	VCP602H.....	62	VCP762D.....	65
VCP436H.....	60	VCP603H.....	62	VCP764D.....	62
VCP437H.....	60	VCP607H.....	70	VCP765D.....	62
VCP441H.....	64	VCP608H.....	70	VCP771D.....	57
VCP442H.....	64	VCP615H.....	70	VCP772D.....	57
VCP443H.....	64	VCP616H.....	70	VCP773D.....	58
VCP451H.....	64	VCP617H.....	70	VCP774D.....	58
VCP452H.....	64	VCP618H.....	70	VCP775D.....	58
VCP458H.....	64	VCP634T.....	69	VCP784D.....	58
VCP459H.....	64	VCP635T.....	70	VCP785D.....	58
VCP460H.....	64	VCP636T.....	70	VCP787D.....	64
VCP463G.....	67	VCP645T.....	70	VCP787Z.....	64
VCP464G.....	67	VCP646T.....	70	VCP789D.....	59
VCP466H.....	65	VCP649G.....	62	VCP790D.....	64
VCP467H.....	65	VCP656G.....	68	VCP801D.....	59
VCP468H.....	65	VCP662H.....	67	VCP822G.....	69
VCP472H.....	65	VCP663H.....	67	VCP823G.....	69
VCP474H.....	65	VCP682G.....	68	VCP823H.....	69
VCP478H.....	65	VCP682H.....	68	VCP824G.....	69
VCP480H.....	65	VCP683G.....	68	VCP834G.....	68
VCP486H.....	65	VCP683H.....	68	VCP835G.....	68
VCP490G.....	67	VCP683H.....	68	VCP838D.....	59
VCP493G.....	67	VCP694H.....	64	VCP839D.....	59
VCP493H.....	67	VCP695H.....	64	VCP840D.....	60
VCP494G.....	67	VCP698H.....	66	VCP841D.....	60
VCP494H.....	67	VCP699H.....	66	VCP844G.....	69
VCP495G.....	67	VCP700D.....	60	VCP845G.....	69
VCP495H.....	67	VCP701D.....	60	VCP849G.....	62
VCP496G.....	67	VCP702D.....	60	VCP864D.....	58
VCP496H.....	67	VCP706T.....	64	VCP868H.....	65
VCP497G.....	67	VCP708T.....	64	VCP869H.....	65
VCP497H.....	67	VCP710T.....	65	VCP870H.....	65
VCP500G.....	67	VCP711T.....	65	VCP871H.....	65
VCP500H.....	67	VCP713D.....	57	VCP880T.....	62
VCP502H.....	64	VCP714D.....	57	VCP885G.....	69
VCP503G.....	68	VCP716T.....	66	VCP916H.....	63
VCP504G.....	68	VCP717T.....	66	VCP917H.....	63
VCP506G.....	68	VCP718T.....	66	VCP931H.....	66
VCP507G.....	68	VCP719T.....	66	VCP935H.....	68
VCP511G.....	68	VCP723D.....	61	VCP936H.....	68
VCP516H.....	63	VCP724D.....	62	VCP940H.....	66
		VCP726D.....	58		

VCP944H.....	59	MCP316H.....	72	PDP149H.....	80
VCP945H.....	59	MCP317H.....	72	PDP259H.....	81
VCP946H.....	60	MCP333H.....	72	PDP303H.....	80
VCP947H.....	60	MCP339H.....	73	PDP304H.....	80
VCP956H.....	61	MCP340H.....	73	PDP305H.....	80
VCP957H.....	61	MCP344H.....	73	PDP310H.....	80
VCP958H.....	61	MCP345H.....	73	PDP311H.....	80
VCP959H.....	61	MCP346H.....	73	PDP315H.....	81
VCP977H.....	62	MCP347H.....	73	PDP316H.....	81
VCP978H.....	62	MCP358H.....	73	PDP317H.....	81
VCP979H.....	61	MCP377H.....	73	PDP320H.....	80
VCP980H.....	61	MCP415H.....	72	PDP332H.....	81
VCPB259H.....	55	MCP416H.....	72	PDP333H.....	81
VCPB260H.....	55	MCP417H.....	72	PDP334H.....	81
VCPB269H.....	55	MCP426H.....	74	PDP338H.....	81
VCPB340H.....	55	MCP427H.....	74	PDP339H.....	81
VCPB341H.....	55	MCP463G.....	74	PDP340H.....	81
VCPB416H.....	55	MCP489G.....	74	PDP341H.....	81
VCPB417H.....	55	MCP490G.....	74	PDP346H.....	81
VCPB41D.....	55	MCP492G.....	74	PDP347H.....	81
VCPB650G.....	56	MCP493G.....	74	PDP352H.....	81
VCPB724D.....	56	MCP494G.....	74	PDP353H.....	81
VCPB725D.....	56	MCP495G.....	74	PDP357H.....	81
VCPB726D.....	55	MCP496G.....	74	PDP358T.....	81
VCPB727D.....	55	MCP497G.....	74	PDP359T.....	81
VCPB739D.....	55	MCP513G.....	74	PDP370T.....	82
VCPB740D.....	55	MCP523H.....	74	PDP371T.....	82
VCPB741D.....	55	MCP604H.....	73	PDP416H.....	81
VCPB765D.....	56	MCP605H.....	73	PDP420G.....	82
VCPB774D.....	55	MCP682H.....	74	PDP421H.....	83
VCPB839D.....	55	MCP822G.....	75	PDP422H.....	83
VCPB840D.....	55	MCP823G.....	75	PDP423H.....	83
VCPB841D.....	55	MCP834G.....	75	PDP432H.....	80
VCPB864D.....	55	MCP835G.....	75	PDP441H.....	83
VCPB944H.....	55	MCP844G.....	75	PDP442H.....	83
VCPB945H.....	55	MCP845G.....	75	PDP443H.....	83
VCPB946H.....	55	MCP935H.....	74	PDP463G.....	84
VCPB947H.....	55	MCP936H.....	74	PDP466H.....	84
VCPB958H.....	56	MCP944H.....	73	PDP467H.....	84
VCPB977H.....	56	MCP945H.....	73	PDP468H.....	84
VCPB978H.....	56	MCP946H.....	73	PDP487G.....	84
VCPP31D.....	60	MCPB316H.....	72	PDP489G.....	84
VCPP40D.....	60	MCPB346H.....	72	PDP493G.....	84
VCPP41D.....	60	MCPB416H.....	72	PDP494G.....	84
VCPP42D.....	59	MCPB944H.....	72	PDP495G.....	85
VCPP71D.....	60	MCPB945H.....	72	PDP496G.....	85
VCPP80D.....	59	MCPB946H.....	72	PDP497G.....	85
VCPP81D.....	59	PDS® Plus Antibacterial		PDP503G.....	85
VCPPB31D.....	55	(polydioxanone) Suture.....76-85		PDP504G.....	85
MONOCRYL® Plus Antibacterial		PDP114H.....	80	PDP507G.....	85
(poligcaprone 25) Suture.....71-75		PDP117H.....	79	PDP513G.....	85
MCP213H.....	72	PDP123H.....	79	PDP534T.....	84
MCP214H.....	72	PDP126H.....	79	PDP535H.....	84
MCP266H.....	74	PDP127H.....	79	PDP569T.....	84
MCP304H.....	72	PDP137H.....	83	PDP582G.....	82
MCP315H.....	72	PDP148H.....	80	PDP662H.....	84
				PDP663H.....	84

PDP682G	85	ESS15	227	SXMP1B116	44
PDP683G	85	EX10G	227	SXMP1B117	44
PDP695H	83	EZ10G	227	SXMP1B118	44
PDP739D	81	JK10G	227	SXMP1B119	44
PDP762D	83	KA200	228	SXMP1B120	45
PDP776D	82	LS3	230	SXMP1B408	43
PDP779D	80	LS4	230	SXMP1B409	43
PDP789D	80	MR20	230	SXMP1B410	43
PDP800G	80	MR4	230	SXMP1B411	43
PDP801D	80	PCP20	229	SXMP1B412	43
PDP802D	80	PCP30	229	SXMP1B413	43
PDP823G	85	PCP40	229	SXMP1B414	43
PDP824G	85	PCP50	229	SXMP1B415	43
PDP833G	85	PCR1	231	SXMP1B416	43
PDP844G	85	PCR2	231	SXMP1B417	43
PDP845G	85	PMR35	228	SXMP1B419	44
PDP877G	82	PMW35	228	SXMP1B420	44
PDP879G	82	PRR35	228	SXMP1B421	44
PDP880G	82	PRW35	228	SXMP1B424	43
PDP881G	82	PSX	229	SXMP1B425	43
PDP968H	83	PXR35	229	SXMP1B427	43
PDP969H	83	PXW35	229	SXMP1B428	43
PDP970H	83	SPR1	231	SXMP1B429	43
PDP990G	82	SRNHI	227	SXMP1B433	43
PDP991G	82	SRNHIL	227	SXMP1B434	43
PDPB316H	77	SW100	228	SXMP1B435	43
PDPB317H	77	SW110	228	SXMP1B437	43
PDPB332H	77	SW112	228	STRATAFIX™ Spiral PDO	
PDPB333H	77	SW120	228	Knotless Tissue Control Device..... 38-40	
PDPB334H	77	SW122	228	SXPD2B200	40
PDPB339H	77	SW222	228	SXPD2B201	40
PDPB340	77	W31G	229	SXPD2B202	40
PDPB345	77	XC200	228	SXPD2B400	39
PDPB346	77	ZK10G	227	SXPD2B401	38
PDPB347	77	STRATAFIX™ Knotless Tissue		SXPD2B402	38
PDPB358T	78	Control Devices 36-52		SXPD2B403	38
PDPB359T	78	STRATAFIX™ Spiral MONOCRYL™ Plus		SXPD2B404	38
PDPB370	78	Knotless Tissue Control Device..... 43-45		SXPD2B405	39
PDPB371	78	SMMP1B438	43	SXPD2B406	38
PDPB879G	78	SXMP1B101	44	SXPD2B407	38
PDPB880G	77	SXMP1B102	44	SXPD2B408	38
PDPB990G	77	SXMP1B103	44	SXPD2B409	38
PDPB991G	77	SXMP1B104	44	SXPD2B410	38
		SXMP1B105	44	SXPD2B411	38
		SXMP1B106	44	SXPD2B412	38
		SXMP1B107	44	SXPD2B413	38
		SXMP1B108	44	SXPD2B414	38
		SXMP1B109	44	SXPD2B415	40
		SXMP1B110	45	SXPD2B416	40
		SXMP1B111	45	SXPD2B417	40
		SXMP1B113	45	SXPD2B418	40
		SXMP1B114	44	SXPD2B419	40
		SXMP1B115	44	SXPD2B420	40
				SXPD2B421	40
				SXPD2B422	40

Primary Suture Storage Carts, Carrying Racks, and Other Wound Closure Devices. 226-231

450G	229
CART	230
E705R	227
E705RL	227
EC1	227
EC12	227
EJ10G	227

SXPD2B423 40
 SXPD2B424 39
 SXPD2B425 39
 SXPD2B426 39

**STRATAFIX™ Spiral PDS™ Plus
 Knotless Tissue Control Device..... 45-48**

SXMP1B438 45
 SXPP1B101 47
 SXPP1B102 47
 SXPP1B103 48
 SXPP1B104 48
 SXPP1B105 48
 SXPP1B106 48
 SXPP1B107 48
 SXPP1B108 48
 SXPP1B109 48
 SXPP1B110 48
 SXPP1B111 48
 SXPP1B112 48
 SXPP1B113 48
 SXPP1B114 48
 SXPP1B116 48
 SXPP1B117 48
 SXPP1B118 48
 SXPP1B120 47
 SXPP1B201 47
 SXPP1B202 47
 SXPP1B203 47
 SXPP1B204 47
 SXPP1B205 46
 SXPP1B401 46
 SXPP1B402 46
 SXPP1B403 47
 SXPP1B404 47
 SXPP1B405 45
 SXPP1B406 46
 SXPP1B407 46
 SXPP1B408 46
 SXPP1B409 46
 SXPP1B410 46
 SXPP1B411 46
 SXPP1B412 46
 SXPP1B413 45
 SXPP1B414 45
 SXPP1B415 45
 SXPP1B416 45
 SXPP1B417 45
 SXPP1B418 47
 SXPP1B419 46
 SXPP1B420 45
 SXPP1B421 45
 SXPP1B422 45
 SXPP1B424 46
 SXPP1B425 45
 SXPP1B426 45
 SXPP1B427 45

SXPP1B428 45
 SXPP1B429 46
 SXPP1B430 46
 SXPP1B431 46
 SXPP1B432 45
 SXPP1B433 45
 SXPP1B450 46
 SXPP1B452 46
 SXPP1B453 46
 SXPP1B454 46
 SXPP1B455 46
 SXPP1B456 46
 SXPP1B457 47
 SXPP1B458 47

**STRATAFIX™ Spiral PGA-PCL
 Knotless Tissue Control Device..... 41-42**

SXMD2B150 41
 SXMD2B400 41
 SXMD2B401 41
 SXMD2B402 41
 SXMD2B403 41
 SXMD2B404 41
 SXMD2B405 41
 SXMD2B406 41
 SXMD2B407 41
 SXMD2B408 41
 SXMD2B409 41
 SXMD2B410 41
 SXMD2B411 42
 SXMD2B412 42
 SXMD2B413 42
 SXMD2B414 42

**STRATAFIX™ Spiral Polypropylene
 Knotless Tissue Control Device..... 42**

SXPL1B400 42
 SXPL2B401 42
 SXPL2B400 42

**STRATAFIX™ Symmetric PDS™ Plus
 Knotless Tissue Control Device..... 48-52**

SXPP1A100 52
 SXPP1A101 52
 SXPP1A200 51
 SXPP1A201 51
 SXPP1A202 51
 SXPP1A203 51
 SXPP1A204 51
 SXPP1A205 51
 SXPP1A300 48
 SXPP1A301 48
 SXPP1A302 48
 SXPP1A303 48
 SXPP1A304 48
 SXPP1A305 48
 SXPP1A306 48

SXPP1A400 50
 SXPP1A401 50
 SXPP1A402 50
 SXPP1A403 50
 SXPP1A404 50
 SXPP1A405 50
 SXPP1A406 50
 SXPP1A407 49
 SXPP1A408 49
 SXPP1A409 49
 SXPP1A410 49
 SXPP1A412 50
 SXPP1A413 49
 SXPP1A414 49
 SXPP1A415 49
 SXPP1A416 50
 SXPP1A417 50
 SXPP1A418 50
 SXPP1A419 49
 SXPP1A420 50
 SXPP1A421 49
 SXPP1A422 49
 SXPP1A423 50
 SXPP1A424 51
 SXPP1A425 50
 SXPP1A426 51
 SXPP1A427 49
 SXPP1A428 49
 SXPP1A429 50
 SXPP1A430 49
 SXPP1A431 49
 SXPP1A432 50
 SXPP1A433 50
 SXPP1A434 50
 SXPP1A435 50
 SXPP1A436 50
 SXPP1A437 51
 SXPP1A438 50
 SXPP1A439 50
 SXPP1A440 50
 SXPP1A441 50
 SXPP1A442 50
 SXPP1A443 50
 SXPP1A444 50
 SXPP1A445 50
 SXPP1A446 49
 SXPP1A447 50

Suture Chart 232-234

Absorbable Sutures Chart 233

Nonabsorbable Sutures Chart 234

Sutures: Absorbable 86-148**Coated VICRYL®****(polyglactin 910) Suture 91-110**

J103G	109	J322H	97	J437H	97
J103T	109	J323H	97	J441H	101
J104G	109	J327H	93	J442H	101
J104T	109	J328H	93	J443H	101
J105G	109	J329H	93	J446G	107
J105T	109	J332H	94	J451H	101
J106G	109	J333H	94	J455H	101
J106T	109	J334H	94	J458H	101
J109T	109	J335H	94	J459H	101
J110T	109	J337H	96	J460H	101
J111T	109	J338H	96	J461H	101
J112T	109	J339H	96	J463G	104
J113T	110	J340H	96	J464G	105
J195H	102	J341H	96	J466H	102
J196H	102	J344H	96	J467H	102
J204G	110	J345H	96	J468H	102
J205G	110	J346H	96	J471H	102
J206G	110	J347H	96	J472H	102
J207G	110	J351H	97	J473H	102
J208G	110	J352H	97	J474H	102
J212H	93	J353H	97	J479H	102
J213H	93	J356H	97	J480H	102
J214H	93	J357H	97	J485H	102
J215H	93	J358H	97	J486H	102
J218H	94	J359H	97	J488G	104
J219H	94	J363H	98	J489G	104
J232H	94	J364H	98	J490G	104
J258H	96	J365H	98	J492G	104
J259H	96	J369H	98	J492H	104
J260H	96	J370H	98	J493G	104
J261H	96	J371H	98	J493H	104
J266H	102	J375H	99	J494G	104
J267H	102	J376H	99	J494H	105
J268H	102	J381H	99	J495G	105
J269H	94	J385H	100	J495H	105
J270H	94	J386H	100	J496G	105
J275H	97	J391H	100	J496H	105
J280H	97	J392H	100	J497G	105
J281H	97	J393H	100	J497H	105
J283G	110	J401G	93	J499G	105
J284G	110	J405G	93	J500G	105
J285G	110	J408G	99	J500H	105
J286G	110	J409G	99	J502H	104
J287G	110	J415H	95	J503G	105
J303H	93	J416H	95	J504G	105
J304H	93	J417H	95	J506G	106
J305H	93	J418H	95	J507G	106
J306H	93	J421H	100	J510G	105
J310H	94	J422H	100	J511G	105
J311H	94	J423H	100	J516H	99
J315H	95	J426H	105	J517H	99
J316H	95	J427H	105	J518H	99
J317H	95	J428H	105	J519H	99
		J430T	96	J523H	104
		J433H	93	J527H	95
		J434H	93	J531H	103
		J436H	97	J532H	103

J533H	102	J662H	104	J789D	94
J534H	102	J663H	104	J790D	101
J535H	102	J670G	108	J801D	95
J544G	108	J671G	108	J808T	106
J546G	108	J682G	105	J813G	107
J547G	108	J682H	105	J814G	107
J548G	108	J683G	105	J822G	107
J551G	104	J683H	105	J823G	107
J552G	108	J693H	101	J823H	107
J553G	108	J694H	101	J824G	107
J555G	107	J695H	101	J833G	106
J556G	107	J698H	103	J834G	106
J562G	107	J699H	103	J835G	106
J566G	108	J700D	97	J838D	96
J568H	103	J701D	97	J839D	96
J569H	103	J702D	97	J840D	96
J570G	108	J707T	101	J841D	96
J571G	108	J708T	101	J843G	106
J574G	108	J710T	102	J844G	106
J575G	108	J711T	102	J845G	106
J576G	108	J713D	93	J849G	98
J578G	108	J714D	93	J864D	95
J579G	108	J716T	103	J868H	102
J581G	98	J717T	103	J869H	102
J582G	98	J718T	103	J870H	102
J583G	98	J719T	103	J871H	102
J588H	101	J723D	98	J880T	98
J589H	101	J724D	98	J885G	110
J590G	108	J726D	94	J904T	109
J591G	108	J727D	94	J905T	109
J593G	105	J738D	96	J906G	109
J594G	105	J739D	96	J907G	110
J596H	104	J740D	96	J910T	109
J599G	108	J741D	96	J911T	109
J602H	99	J743D	93	J912G	109
J603H	99	J748T	101	J913G	110
J606H	109	J749T	101	J916H	99
J607H	109	J750T	102	J917H	99
J608H	109	J751D	97	J935H	105
J610H	104	J752D	97	J936H	105
J614H	109	J753D	97	J940H	103
J615H	109	J754T	102	J944H	96
J616H	109	J756T	103	J945H	96
J617H	110	J757T	103	J946H	96
J618H	110	J761D	102	J947H	96
J633H	109	J762D	102	J956H	97
J634H	109	J764D	98	J957H	97
J635H	109	J765D	98	J958H	97
J636H	109	J771D	94	J959H	97
J643H	109	J772D	94	J974G	108
J644H	109	J773D	95	J975G	108
J645H	109	J774D	95	J977H	98
J646H	109	J775D	95	J978H	98
J649G	98	J783G	95	J979H	98
J650G	98	J784G	95	J980H	98
J656G	106	J785G	95	JB259	92
J657G	106	J787D	101	JB260	92

JB340.....	92	Y378H.....	117	Y935H.....	119
JB341.....	92	Y397H.....	117	Y936H.....	119
JB724.....	92	Y398H.....	117	Y944H.....	116
JB725.....	92	Y399H.....	117	Y945H.....	116
JB839.....	92	Y415H.....	115	Y946H.....	116
JB840.....	92	Y416H.....	115	Y947H.....	117
JB841.....	92	Y417H.....	115	Y967H.....	118
JB864.....	92	Y426H.....	119	YB257.....	114
JB944.....	92	Y427H.....	119	YB258.....	114
JB945.....	92	Y432H.....	114	YB416.....	114
JB946.....	92	Y433H.....	114	YB944.....	114
JB947.....	92	Y434H.....	114	YB946.....	114
JB977.....	92	Y436H.....	117	YB947.....	114
JB978.....	92	Y437H.....	117	YB978.....	114
JJ31G.....	96	Y463G.....	119	YY30G.....	117
JJ40G.....	96	Y464G.....	119	YY31G.....	117
JJ41G.....	96	Y489G.....	119	YY70G.....	117
JJ42G.....	96	Y490G.....	119	YY71G.....	117
JJ80G.....	95	Y492G.....	119		
JJ81G.....	95	Y493G.....	119		
		Y494G.....	119		
		Y495G.....	119		
		Y496G.....	119		
		Y497G.....	119		
		Y510G.....	119		
		Y513G.....	119		
		Y523H.....	118		
		Y527H.....	115		
		Y604H.....	117		
		Y605H.....	117		
		Y606H.....	117		
		Y682H.....	119		
		Y701D.....	117		
		Y732H.....	115		
		Y738D.....	116		
		Y739D.....	116		
		Y740D.....	116		
		Y752D.....	117		
		Y762H.....	116		
		Y801D.....	116		
		Y814G.....	120		
		Y817H.....	118		
		Y822G.....	120		
		Y823G.....	120		
		Y833G.....	119		
		Y834G.....	119		
		Y835G.....	119		
		Y844G.....	120		
		Y845G.....	120		
		Y872D.....	115		
		Y908G.....	120		
		Y909G.....	120		
		Y910G.....	120		
		Y911G.....	120		
		Y917H.....	118		
		Y931H.....	118		
		Y932H.....	118		
				PDS® II (polydioxanone) Suture.....	121-129
				Z00JG.....	126
				Z114H.....	123
				Z117H.....	122
				Z123H.....	122
				Z126H.....	122
				Z127H.....	122
				Z135H.....	122
				Z137H.....	125
				Z148H.....	123
				Z149H.....	123
				Z155H.....	122
				Z195T.....	127
				Z259H.....	124
				Z303H.....	123
				Z304H.....	123
				Z305H.....	123
				Z310H.....	123
				Z311H.....	123
				Z315H.....	124
				Z316H.....	124
				Z317H.....	124
				Z320H.....	123
				Z332H.....	123
				Z333H.....	123
				Z334H.....	123
				Z338H.....	124
				Z339H.....	124
				Z340H.....	124
				Z341H.....	124
				Z346H.....	124
				Z347H.....	124
				Z352H.....	124
				Z353H.....	124
				Z357H.....	124
				Z358T.....	124
				Z359T.....	124
				Z370T.....	125

MONOCRYL®**(poliglecaprone 25) Suture 113-120**

Z371T.....	125	Z877G.....	125	4923H.....	133
Z397H.....	126	Z879G.....	125	4924H.....	133
Z398H.....	126	Z880G.....	125	4925H.....	133
Z416H.....	124	Z881G.....	125	4995H.....	132
Z420G.....	125	Z968H.....	127	5589T.....	133
Z421H.....	126	Z969H.....	127	598H.....	139
Z422H.....	126	Z970H.....	127	634G.....	137
Z423H.....	126	Z990G.....	125	635H.....	137
Z432H.....	123	Z991G.....	125	636H.....	137
Z441H.....	126	Z997G.....	126	637H.....	137
Z442H.....	126	ZB340.....	122	654H.....	138
Z443H.....	126	ZB345.....	122	687G.....	139
Z451G.....	126	ZB346.....	122	723G.....	137
Z451H.....	126	ZB347.....	122	724G.....	137
Z452G.....	126	ZB370.....	122	744G.....	137
Z452H.....	126	ZB371.....	122	751G.....	137
Z463G.....	128	ZB740.....	122	752G.....	137
Z466H.....	127			790G.....	139
Z467H.....	127	Surgical Gut Suture - Chromic.....	130-142	792G.....	139
Z468H.....	127	1621H.....	139	793G.....	139
Z487G.....	128	1635G.....	139	794G.....	139
Z489G.....	128	1636G.....	139	795G.....	139
Z493G.....	128	1637G.....	139	796G.....	139
Z494G.....	128	1638G.....	139	798G.....	140
Z495G.....	128	1641G.....	140	800H.....	134
Z496G.....	128	1642G.....	140	801H.....	134
Z497G.....	128	1643G.....	140	802H.....	134
Z503G.....	128	1652G.....	140	803H.....	134
Z504G.....	128	1653G.....	140	805H.....	138
Z507G.....	128	1654G.....	139	810H.....	133
Z513G.....	128	1658G.....	140	811H.....	133
Z534T.....	127	1659G.....	140	812H.....	133
Z535H.....	127	1731G.....	141	813H.....	133
Z569T.....	127	1744G.....	141	814H.....	133
Z583G.....	125	1745G.....	141	816H.....	138
Z662H.....	127	1754G.....	141	817H.....	138
Z663H.....	127	1757G.....	141	832H.....	138
Z682G.....	128	1766G.....	141	833H.....	138
Z683G.....	128	1771G.....	141	837H.....	136
Z690G.....	125	1791G.....	141	838H.....	136
Z695H.....	126	1792G.....	141	839H.....	136
Z739D.....	124	1793G.....	141	862H.....	135
Z740D.....	124	1797G.....	141	863H.....	135
Z741D.....	124	1798G.....	141	864H.....	135
Z761D.....	127	1816G.....	140	865H.....	135
Z762D.....	127	2643G.....	140	882H.....	132
Z765D.....	125	2644G.....	140	883H.....	132
Z773D.....	124	3802T.....	134	884H.....	132
Z774D.....	124	3803T.....	134	885H.....	132
Z775D.....	124	3811T.....	133	886H.....	137
Z776D.....	125	3812T.....	133	887H.....	137
Z807T.....	129	3813G.....	133	888H.....	137
Z823G.....	129	44T.....	135	892H.....	132
Z824G.....	129	47T.....	131	893H.....	132
Z833G.....	129	4887G.....	131	903H.....	135
Z844G.....	129	48G.....	131	904H.....	135
Z845G.....	129	4915H.....	134	905H.....	135

912H.....	134	SG13T	142	H821H.....	145
913H.....	134	SG14T	142	H822H.....	145
914H.....	134	SG15T	142	H823H.....	145
915H.....	134	U202H.....	131	L101G.....	148
922H.....	133	U203H.....	131	L102G.....	148
923H.....	133	U204H.....	131	L103G.....	148
924H.....	133	U213H.....	137	L104G.....	148
925H.....	133	U233H.....	136	N862H.....	144
942H.....	136	U245H.....	135	N863H.....	144
943H.....	136	U246H.....	135	N864H.....	144
944H.....	136	U247H.....	135	S102H.....	148
945H.....	136	U255H.....	136	S103H.....	148
975H.....	138	U256H.....	136	S104H.....	148
977H.....	138			S105H.....	148
982H.....	138	Surgical Gut Suture - Plain.....143-148		SG2T.....	148
983H.....	138	1626G.....	147	U207H.....	144
990H.....	134	1627G.....	147		
991H.....	134	1630H.....	147	VICRYL RAPIDE™	
CC01G.....	133	1632G.....	147	(polyglactin 910) Suture..... 87-90	
CC02G.....	133	1633G.....	147	VR1695.....	88
CC03G.....	132	1640H.....	147	VR214.....	88
CC04G.....	132	1644G.....	147	VR2218.....	88
CC38G.....	133	1657G.....	146	VR416.....	88
CC40G.....	133	1735G.....	148	VR417.....	88
CC41G.....	133	1772G.....	148	VR426.....	89
G121H.....	132	1775G.....	148	VR490.....	89
G122H.....	132	1824H.....	146	VR493.....	89
G123H.....	132	1828H.....	146	VR494.....	89
G124H.....	132	1915G.....	147	VR496.....	89
G126H.....	133	52T.....	144	VR497.....	89
G127H.....	133	53T.....	144	VR834.....	90
G128H.....	133	5503T.....	144	VR835.....	90
G133H.....	134	612G.....	145	VR844.....	90
G152H.....	136	658G.....	147	VR845.....	90
G172H.....	132	665H.....	146	VR916.....	88
G173H.....	132	686G.....	147	VR917.....	88
G181H.....	132	742G.....	145	VR931.....	89
G182H.....	132	770G.....	146	VR932.....	89
GL30G.....	135	772G.....	146	VR935.....	89
GL31G.....	135	773G.....	146	VR944.....	88
K895H.....	131	774G.....	146	VR945.....	88
L111G.....	142	775G.....	146		
L112G.....	142	842H.....	144	VICRYL® (polyglactin 910) Suture..... 111-112	
L113G.....	142	843H.....	144	V402G.....	112
L114G.....	142	844H.....	144	V448G.....	112
L115G.....	142	846H.....	145	V449G.....	112
M760G.....	135	849G.....	145	V450G.....	112
N877H.....	135	852H.....	144	V549G.....	112
N878H.....	135	853H.....	144	V950G.....	112
N879H.....	135	872H.....	144	V960G.....	112
S111H.....	142	873H.....	144	V966G.....	112
S112H.....	142	898H.....	145	V980G.....	112
S113H.....	142	CP53G.....	144		
S114H.....	142	G321H.....	144		
S115H.....	142	G322H.....	144		
S116H.....	142	H810H.....	145		
S117H.....	142	H820G.....	145		

Sutures: Nonabsorbable149-222**ETHIBOND EXCEL®****Polyester Suture.....155-166**

10X42.....	157
10X52.....	162
10X53.....	162
10X66.....	163
10X72.....	163
10X82.....	159
10X86.....	159
B409T.....	164
B499T.....	164
CX12D.....	158
CX13D.....	158
CX14D.....	158
CX21D.....	159
CX22D.....	159
CX26D.....	158
CX27D.....	158
CX30D.....	160
CX31D.....	160
CX41D.....	157
CX42D.....	157
CX45D.....	158
CX46D.....	158
CX74T.....	164
CX75T.....	164
MB46G.....	161
MB47G.....	166
MB66G.....	163
MPX17.....	162
MPX37.....	162
MX523.....	157
MX551.....	157
MX552.....	157
MX553.....	157
MX563.....	159
MX69G.....	163
MX763.....	158
MX823.....	158
MX833.....	158
MX927.....	161
MX944.....	161
MX977.....	162
PX17H.....	162
PX514.....	159
PX54H.....	162
PX62H.....	163
PX64H.....	163
PX65H.....	163
PX82H.....	157
PX83H.....	157
PX937.....	162
PX977.....	162
PXX22.....	161
PXX41.....	157
PXX43.....	157
PXX50.....	162
PXX52.....	162
PXX76.....	159
PXX77.....	163
PXX79.....	163
PXX80.....	157
PXX82.....	157
PXX85.....	159
PXX86.....	159
PXX87.....	159
PXX92.....	158
SX17H.....	162
SX54H.....	162
SXP02.....	159
SXP07.....	162
SXP22.....	161
SXP52.....	162
SXP54.....	162
SXP76.....	159
SXP77.....	163
SXP82.....	157
SXP86.....	159
SXX43.....	157
SXX50.....	162
SXX53.....	162
SXX54.....	162
SXX86.....	159
X114H.....	156
X182H.....	166
X183H.....	166
X184H.....	166
X185H.....	166
X186H.....	166
X190G.....	160
X204H.....	156
X205H.....	156
X303H.....	166
X304H.....	166
X305H.....	166
X306H.....	166
X406T.....	164
X411H.....	158
X412H.....	158
X423H.....	159
X424H.....	159
X425H.....	159
X444H.....	161
X496T.....	164
X513H.....	159
X517H.....	164
X518H.....	164
X519H.....	164
X521H.....	159
X522H.....	159
X523H.....	159
X524H.....	159
X538H.....	165
X548H.....	165
X550H.....	157
X551H.....	157
X556H.....	157
X557H.....	157
X558H.....	157
X559H.....	157
X562H.....	159
X563H.....	159
X581H.....	156
X582H.....	157
X588H.....	156
X622H.....	165
X663H.....	164
X664H.....	164
X681G.....	165
X683G.....	163
X692G.....	165
X695G.....	165
X698G.....	165
X706H.....	156
X717H.....	156
X726H.....	156
X727H.....	156
X728H.....	160
X761H.....	158
X762H.....	158
X763H.....	158
X806H.....	162
X807H.....	162
X831H.....	158
X832H.....	158
X833H.....	158
X834H.....	158
X843H.....	159
X844H.....	159
X864H.....	160
X865H.....	160
X870H.....	157
X871H.....	157
X872H.....	157
X873H.....	157
X881H.....	156
X889H.....	156
X890H.....	156
X905H.....	163
X916H.....	162
X917H.....	162
X925H.....	161
X926H.....	161
X927H.....	161
X934H.....	162
X935H.....	162
X936H.....	162
X937H.....	162
X942H.....	160

X951H.....	161	1993G.....	174	7711G.....	175
X963H.....	161	1994G.....	174	7717G.....	175
X964H.....	161	1995G.....	174	7718G.....	175
X966H.....	163	2800G.....	168	7731G.....	174
X975H.....	162	2808G.....	168	7756G.....	175
X976H.....	162	2809G.....	168	7757G.....	175
X977H.....	162	2810G.....	168	7760G.....	175
X987H.....	162	2813G.....	168	7770G.....	175
X997G.....	163	2815G.....	168	824G.....	169
		2819G.....	168	825G.....	169
ETHILON® Nylon Suture.....	167-176	2820G.....	168	830H.....	169
1603G.....	173	2822G.....	168	9000G.....	176
1611G.....	172	2829G.....	168	9001G.....	176
1629H.....	170	2830G.....	168	9002G.....	176
1647G.....	172	2849G.....	169	9003G.....	176
1660G.....	173	2850G.....	168	9005G.....	176
1661G.....	173	2860G.....	168	9006G.....	176
1662G.....	173	2870G.....	168	9007G.....	176
1663G.....	173	2881G.....	168	9010G.....	176
1663H.....	173	2890G.....	174	9011G.....	176
1665G.....	172	2892G.....	169	9013G.....	176
1666G.....	172	2894G.....	169	9016G.....	176
1666H.....	172	2897G.....	169	9030G.....	176
1667G.....	172	2898G.....	169	9031G.....	176
1667H.....	172	460T.....	171	9032G.....	176
1668G.....	172	470G.....	171	9033G.....	176
1668H.....	172	489T.....	171	9036G.....	176
1669H.....	172	490T.....	171	9038G.....	176
1670H.....	170	562H.....	170	9060G.....	175
1671H.....	170	585H.....	173	9061G.....	176
1673H.....	170	593H.....	172	9065G.....	175
1674H.....	170	626H.....	171	G666G.....	172
1683H.....	171	627H.....	171	G667G.....	172
1691H.....	170	628H.....	171	G695G.....	172
1696G.....	172	642G.....	171	G696G.....	172
1697H.....	171	661G.....	170	G697G.....	172
1698G.....	172	661H.....	170	G698G.....	172
1714G.....	174	662G.....	170	L880G.....	169
1716G.....	175	662H.....	170	L886T.....	168
1740G.....	175	663G.....	170		
1845G.....	174	663H.....	170	MERSILENE® Polyester	
1854G.....	173	664G.....	170	Fiber Suture.....	177-181
1855G.....	173	664H.....	170	1722G.....	180
1856G.....	173	667G.....	169	1760G.....	180
1864G.....	173	668G.....	170	1761G.....	180
1865G.....	173	669H.....	170	1762G.....	180
1866G.....	173	689G.....	172	1764G.....	180
1893G.....	174	690G.....	172	1778G.....	180
1894G.....	174	691G.....	172	1779G.....	180
1895G.....	174	697G.....	172	740G.....	179
1945G.....	174	697H.....	172	746G.....	179
1954G.....	173	698G.....	172	MR21T.....	178
1955G.....	173	698H.....	172	R424H.....	178
1956G.....	173	699G.....	172	R495T.....	179
1964G.....	173	699H.....	172	R633H.....	178
1965G.....	173	749G.....	174	R647H.....	178
1966G.....	173	7707G.....	175	R665H.....	178

R690G	179	1780G	194	A184H	196
R691G	179	1925G	194	A185H	196
R749G	180	1984G	194	A186H	196
R750G	180	413H	188	A187H	196
R756G	180	423H	189	A302H	196
R770G	180	424H	189	A303H	196
R832H	178	443H	192	A304H	196
R833H	178	482T	192	A305H	196
RS22	181	485T	192	A306H	196
RS23	181	486H	194	A307H	196
		580H	194	CO03D	188
NUROLON® Nylon Suture	182-185	583H	193	CO04D	188
5424H	184	590H	194	CO125	189
5425H	184	591H	192	CO12D	189
5665G	185	621H	193	CO135	189
5666G	185	622H	193	CO13D	189
C502D	183	623H	193	CO145	189
C503D	183	624H	193	CO14D	189
C512D	183	629G	191	CO16D	189
C513D	183	632G	192	CO17D	189
C514D	183	639G	193	CO18D	189
C520D	184	640G	193	CO41D	188
C521D	184	641G	193	CO45D	189
C522D	184	673H	194	CO53D	188
C526D	183	675H	191	CO54D	188
C527D	183	677G	191	CO84D	188
C541D	183	678G	191	CO21D	189
C542D	183	678H	191	CO22D	189
C545D	183	679H	191	CO26D	188
C546D	183	680H	191	CO27D	188
C550D	184	682G	191	K533H	189
C551D	184	683G	191	K534H	189
C553D	183	684G	191	K570H	188
C554D	183	684H	191	K571H	188
C563D	184	685G	191	K572H	188
C566D	184	685H	191	K580H	187
C570D	184	706G	187	K581H	187
C572D	185	707G	187	K801H	187
C575D	185	709H	187	K802H	187
C576D	185	714G	191	K809H	187
C584D	183	734G	192	K831H	189
C595T	185	735G	191	K832H	189
N104T	183	736G	191	K833H	189
N124T	183	737G	192	K834H	189
		765G	192	K835H	189
PERMA-HAND® Silk Suture	186-196	766G	192	K843H	190
1588H	194	768G	192	K844H	190
1639G	193	769G	192	K845H	190
1675H	193	7729G	191	K846H	190
1676G	193	7733G	194	K852H	190
1677G	193	7734G	191	K870H	188
1678H	193	780G	193	K871H	188
1679H	193	783G	193	K872H	188
1684G	193	786G	193	K873H	188
1687G	194	789G	193	K880H	187
1689H	191	A182H	196	K881H	187
1732G	194	A183H	196	K882H	187

K883H	187	7761G	216	8621H	213
K889H	187	7796G	216	8622H	213
K890H	187	788G	210	8623H	213
K891H	187	8190G	210	8629H	213
K945G	190	8204H	207	8630G	215
K952H	190	8205H	207	8631G	215
K964H	190	8206H	206	8632G	215
LA53G	195	8301H	203	8634G	215
LA54G	195	8304H	203	8635G	215
LA55G	195	8305H	203	8636G	215
M104T	188	8307H	205	8648G	214
M13T	189	8318H	205	8649H	213
M3T	188	8321H	205	8655G	214
M72T	190	8325H	205	8656G	214
N266H	188	8356H	207	8661G	212
N267H	187	8357H	207	8663G	214
N271H	190	8381H	206	8665G	212
N272H	190	8410H	208	8675H	213
Q706H	187	8411H	208	8680G	214
Q809H	187	8412H	208	8681G	214
SA10G	195	8417H	208	8682G	214
SA11G	195	8418H	208	8682H	214
SA5H	195	8422H	209	8683G	212
SA63H	195	8423H	209	8684G	213
SA64H	195	8424H	209	8685H	213
SA65H	195	8425H	209	8686G	214
SA66G	195	8433H	209	8687H	214
SA6H	195	8434H	209	8689H	213
SA73H	195	8435H	209	8690H	213
SA74H	195	8444H	212	8695G	214
SA75H	195	8454H	209	8696G	214
SA76G	195	8455H	209	8697G	214
SA77G	195	8519H	206	8698G	214
SA7H	195	8521H	208	8699G	214
SA83H	195	8522H	208	8701H	203
SA84H	195	8523H	208	8702H	203
SA85H	195	8526H	208	8703H	203
SA86G	195	8533H	208	8704H	211
SA87G	195	8534H	208	8705H	210
SA8H	195	8554H	206	8706H	205
SA9G	195	8555H	207	8707H	211
		8556H	207	8708H	204
		8557H	207	8709H	203
		8558H	207	8710H	206
		8559H	207	8711H	206
		8571G	213	8713H	206
		8580H	206	8714H	206
		8581H	206	8717H	205
		8592G	212	8718H	205
		8603G	214	8719H	211
		8604G	214	8720H	205
		8605G	214	8721H	211
		8606G	214	8725H	205
		8610H	204	8726H	205
		8617G	215	8730H	198
		8618G	215	8732H	199
		8619G	215	8734H	199
PROLENE® Polypropylene Suture ... 197-216					
1710G	216				
1713G	213				
1754G	215				
1755G	216				
1757G	215				
1771G	215				
2775G	198				
2790G	198				
2793G	198				
2794G	198				
2795G	198				
3844G	210				
3846T	213				
7740G	215				

8735H.....	200	C841G.....	207	EPM8731.....	200
8740H.....	199	C845G.....	208	EPM8732.....	199
8741H.....	199	EP8610H.....	204	EPM8733.....	200
8746H.....	201	EP8706H.....	205	EPM8734.....	199
8747H.....	201	EP8707H.....	211	EPM8735.....	200
8752H.....	201	EP8708H.....	204	EPM8736.....	199
8753H.....	201	EP8709H.....	203	EPM8737.....	200
8754H.....	201	EP8711H.....	206	EPM8738.....	198
8755H.....	202	EP8714H.....	206	EPM8739.....	199
8756H.....	207	EP8718H.....	205	EPM8741.....	199
8757H.....	207	EP8719H.....	211	EPM8742.....	199
8762H.....	208	EP8726H.....	205	EPM8743.....	199
8766H.....	200	EP8730H.....	198	EPM8745.....	202
8776H.....	204	EP8732H.....	199	EPM8746.....	201
8800H.....	204	EP8734H.....	199	EPM8747.....	202
8805H.....	203	EP8735H.....	200	EPM8751.....	202
8806H.....	203	EP8740H.....	199	EPM8752.....	201
8807H.....	211	EP8741H.....	199	EPM8753.....	201
8824G.....	209	EP8746H.....	201	EPM8754.....	201
8825G.....	209	EP8747H.....	201	EPM8755.....	202
8831H.....	208	EP8752H.....	201	EPM8765.....	200
8832H.....	208	EP8753H.....	201	EPM8766.....	200
8833H.....	208	EP8754H.....	201	EPM8767.....	200
8834H.....	208	EP8755H.....	202	EPM8768.....	200
8842H.....	209	EP8766H.....	200	EPM8775.....	200
8843H.....	209	EP8776H.....	204	EPM8776.....	204
8845G.....	210	EP8800H.....	204	EPM8785.....	200
8853H.....	209	EP8805H.....	203	EPM8800.....	204
8860H.....	206	EP8806H.....	203	EPM8805.....	203
8861H.....	206	EP8807H.....	211	EPM8806.....	203
8862H.....	206	EP8889H.....	204	EPM8807.....	211
8870H.....	207	EPH7477H.....	205	EPM8812SL.....	212
8871H.....	207	EPH8710H.....	206	HS6822H.....	208
8872H.....	207	EPH8716H.....	206	HS6855H.....	208
8881H.....	206	EPH8720H.....	205	HS6856H.....	207
8889H.....	204	EPH8725H.....	205	HS6857H.....	207
8890H.....	204	EPH8890H.....	204	HS6861H.....	207
8925H.....	211	EPH9702H.....	203	HS8556H.....	207
8926H.....	211	EPM8203.....	203	M8203.....	203
8934H.....	212	EPM8205.....	203	M8205.....	203
8935H.....	212	EPM8206.....	205	M8206.....	205
8936H.....	212	EPM8610.....	204	M8220.....	205
8937H.....	212	EPM8701.....	203	M8301.....	203
8975H.....	212	EPM8702.....	203	M8303.....	203
8976H.....	212	EPM8703.....	203	M8304.....	203
8977H.....	212	EPM8704.....	211	M8305.....	203
8978H.....	212	EPM8704SL.....	212	M8307.....	205
9090G.....	216	EPM8705.....	210	M8318.....	205
9091G.....	215	EPM8706.....	205	M8321.....	205
9556H.....	207	EPM8707.....	211	M8325.....	205
9702H.....	203	EPM8709.....	203	M8556.....	207
9706H.....	206	EPM8711.....	206	M8557.....	207
C811G.....	209	EPM8718.....	205	M8600.....	204
C821G.....	209	EPM8719.....	211	M8610.....	204
C822G.....	209	EPM8726.....	205	M8701.....	203
C826G.....	208	EPM8727.....	204	M8702.....	203
C827G.....	208	EPM8730.....	198	M8703.....	203

M8704	211	3307H	152	DS20	222
M8705	210	3321H	152	DS22	222
M8706	205	3521H	153	DS23	222
M8707	211	3522H	153	DS24	222
M8709	203	3523H	153	DS25	222
M8710	206	3526H	153	DS26	222
M8711	206	3555H	153	DS28	222
M8717	205	3556H	153	DS30	222
M8718	205	3557H	153	DS32	222
M8719	211	3558H	153	DS35	222
M8720	205	3580H	153	M400G	219
M8721	211	3581H	153	M414G	218
M8725	205	3701H	152	M435G	218
M8726	205	3702H	152	M436G	220
M8727	204	3703H	152	M645G	221
M8730	198	3706H	152	M649G	218
M8731	200	3707H	154	M650G	218
M8732	199	3709H	152	M651G	218
M8733	200	3710H	153	M652G	220
M8734	199	3711H	153	M653G	220
M8735	200	3717H	152	M654G	220
M8736	199	3720H	152	M655G	220
M8737	200	3725H	152	M656G	218
M8738	198	3726H	152	M657G	220
M8739	199	3735H	151	M660G	218
M8741	199	3762H	153	M682G	218
M8742	199	3776H	152	SS26G	221
M8743	199	3805H	152	SS28G	221
M8745	202	3806H	152	SS30G	221
M8746	201	3842H	154	TPW10	220
M8747	202	3975H	154	TPW20	220
M8752	201	3976H	154	TPW30	221
M8753	201	M3702	152	TPW32	221
M8754	201	M3703	152	TPW40	220
M8755	202	M3706	152	TPW42	220
M8764	200	M3709	152	TPW50	221
M8765	200	M3720	152	TPW52	221
M8766	200	M3725	152	TPW62	221
M8767	200	M3735	151	TPW72	221
M8768	200	M3737	151	TPW80	221
M8775	200	M3739	151	TPW90	221
M8776	204	M3743	151	TPW92	221
M8785	200	M3747	151		
M8800	204	M3766	151		
M8805	203	M3767	151		
M8806	203	M3768	151		
M8807	211	M3770	151		
P8522	208	M3805	152		
P8935	212				
PB20H	198	Surgical Stainless Steel Suture217-222			
PB21H	198	467G	219	True-to-Scale Needle Diagram . . . 29-33	
PB30H	198	497G	219	Blunt Point	33
		603G	219	Cutting Edge	30-32
		607G	219	Endoscopy	33
		618G	219	TAPERCUT®	33
		631G	219	Taperpoint	32
		DS18	222	Trocar	33
PRONOVA® Poly					
(Hexafluoropropylene - VDF) Suture .. 150-154					
3304H	152				
3305H	152				

Place an Order

Contact your Ethicon Representative
or call 1-877-ETHICON (384-4266) to order.

More Information

To find the right product for your wound
closure needs, visit Ethicon.com or the
Wound Closure Resource Center.